# RIA 17 $\alpha$ - Hydroxyprogesterone

**REF** IM1452

## TABLE OF CONTENTS

English . . . . .	2	Polski (PL) . . . . .	32
Français (FR) . . . . .	5	Čeština (CZ) . . . . .	35
Deutsch (DE) . . . . .	8	Slovenčina (SK) . . . . .	38
Italiano (IT) . . . . .	11	Türkçe (TR) . . . . .	41
Español (ES) . . . . .	14	Русский (RU) . . . . .	44
Português Portugal (PT-PT) . . . . .	17	Србија (SR) . . . . .	47
Ελληνικά (EL) . . . . .	20	Português Brasil (PT-BR) . . . . .	50
中文 (ZH-CN) . . . . .	23	APPENDIX . . . . .	53
Lietuviškai (LT) . . . . .	26	REFERENCES . . . . .	56
Magyar (HU) . . . . .	29		

---

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# RIA 17 $\alpha$ - Hydroxyprogesterone

**REF** IM1452

## RADIOIMMUNOASSAY FOR THE IN VITRO DETERMINATION OF 17 $\alpha$ -HYDROXYPROGESTERONE IN HUMAN SERUM AND PLASMA

For *in vitro* diagnostic use.

### PRINCIPLE

The radioimmunoassay of 17 $\alpha$ -Hydroxyprogesterone (17 OH-P) is a competition assay. The same kit may be employed for the measurement of 17 OH-P:

either directly in serum or EDTA plasma,

or after the extraction of the sample with ether and evaporation of the solvent, followed by re-suspension in the zero calibrator.

Serum or plasma samples, extracts, controls and calibrators are incubated with 125I-labeled 17 OH-P as a tracer, in antibody-coated tubes. After incubation, the contents of the tubes are rinsed so as to remove unbound 125I-labelled 17 OH-P. The bound radioactivity is then determined in a gamma counter. A calibration curve is established and unknown values are determined by interpolation from the standard curve.

### WARNING AND PRECAUTIONS

#### General remarks:

- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be established with each assay.
- It is recommended to perform the immunoassay in duplicate.
- Each tube must be used only once.

#### Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

#### Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide may react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

#### Materials of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HBsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All serum and plasma samples should be handled as if capable of transmitting hepatitis or AIDS and waste should be discarded according to the country rules.

#### Ethyl ether

Ethyl ether is a volatile and highly flammable organic solvent. Extraction and evaporation must be done in a ventilated hood. Avoid any contact with a flame and do not pipette reagents by mouth.

### GHS HAZARD CLASSIFICATION

Wash Solution (20X) DANGER

H360

May damage fertility or the unborn child.

P201

Obtain special instructions before use.

P280

Wear protective gloves, protective clothing and eye/face protection.

P308+P313

IF exposed or concerned:  
Get medical advice/attention.  
Boric Acid 0.1 - 0.3%  
Sodium Borate Decahydrate 0.1 - 0.3%

**SDS**

Safety Data Sheet is available at [techdocs.beckmancoulter.com](http://techdocs.beckmancoulter.com)

### SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in dry tubes or in tubes with EDTA.
- Separate serum or plasma from cells by centrifugation.
- Serum and plasma samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage keep frozen (at < -18°C, 1 year maximum) after aliquoting so as to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- If samples have concentrations greater than the highest calibrator, they must be diluted into the zero calibrator.

Serum and EDTA plasma values for 30 samples (serum values ranging from 0.24 to 1.88 ng/mL) were compared using the 17 OH-P RIA kit. Results are as follows:

$$[\text{EDTA-plasma}] = 0.9026[\text{serum}] + 0.0283,$$

$$R = 0.9747$$

### MATERIALS PROVIDED

All reagents in the kit are stable until the expiration date indicated on the kit label, when stored at 2-8°C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Storage conditions for reagents after reconstitution or dilution are indicated in paragraph Procedure.

**Anti-17 OH-P antibody-coated tubes: 2 x 50 tubes** (ready-to-use)

**125I-labeled 17 $\alpha$ -Hydroxyprogesterone: one vial (45 mL)** (ready-to-use)

The vial contains <640 kBq, at the date of manufacture, of 125I-labeled 17 OH-P in buffer containing proteins and a dye.

Note: Occasional presence of clotted particles in the tracer does not affect assay performance.

**Calibrators: one 2 mL vial, five 0.5 mL vials** (ready-to-use)

The calibrator vials contain from 0 to approximately 50 ng/mL (0 to approximately 151 nmol/L) of 17 OH-P in human serum with sodium azide (<0.1%). The exact concentration is indicated on each vial label. Calibrators are verified to an internal reference standard.

The zero calibrator (5 mL) may be ordered separately, too (cat. #B23373).

### Control samples: two vials (lyophilized)

The control vial contains 17 OH-P in human serum with sodium azide (<0.1%). The volume for reconstitution is indicated on the vial label. The expected values are in the concentration range indicated on a supplement.

### Wash solution (20x): one 50 mL vial

Concentrated solution has to be diluted before use.

## MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

For the assay:

- Precision micropipet (25 µL and 400 µL).
- Semi-automatic pipets (25 µL, 400 µL and 2 mL).
- Vortex type mixer.
- Horizontal or orbital shaker.
- Aspiration system.
- Gamma counter set for <sup>125</sup>I

For the extraction step (optional):

- Precision micropipette (200 µL)
- Glass pipets (2.5 mL, 5 mL, 10 mL).
- Glass vials (from 6 mL to 15 mL).
- Glass tubes for recovery of ether phase.
- Evaporator (speedvac type) or 37°C water bath.
- Analytical grade ethyl ether.

## PROCEDURE

### Preparation of reagents

Let all the reagents come to room temperature.

### Reconstitution of control samples

The content of the vials is reconstituted with the volume of distilled water indicated on the label. Wait for 10 min following reconstitution and mix gently to avoid foaming before dispensing. Store the reconstituted solutions at 2-8°C until the expiry date of the kit.

### Preparation of wash solution

Pour the content of the vial into 950 mL of distilled water and homogenize. The diluted solution may be stored at 2-8°C until the expiry date of the kit.

### Extraction of samples (optional; see § Limitations)

Note: The extraction must be done in clean glass vials or tubes, pre-rinsed with ethyl ether.

Samples only are extracted before assay; do not extract calibrators.

- Bring samples to room temperature and mix well before starting extraction.
- Number one vial for each sample.
- Place 200 µL of each sample into corresponding vial.
- Add 5 mL of ethyl ether to vials. Stopper carefully.
- Vortex vials vigorously (2 x 1 minute).
- Let stay vials on the table for minimum 5 minutes.
- Take off carefully 2.5 mL of organic phase without contaminating with aqueous phase and place it into numbered glass tubes.
- Evaporate ether phase completely with either evaporator (e.g. speedvac type) or by placing tubes into 37°C water bath.

Note 1: The tubes must be firmly attached to test tubes rack, since after evaporation of the ether; they will be lighter and tend to float off.

At this stage it is possible to stopper the tubes containing the dry extracts and to store them for up to 7 days in the cold (2-8°C) prior to continuing the assay.

- Re-dissolve dry ether extracts in 200 µL of the zero calibrator. Vortex vigorously (30 sec, 2000 rpm), wait for 15 minutes, vortex again.

Note 2: If a larger volume is necessary, instead of taking off 2.5 mL of organic phase, keep vials at <-18°C until aqueous phase freezes, take off carefully 5 mL of organic phase without contaminating with aqueous phase, evaporate it and re-dissolve it in 400 µL of the zero calibrator.

## Immunoassay procedure

Bring all reagents to room temperature before pipeting.

Step 1 Additions	Step 2 Incubation	Step 3 Counting
To antibody coated tubes add successively:  25 µL of calibrator, control, sample or 50 µL of extract and 400 µL of tracer.* Mix.	Incubate 120 minutes at 18 - 25°C with shaking (≥280 rpm).	Aspirate carefully the content of tubes (except the 2 tubes «total cpm»)  Wash twice with 2 mL of wash solution. Aspirate.  Count bound cpm (B) and total cpm (T) for 1 min.

\*Add 400 µL of tracer to 2 additional tubes to obtain total cpm.

## RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of 17 OH-P concentrations in samples measured at the same time as the calibrator.

### Standard curve

The results in the quality control department were calculated using *spline* curve fit with *B/T* or *B/B<sub>0</sub>* on the logit vertical axis and analyte concentration of the calibrators on the logarithmic horizontal axis (ng/mL).

Other data reduction methods may give slightly different results.

Total activity: 179,362 cpm				
Calibrators	17OH-P (ng/mL)	cpm (n=3)	B/T (%)	B/B <sub>0</sub> (%)
0	0	55,975	31.21	100.0
1	0.12	44,882	25.02	80.18
2	0.40	31,190	17.39	55.72
3	1.90	12,557	7.00	22.43
4	12.5	2,515	1.40	4.49
5	50.0	895	0.50	1.60

(Example of standard curve, do not use for calculation)

### Samples

For each sample, locate the B/T or B/B<sub>0</sub> on the vertical axis and read off the corresponding 17 OH-P concentration on the horizontal axis.

To convert ng/mL into nmol/L, multiply results by 3.026.

## EXPECTED VALUES

It is suggested that each laboratory establishes its own normal values. The following values obtained with healthy subjects are indicative only.

Adults	N	Median	Min.	Max.	2.5 <sup>th</sup>	97.5 <sup>th</sup>
					Percentile	Percentile
(ng/mL)						
Men	55	0.93	0.53	2.22	0.55	1.99
Women						
Follicular Phase	111	0.48	0.13	1.67	0.21	1.45
Luteal Phase	112	1.52	0.42	3.20	0.61	2.88
Preovulatory Peak	22	1.39	0.54	2.04	0.55	2.01
Contraception	30	0.37	0.17	1.48	0.18	1.47
Women after menopause	38	0.38	0.13	0.92	0.16	0.79
Pregnancy 1 <sup>st</sup> trimester	45	2.03	0.78	5.75	0.93	3.82
Pregnancy 2 <sup>nd</sup> trimester	44	2.23	0.60	6.86	1.23	3.70

Detail information about expected values for children (sorted according to age and sex) can be found in the data sheet "APPENDIX".

## QUALITY CONTROL

Good laboratory practices imply that control samples must be used regularly to ensure the quality of the results obtained. These samples

must be processed exactly the same way as the assay samples, and it is recommended to analyze their results using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: [imunochem@beckman.com](mailto:imunochem@beckman.com)

## PERFORMANCE CHARACTERISTICS

(For more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

### Sensitivity

**Limit of detection (LoD):** 0.07 ng/mL

The LoD for 17 OH-P is 0.07 ng/mL, determined consistent with guidelines in CLSI document EP17-A2 [1] based on the proportions of false positives ( $\alpha$ ) less than 5% and false negatives ( $\beta$ ) less than 5%; using determinations, with 120 blank and 120 low level samples; and Limit of Blank (LoB) of 0.04 ng/mL.

### Specificity

The antibody used in the immunoassay is specific for 17- $\alpha$ -hydroxyprogesterone.

### Precision

#### Repeatability and within-laboratory precision

The precision of the assay was determined consistent with guidelines in CLSI document EP05-A3 [2]. For repeatability the coefficients of variation were found below or equal to 10.5% for serum samples. For within laboratory precision the coefficients of variation were found below or equal to 12.8% for serum samples.

### Accuracy

#### Linearity

The assay demonstrated to be linear from 0.12 to 50.88 ng/mL using serum samples (determined consistent with guidelines in CLSI document EP06-A [3]).

### Dilution test

High-concentration serum samples were serially diluted with zero calibrator. The recovery percentages ranged from 95.7% to 119%.

### Recovery test

Serum samples were spiked with known quantities of 17 OH-P. The recovery percentages ranged from 100% to 119%.

**Measurement range** (from LoD to the highest calibrator):

0.07 ng/mL to approximately 50 ng/mL

## LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly.

Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

For more details, see the data sheet "APPENDIX".

Very young infants (particularly less than 6 months old) have very high levels of 17OH-pregnenolone sulfate. In that particular case, extraction with ether prior to the 17 OH-P assay is recommended.

In spite of low cross reactivity, one can observe falsely high levels of 17 OH-P after regular administration of spironolactone. Spironolactone therapy must therefore be interrupted at least three weeks prior to the assay to ensure absence of interference.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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#### Interference

Serum samples containing 17- $\alpha$ -hydroxyprogesterone concentrations of app. 0.5 and 3 ng/mL were spiked with multiple concentrations of the substances below and assayed using 17 OH-P RIA kit. Values were calculated as described in CLSI EP07, 3rd ed. [4] and EP37, 1st ed. [5]. Interference was determined by testing controls (no interfering substance added) and matched test samples (with interfering substance added). No interference (defined as a shift in dose >15%) was found for addition of interferent up to concentration stated in the table below.

Interferent	Test concentration
Hemoglobin	10,191 $\mu$ g/mL
Conjugated bilirubin	435.9 $\mu$ g/mL
Unconjugated bilirubin	130.6 $\mu$ g/mL
Biotin	1,705 ng/mL
Ascorbic acid	65.67 $\mu$ g/mL
Acetylsalicylic acid	64.48 $\mu$ g/mL
Ibuprofen	211.5 $\mu$ g/mL
Cholesterol	8.13 mg/mL
Heparin	6,227 ng/mL
Prednisone	282.6 ng/mL
Prednisolone	1,282 ng/mL
Protein ( $\gamma$ -globulin)	136.8 mg/mL
Rheumatoid factor	19.76 IU/mL
TAG	12.32 mg/mL

#### Specificity

Data on cross-reactivity with several hormones are presented in the following table:

$$\text{Cross-reactivity (\%)} = \frac{17 \text{ OH-P concentration}}{\text{Compound concentration}} \times 100 \text{ at } 50\% \text{ binding of the zero calibrator}$$

COMPOUND	% CROSS-REACTIVITY
17- $\alpha$ -hydroxyprogesterone	100.0
11-deoxycortisol	2.24
17- $\alpha$ -hydroxypregnenolone sulfate	1.78
17- $\alpha$ -hydroxypregnenolone	1.00
Progesterone	0.28
11-deoxycorticosterone	ND
Aldosterone	ND
Androstenedione	ND
Corticosterone	ND
Cortisone	ND
DHEA	ND
Estradiol	ND
Estriol	ND
Estrone	ND
Etiocolanalone	ND
Pregnenolone	ND
Spironolactone	ND
Testosterone	ND

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Five serum samples, five EDTA-plasma samples and five extracts were assayed for 20 days, 2 runs per day, triplicates per run. Assays were performed by two lab technicians, by two reagent lots. There were 120 individual measurements per sample with no invalid results.

Serum	Mean (ng/mL)	Repeatability		Within laboratory precision	
		SD (ng/mL)	C.V. (%)	SD (ng/mL)	C.V. (%)
S1	29.03	3.05	10.50	3.73	12.84
S2	15.88	1.19	7.50	1.44	9.05
S3	11.06	0.60	5.41	0.84	7.58
S4	1.00	0.04	4.28	0.07	7.13
S5	0.20	0.02	8.38	0.02	10.30

EDTA-plasma	Mean (ng/mL)	Repeatability		Within laboratory precision	
		SD (ng/mL)	C.V. (%)	SD (ng/mL)	C.V. (%)
P1	26.02	3.36	12.91	4.41	16.94
P2	19.49	1.54	7.91	2.74	14.05
P3	8.80	0.59	6.69	0.85	9.64
P4	0.69	0.05	7.35	0.07	9.98
P5	0.25	0.02	7.54	0.03	11.19

Extract	Mean (ng/mL)	Repeatability		Within laboratory precision	
		SD (ng/mL)	C.V. (%)	SD (ng/mL)	C.V. (%)
E1	30.49	2.63	8.62	5.83	19.13
E2	16.52	1.15	6.96	2.59	15.70
E3	10.73	0.64	5.93	1.45	13.54
E4	0.74	0.04	5.65	0.08	10.18
E5	0.26	0.02	7.11	0.03	11.45

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The assay demonstrated to be linear from 0.15 to 62.04 ng/mL using EDTA-plasma samples (determined consistent with guidelines in CLSI document EP06-A [3]).

The assay demonstrated to be linear from 0.16 to 64.39 ng/mL using extracts (determined consistent with guidelines in CLSI document EP06-A [3]).

##### Dilution test

Three serum / EDTA-plasma samples / extracts were diluted into the zero calibrator and assayed according to assay procedure.

Serum	Dilution factor	17 $\alpha$ -hydroxyprogesterone (ng/mL)		Ratio (%)
		Measured	Expected	Measured/Expected
S1	-	30.00	-	-
	1:2	14.36	15.00	95.73
	1:4	7.45	7.50	99.33
	1:8	4.45	3.75	118.7
	1:16	2.22	1.88	118.4
S2	-	10.78	-	-
	1:2	5.91	5.39	109.6
	1:4	3.00	2.70	111.3
	1:8	1.48	1.35	109.8
	1:16	0.68	0.67	100.9
S3	-	5.83	-	-
	1:2	3.03	2.92	103.9
	1:4	1.61	1.46	110.5
	1:8	0.74	0.73	101.5
	1:16	0.36	0.36	98.80

EDTA-plasma	Dilution factor	17 $\alpha$ -hydroxyprogesterone (ng/mL)		Ratio (%)
		Measured	Expected	Measured/Expected
P1	-	14.73	-	-
	1:2	7.57	7.37	102.8
	1:4	3.93	3.68	106.7
	1:8	2.08	1.84	113.0
P2	1:16	0.99	0.92	107.5
	-	15.04	-	-
	1:2	8.64	7.52	114.9
	1:4	4.29	3.76	114.1
P3	1:8	2.20	1.88	117.0
	1:16	1.08	0.94	114.9
	-	9.02	-	-
	1:2	5.01	4.51	111.1
	1:4	2.59	2.26	114.9
	1:8	1.29	1.13	114.4
	1:16	0.61	0.56	108.2

EDTA-plasma	Endogen. conc. (ng/mL)	Added conc. (ng/mL)	Expected conc. (ng/mL)	Measured conc. (ng/mL)	Ratio (%) Measured/Expected
P1	0.45	0.25	0.70	0.78	112.1
	0.44	0.48	0.92	1.10	119.2
	0.44	1.92	2.37	2.53	107.0
P2	0.58	0.36	0.95	1.03	108.8
	0.57	0.71	1.27	1.51	118.6
	0.57	2.38	2.95	3.39	114.8
P3	1.80	0.93	2.72	2.72	99.92
	1.87	1.92	3.79	3.75	98.98
	1.76	4.55	6.31	5.57	88.29

Three samples were spiked by calibrator with higher concentration. Volume of added calibrator was up to 10% of volume of samples. Samples were extracted and assayed according to the procedure of the kit.

Extract	Endogen. conc. (ng/mL)	Added conc. (ng/mL)	Expected conc. (ng/mL)	Measured conc. (ng/mL)	Ratio (%) Measured/Expected
E1	1.13	0.66	1.78	1.75	98.15
	1.10	1.61	2.72	3.09	113.8
	1.07	3.13	4.19	4.90	116.8
E2	0.76	0.66	1.42	1.54	108.6
	0.75	1.30	2.05	2.29	111.8
	0.73	2.83	3.56	4.21	118.4
E3	0.38	0.29	0.67	0.74	111.1
	0.38	0.82	1.20	1.28	106.4
	0.38	1.92	2.30	2.57	111.8

#### Comparison of direct and indirect procedure

30 serum samples were assayed by direct and indirect procedure (sample values in direct procedure ranging from 0.23 to 2.22 ng/mL) using the 17 OH-P RIA kit. Results are as follows:

$$[\text{indirect procedure}] = 0.8448[\text{direct procedure}] + 0.0088,$$

$$R = 0.9633$$

#### Expected data for children

Results are sorted according to the age and sex.

Children	N	Median	Min.	Max.	2.5 <sup>th</sup>	97.5 <sup>th</sup>
					percentile	percentile
(ng/mL)						
0-2 months (after extraction)	40	1.41	0.23	3.08	0.42	2.91
3-5 months (after extraction)	20	0.67	0.33	1.69	0.33	1.68
6-23 months	26	0.72	0.13	2.51	0.14	2.35
Boys 2-11 years	38	0.31	0.09	2.21	0.14	1.41
Girls 2-9 years	32	0.50	0.16	2.41	0.19	1.63
Boys 12-15 years	17	0.67	0.30	2.18	0.32	2.10
Girls 10-15 years	25	0.86	0.33	3.21	0.42	2.64

#### <sup>125</sup>I Characteristics

$$T_{1/2} (^{125}\text{I}) = 1443 \text{ h} = 60.14 \text{ d}$$

<sup>125</sup> I	E (MeV)	%
Y	0.035	
X	0.027	114
	0.032	25

Extract	Dilution factor	17 $\alpha$ -hydroxyprogesterone (ng/mL)		Ratio (%)
		Measured	Expected	Measured/Expected
E1	-	22.50	-	-
	1:2	10.86	11.25	96.53
	1:4	6.16	5.63	109.5
	1:8	3.01	2.81	107.0
	1:16	1.44	1.41	102.4
E2	-	15.26	-	-
	1:2	7.33	7.63	96.07
	1:4	3.75	3.82	98.30
	1:8	1.82	1.91	95.41
	1:16	0.88	0.95	92.27
E3	-	10.11	-	-
	1:2	4.88	5.06	96.54
	1:4	2.54	2.53	100.5
	1:8	1.29	1.26	102.1
	1:16	0.64	0.63	101.3

#### Recovery test

Three serum / EDTA-plasma samples were spiked by calibrator with higher concentration. Volume of added calibrator was up to 10% of volume of samples. Samples were then assayed.

Serum	Endogen. conc. (ng/mL)	Added conc. (ng/mL)	Expected conc. (ng/mL)	Measured conc. (ng/mL)	Ratio (%) Measured/Expected
S1	0.76	0.36	1.12	1.13	100.8
	0.74	0.60	1.34	1.38	103.1
	0.74	2.38	3.12	3.24	103.7
S2	1.32	0.60	1.92	2.08	108.4
	1.34	1.92	3.26	3.43	105.2
	1.29	3.70	4.99	5.95	119.2
S3	1.82	0.82	2.64	2.64	99.99
	1.88	1.92	3.80	4.44	116.9
	1.77	4.55	6.32	7.06	111.7

## Symbols Key

**REF** Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferens / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Úrün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

**IVD** In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / За ин витро диагностика / 體外診斷

**CONTENTS** Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 組成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

 Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 製造商 / Gamintojas / Gyártó: / Producent / Výrobce / Výrobca / 제조 / Üretici / Изготовлено / Произведено от / 製造商

 Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Räcker till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

**CE** CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Mercado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

**SDS** Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

 Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문의 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明

 Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度範圍 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sıcaklık aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

 Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Upozornění / Upozornenie / 주의 / Dikkat / Внимание / 注意

 Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirație / 만료 날짜 / Son Kullanma Tarihi / Срок годности / Срок на годност / 到期日

**LOT** Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer serii / Číslo šarže / 로트 번호 / Lot Numarası / Номер партии / Номер на партида / 批號

 Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Datum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biyolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktiv / Ραδιενεργό / 放射性 / Radioaktyvioji medžiaga / Radioaktiv / Radioaktywny / Radioaktivní / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

**DANGER**

DANGER / DANGER / GEFAHR / PERICOLO / PELIGRO / PERIGO / FARA / ΚΙΝΔΥΝΟΣ / 危險 / PAVOJUS / VESZÉLY / NIEBEZPIECZENSTWO / NEBEZPEČÍ / NEBEZPEČENSTVO / 위험 / TEHLÍKE / ОПАСНО / ОПАСНОСТ / 危險

**Ag**<sup>125I</sup>

**Ab**<sup>125I</sup>

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / Tracer / Ανιχνευτής / 追踪剂 / Atsekamoji medžiaga / Nyomjelző / Znacznik / Radioindikátor / Indikátor (tracer) / 트레이서 / Tracer'lar / метка / Индикатор / 追蹤劑

**CAL**

**CAL 0**

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정 물질 / Kalibratör / Калибратор / Калибратор / 校正液

**CTRL**

Control / Contrôle / Kontrolle / Controllo / Control / Controllo / Kontrolle / Μάρτυρας / 质控品 / Kontrolliné / Kontroll / Kontrola / Kontrola / Kontrola / 정도관리 / Kontrol / Контроль / Контролна / 質控品

**TUBE**

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνάκια / 试管 / Mégintüveliai / Csövek / Probówki / Zkumavky / Skúmavky / 튜브 / Tüpler / пробирки / Епруветки / 試管

**IFU**

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用说明

**SOLN|WASH|20x**

Wash Solution Concentrate 20X / Solution de lavage concentrée 20X / Waschlösungskonzentrat 20X / Concentrato di soluzione di lavaggio 20X / Solución de lavado concentrada 20X / Concentrado de solução de lavagem 20X / Tvättlösningkoncentrat 20X / Συμπυκνωμένο διάλυμα πλύσης 20X / 浓缩清洗液 20X / Plovimo tirpalu koncentratas 20X / 20X mosóoldat-koncentrátum / Koncentrat 20X roztworu płuczacego / Koncentrát mycího roztoku 20X / Koncentrát premývacieho roztoku 20X / 농축 세척액(20배) / Yıkama Çözeltisi Konsantresi 20X / Концентрат промывочного раствора 20X / Концентрат на разтвор за промиване 20X / 清洗溶液濃縮 20X

## REFERENCES

1. Approved Guideline - Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures, EP17-A2. June 2012. Clinical and Laboratory Standards Institute.
2. Approved Guideline – Evaluation of Precision of Quantitative Measurement Procedures, EP05-A3. October 2014. Clinical and Laboratory Standards Institute.
3. Approved Guideline - Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach, EP06-A. April 2003. Clinical and Laboratory Standards Institute.
4. Approved Guideline - Interference Testing in Clinical Chemistry, EP07, 3rd.ed. April 2018. Clinical and Laboratory Standards Institute.
5. Approved Guideline - Supplemental Tables for Interference Testing in Clinical Chemistry, EP37, 1st ed. April 2018. Clinical and Laboratory Standards Institute.