AFF IRMA KIT

REF IM1441, IM3285

TABLE OF CONTENTS

English	2	Čeština (CZ)	29
Français (FR)	5	Slovenčina (SK)	32
Deutsch (DE)	8	한국어 (KO)	35
Italiano (IT)	11	Türkçe (TR)	38
Español (ES)	14	Русский (RU)	41
Ελληνικά (EL)	17	中文 (ZH-TW)	44
Lietuviškai (LT)	20	Србија (SR)	46
Magyar (HU)	23	APPENDIX	49
Polski (PL)	26		

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AFP IRMA KIT

REF IM1441, IM3285

IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION OF AFP

IN HUMAN SERUM OR AMNIOTIC FLUID
For *in vitro* diagnostic use

PRINCIPLE

The AFP assay is a two-step "sandwich" type assay in which two monoclonal mouse antibodies, directed against two different epitopes of the molecule, are employed. The samples or calibrators are first incubated in tubes coated with the first monoclonal antibody. The contents of the tubes are then aspirated and the presence of AFP in the sample is revealed by incubation with the second, ¹²⁵I labelled antibody. The contents of the tubes are aspirated and the bound radioactivity is determined in a gamma counter. The AFP concentrations in the samples are obtained by interpolation from the standard curve. The concentration of AFP in the samples is directly proportional to the radioactivity.

Data may be used:

- In oncology,
- In conjunction with hCG and unconjugated estriol concentrations in Down's syndrome (trisomy 21) risk assessment. Such an approach has been verified on the basis of multivariate Gaussian distribution analysis (Wald et al.: Brit Med J 297, 883 – 887; 1988). It's strongly recommended to use validated (CE marked) software designated specifically for trisomy 21 risk evaluation, e.g. software alpha™ (alpha is the trademark of Logical Medical Systems).

WARNING AND PRECAUTIONS

General remarks:

- Do not mix the reagents from kits of different lots.
- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- A standard curve must be included with each assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material are subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Materials of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HBsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All serum samples should be handled as if capable of transmitting hepatitis or AIDS and waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Wash Solution (20X)

DANGER

H360

May damage fertility or the unborn child.

P201

Obtain special instructions before use.

P280

Wear protective gloves, protective clothing and eye/face protection.

P308+P313

IF exposed or concerned:
Get medical advice/attention.
Boric Acid 0.1 - 0.3%
Sodium Borate Decahydrate
0.1 - 0.3%

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

Serum samples

- Collect blood in tubes without additive.
- Separate serum by centrifugation.
- Serum samples may be kept at 2-8°C, if they are to be assayed within 24 hours. For longer storage keep frozen at <-18°C (2 months maximum) or preferably at <-80°C (1 year maximum) in aliquots to avoid repeated freezing and thawing. Thawing of sample must be performed at room temperature.
- If samples have concentrations greater than the highest calibrator, they must be diluted into phosphate buffer.

Amniotic fluid samples

- Collect amniotic fluid in dry tubes without additive.
- Amniotic fluid samples may be kept at 2-8°C, if they are to be assayed within 24 hours. For longer storage keep frozen at <-18°C (2 months maximum) or preferably at <-80°C (1 year maximum) in aliquots to avoid repeated freezing and thawing. Thawing of sample must be performed at room temperature.
- Amniotic fluid samples are diluted 1:100 in phosphate buffer prior to assay.

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Expiry dates printed on component vial labels apply to the long-term storage by manufacturer only, prior to assembling of the kit. Do not take into account.

Storage conditions for reagents after reconstitution or dilution are indicated in paragraph Procedure.

Kit for determination of AFP, 100 tubes (Cat. # IM1441)

Anti-AFP monoclonal antibody-coated tubes: 2x 50 tubes (ready-to-use)

Monoclonal ¹²⁵I-labelled anti-AFP tracer antibody: one 22 mL (ready-to-use).

The vial contains 320 kBq at the date of manufacture of ¹²⁵I-labelled immunoglobulin in buffer containing bovine serum albumin, sodium azide (<0.1%) and a dye.

AFP calibrators: six 0.5 mL vials (ready-to-use).

The calibrator vials contain from 0 to approximately 400 IU/mL AFP in bovine serum with sodium azide (<0.1%). The exact concentrations are indicated

on each vial label. The calibrators were calibrated using the International Standard of Alpha-fetoprotein, Human 72/225. (1 IU WHO = 1.21 ng AFP)

Control samples: two vials (lyophilized)

The vials contain AFP in lyophilised human serum. The expected values are in the concentration range indicated on a supplement.

Phosphate buffer: one 30 mL vial; ready for use.

The vial contains buffer with bovine serum albumin.

Wash solution (20x): one 50 mL vial

Concentrated solution has to be diluted before use.

Kit for determination of AFP, 400 tubes (Cat. # IM3285)

Anti-AFP monoclonal antibody-coated tubes: 8 x 50 tubes (ready-to-use).

Monoclonal ¹²⁵I-labelled anti-AFP tracer antibody: four 22 mL vials. (ready-to-use).

AFP Calibrators: six 0.5 mL vials, (ready-to-use).

Control samples: two vials (lyophilized)

Phosphate buffer: four 30 mL bottles; ready-to-use

Wash solution (20X): two 50 mL vials

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- precision micropipets (50 µL)
- semi-automatic pipets (150 µL, 200 µL, 2 mL)
- vortex type mixer
- horizontal or orbital shaker
- aspiration system
- gamma counter set for 125 iodine

PROCEDURE

Preparation of reagents

Let all the reagents come to room temperature.

Reconstitution of control samples

The content of the vials is reconstituted with the volume of distilled water indicated on the label. Wait for 10 min following reconstitution and mix gently to avoid foaming before dispensing. Store the reconstituted solutions aliquoted at < -18°C, until the expiry date of the kit.

Preparation of wash solution

Pour the content of the vial into 950 mL of distilled water and homogenize. The diluted solution may be stored at 2-8°C until the expiry date of the kit.

Summary of assay procedure

Let all the reagents come to room temperature.

Step 1 1 st incubation	Step 2 Washing
To coated tubes add successively: 50 µL of calibrator or sample and 150 µL of phosphate buffer. Mix. Incubate for 15 min at 18-25°C with shaking at >280 rpm.	Aspirate carefully the contents of each tube. Add 2 mL of wash solution and aspirate carefully.

Step 3 2 nd incubation	Step 4 Washing and counting
Add 200 µL of tracer to all tubes. Mix. Incubate for 30 min at 18-25°C with shaking at >280 rpm.	Aspirate carefully the contents of each tube (except of 2 tubes for total activity T). Wash the tubes with 2 mL of wash solution and then aspirate. Count bound cpm (B) and total cpm (T) for 1 min.

*Add 200 µL of tracer to each of 2 additional tubes to obtain total cpm.

RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of AFP concentrations in samples measured at the same time as the calibrators.

Standard curve

The results in this package insert were calculated using a log-log curve fit ("spline" mode) with determined radioactivity (cpm_{cat} - cpm_{cal0}) on vertical axis and the AFP concentration of the calibrators on the horizontal axis (IU/mL). Other data reduction methods may give slightly different results.

Total activity: 107,924 cpm				
Calibrators	IU/mL	cpm (n=3)	B/T (%)	cpm _{cat} - cpm _{cal0}
0	0	20	0.02	-
1	3	2,096	1.83	2,076
2	10	6,511	5.68	6,491
3	40	22,881	20.0	22,861
4	150	61,576	53.7	61,556
5	400	94,180	82.2	94,160

(Example of standard curve, do not use for calculation)

Samples

Locate the (cpm_{sample} - cpm_{cal0}) value for each serum or diluted amniotic fluid sample on the vertical axis of the standard curve, read-off the AFP concentration of the sample on the horizontal axis in IU/mL. The concentrations of diluted samples must be corrected by the dilution factor.

To convert IU/mL to ng/mL, multiply results by 1.21.

EXPECTED VALUES

Oncology - normal AFP serum concentrations

It is suggested that each laboratory establishes its own normal values. Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

The following values obtained with healthy subjects are indicative only. The serum AFP concentrations shown in the table were obtained from 58 healthy male and 36 healthy female donors.

	Number of samples	Conc. range (IU/mL)	Median (IU/mL)
Male	58	0.46 - 6.41	3.12
Female	36	0.49 - 9.84	1.94

Screening of congenital malformations

Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Laboratoires should establish their own references values according to the method they have chosen for the Down syndrome risk assessment. Out of this context, an isolated AFP data have no diagnostic value by itself.

QUALITY CONTROL

Good laboratory practices imply that control samples must be used regularly to ensure the quality of the results obtained. These samples must be processed exactly in the same way as the assay samples, and it is recommended to analyse their results using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PERFORMANCE CHARACTERISTICS

(for more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Analytical sensitivity: 0.11 IU/mL

Functional sensitivity: 0.2 IU/mL

Specificity

The antibodies used in this kit exhibit no cross-reaction with human serum albumin, transferrin, human IgG, haemoglobin or bilirubin.

Precision

Intra-assay

Samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 6.9% for serum and 3.0% for amniotic fluid.

Inter-assay

Samples were assayed in duplicate in 10 different series. The coefficients of variation were found below or equal to 6.3% for serum and 10.9% for amniotic fluid.

Accuracy**Dilution test**

High-concentration samples were serially diluted in phosphate buffer. The recovery percentages obtained were between 86% and 111% for serum and between 92% to 112% for amniotic fluid.

Recovery test

Low-concentration samples were spiked with known quantities of AFP. The recovery percentages obtained were between 88% and 108% for serum and between 100% to 119% for amniotic fluid.

Measurement range (from analytical sensitivity to highest calibrator): 0.11 to approximately 400 IU/mL.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly.

Do not use hemolyzed, lipemic or icteric samples.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

Hook effect

There is no hook effect when the two-step procedure is used.

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