# ACTIVE® Free Testosterone RIA

**REF** DSL4900

## TABLE OF CONTENTS

English . . . . .	2	Slovenčina (SK) . . . . .	34
Français (FR) . . . . .	5	한국어 (KO) . . . . .	37
Deutsch (DE) . . . . .	8	Türkçe (TR) . . . . .	39
Italiano (IT) . . . . .	11	Русский (RU) . . . . .	42
Español (ES) . . . . .	14	中文 (ZH-TW) . . . . .	45
Português Portugal (PT-PT) . . . . .	17	Србија (SR) . . . . .	47
Svenska (SV) . . . . .	20	Português Brasil (PT-BR) . . . . .	50
Ελληνικά (EL) . . . . .	23	APPENDIX . . . . .	53
中文 (ZH-CN) . . . . .	26	REFERENCES . . . . .	56
Polski (PL) . . . . .	28		
Čeština (CZ) . . . . .	31		

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# ACTIVE® Free Testosterone RIA

**REF** DSL4900

## RADIOIMMUNOASSAY FOR THE QUANTITATIVE MEASUREMENT OF FREE TESTOSTERONE IN HUMAN SERUM

**THIS ASSAY IS INTENDED FOR IN VITRO DIAGNOSTIC USE**

### PRINCIPLE

The radioimmunoassay of free testosterone is a competitive assay. The procedure follows the basic principle of radioimmunoassay where there is competition between a radioactive and a non-radioactive antigen for a fixed number of antibody binding sites [1]. The amount of [<sup>125</sup>I]-labeled testosterone analog bound to the antibody is inversely proportional to the concentration of unlabeled free testosterone present. The separation of the free and bound antigen is achieved by decanting or aspirating the antibody-coated tubes. A standard curve is constructed and unknown free testosterone values are obtained from the curve by interpolation.

*For Summary and Explanation of the Test see APPENDIX.*

### WARNING AND PRECAUTIONS

#### General remarks:

- For *in vitro* diagnostic use.
- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- Do not use any component beyond the expiration date shown on its label.
- A standard curve must be established with each assay.
- It is recommended to perform the assay in duplicate.
- Calibrators and controls should be mixed before use by inverting or swirling gently rather than vortexing.

#### Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material are subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should take place in the presence of radioactive materials.
- No pipeting by mouth.
- Avoid all contact with radioactive materials by using gloves and lab coat.
- All manipulation of radioactive substances should be done in an appropriate location, away from corridors and other busy areas.
- Radioactive materials should be stored in the container provided and in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

#### Sodium azide

Some reagents contain sodium azide as preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal of liquids, flush with a large volume of water to prevent azide build-up [2].

#### Materials of human origin

Patient samples and blood-derived products may be routinely processed with minimum risk using the procedure described. However, handle these products as potentially infectious according to universal precautions and good clinical laboratory practices, regardless of their origin, treatment, or

prior certification. Use an appropriate disinfectant for decontamination. Store and dispose of these materials and their containers in accordance with local regulations and guidelines.

### GHS HAZARD CLASSIFICATION

Not classified as hazardous

**SDS**

Safety Data Sheet is available at [techdocs.beckmancoulter.com](http://techdocs.beckmancoulter.com)

### SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Serum is the recommended sample type.
- Allow serum samples to clot completely before centrifugation.
- Serum samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage keep frozen at < -20°C up to 2 months [3]. It is recommended to prepare aliquots to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- Frozen samples should be thawed and mixed thoroughly by gentle swirling or inversion prior to use.
- Any sample reading above the concentration of the highest calibrator should not be diluted. Dilution will change the proportion of free to bound testosterone.

### MATERIALS PROVIDED

All unopened reagents in the kit are stable until the expiration date indicated on the kit label, when stored at 2-8°C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Storage conditions for opened reagents are indicated in appropriate paragraphs.

#### Free Testosterone Antibody-Coated Tubes: 2 x 50 tubes (ready-to-use)

Plastic tubes with rabbit anti-testosterone immunoglobulin immobilized to the inside wall of each tube.

#### <sup>125</sup>I-labeled Testosterone Analog Tracer (YELLOW): one 22 mL vial (ready-to-use)

At the time of manufacture, the vial contains 185 kBq, (<5 µCi), of <sup>125</sup>I-labeled testosterone analog in buffer with proteins (BSA), sodium azide (<0.1%) and a dye.

#### Calibrators: one 1.0 mL vial labeled 0, and seven 0.5 mL vials labeled 1-7 (ready-to-use)

The calibrator vials contain from 0 to approximately 100 pg/mL (0 to approximately 346.7 pmol/L) of free testosterone in human serum with sodium azide (< 0.1%). The exact concentration is indicated on each vial label. Once opened, store at 2-8°C for up to 3 weeks, or at < -20°C until expiration date of kit. The calibrator values were established using an internal standard.

#### Controls: two 0.5 mL vials labeled 1, 2 (ready-to-use)

The vials contain free testosterone in human serum with sodium azide (<0.1%). The expected values are in the concentration range indicated in a supplement found in the kit. Once opened, store at 2-8°C for up to 3 weeks, or at < -20°C until expiration date of kit.

### MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- 12 x 75 mm plastic or glass test tubes for total counts.
- Test tube rack for 12 x 75 mm tubes.
- Precision micropipet (50 µL).
- Semi-automatic pipet (200 µL, 3.0 mL).
- Deionized water.
- Vortex mixer.
- Temperature-controlled water bath, 37°C ± 2°C.
- A sponge rack or similar device for decantation.
- Absorbent material for blotting tubes.
- Gamma counter set for 125 iodine.

## PROCEDURE

### Preparation of reagents

Let all the reagents come to room temperature and mix them thoroughly by gentle inversion before use.

### Assay procedure

Run calibrators, controls and patient samples in duplicate.

Step 1 Additions	Step 2 Incubation & Washing	Step 3 Washing & Counting
To the bottom of antibody coated tubes successively add:	Cover and incubate 60 min in a water bath at $37 \pm 2^\circ\text{C}$ .	Decant all tubes, (except «total cpm» tubes), by simultaneous inversion with a sponge rack into a radioactive waste receptacle.
50 $\mu\text{L}$ of calibrator, control or sample, immediately add 200 $\mu\text{L}$ of tracer.*	Decant all tubes, (except «total cpm» tubes), by simultaneous inversion with a sponge rack into a radioactive waste receptacle.	Strike the tubes sharply on absorbent material to facilitate complete drainage.
Mix by vortex.	Strike the tubes sharply on absorbent material to facilitate complete drainage.	Drain on absorbent material for a minimum of 2 minutes. Blot droplets from the rims of tubes.
	Add 3.0 mL of deionized water to all tubes, except «total cpm» tubes.	Count bound cpm (B) and total cpm (T) for 1 minute.

\*Add 200  $\mu\text{L}$  of tracer to 2 additional tubes to obtain «total cpm».

## RESULTS

Results are obtained from the standard curve by interpolation. The curve is used for the determination of free testosterone concentrations in samples measured at the same time as the calibrators.

### Standard curve

The results in the quality control department were calculated using *spline* curve fit with  $B/T$  or  $B/B_0$  on the lin vertical axis and analyte concentration of the calibrators on the log horizontal axis (pg/mL).

Other data reduction methods may give slightly different results.

$ED_{50} = 3.710 \text{ pg/mL}$

Total activity: 62,502 cpm				
Calibrators	Free testosterone (pg/mL)	cpm (n=2)	B/T (%)	B/B <sub>0</sub> (%)
0	0	28,056	44.89	100.0
1	0.23	23,889	38.22	85.15
2	0.90	18,675	29.88	66.56
3	2.80	15,016	24.03	53.52
4	9.50	10,466	16.75	37.30
5	26.00	6,682	10.69	23.82
6	50.00	4,466	7.15	15.92
7	100.0	2,989	4.78	10.66

(Example of standard curve, do not use for calculation)

### Samples

For each sample locate the ratio  $B/T$  or  $B/B_0$  on the vertical axis of the standard curve and read off the corresponding free testosterone concentration of the sample on the horizontal axis in pg/mL.

To convert concentrations from pg/mL to pmol/L, multiply results by 3.47.

## EXPECTED VALUES

Each laboratory should establish its own reference ranges. The following values were obtained using serum samples from apparently normal healthy adults.

Men (Age in years)	n	Median (pg/mL)	Min (pg/mL)	Max (pg/mL)	2.5 <sup>th</sup> Percentile (pg/mL)	97.5 <sup>th</sup> Percentile (pg/mL)
20 - 29	30	15.40	7.54	27.79	8.68	25.09
30 - 39	30	14.94	6.77	22.36	8.85	21.40
40 - 49	30	11.48	7.37	20.47	7.56	18.64
50 and above	30	9.05	4.52	15.05	5.72	14.21
All	120	12.35	4.52	27.79	6.76	22.76

Women	n	Median (pg/mL)	Min (pg/mL)	Max (pg/mL)	2.5 <sup>th</sup> Percentile (pg/mL)	97.5 <sup>th</sup> Percentile (pg/mL)
Follicular Phase	32	1.48	0.59	4.28	0.64	3.41
Luteal Phase	32	1.44	0.59	3.31	0.60	2.95
Preovulatory Peak	21	1.51	0.81	3.79	0.90	3.79
Postmenopausal	13	1.17	0.33	1.97	0.36	1.85
Contraceptives	19	1.45	0.49	2.81	0.57	2.70

Random women	n	Median (pg/mL)	Min (pg/mL)	Max (pg/mL)	2.5 <sup>th</sup> Percentile (pg/mL)	97.5 <sup>th</sup> Percentile (pg/mL)
All	120	1.29	0.33	4.28	0.49	2.87

Detail information about expected values for children (sorted according to age and sex) can be found in the data sheet "APPENDIX".

## QUALITY CONTROL

Good laboratory practices require that control samples be used regularly to ensure the quality of the results obtained. These controls must be processed in exactly the same way as the patient samples, and it is recommended that their results be analyzed using appropriate statistical methods.

Failure to obtain the appropriate values for controls may indicate imprecise manipulations, improper sample handling or deterioration of reagents. In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or [imunochem@beckman.com](mailto:imunochem@beckman.com).

In the US, contact the Beckman Coulter technical support at 1-800-854-3633; or by email at: [immunoassay@beckman.com](mailto:immunoassay@beckman.com)

## PERFORMANCE CHARACTERISTICS

(For more details, see APPENDIX)

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

### Sensitivity

Limit of detection (LoD): 0.13 pg/mL

The LoD for free testosterone is 0.13 pg/mL, determined consistent with guidelines in CLSI document EP17-A2 [4] based on the proportions of false positives ( $\alpha$ ) less than 5% and false negatives ( $\beta$ ) less than 5%; using determinations, with 120 blank and 120 low level samples; and Limit of Blank (LoB) of 0.05 pg/mL.

Limit of quantitation (LoQ): 0.37 pg/mL

The LoQ for Free Testosterone is 0.37 pg/mL, determined consistent with guidelines in CLSI document EP17-A2 [4] as the lowest concentration where criterion for maximum total error was achieved.

### Specificity

The antibody used in the immunoassay is highly specific for free testosterone. Low cross reactivities were obtained against several related molecules.

### Precision

#### Repeatability and within-laboratory precision

The precision of the assay was determined consistent with guidelines in CLSI document EP05-A3 [5]. For repeatability the coefficients of variation were found below or equal to 13.2% for serum samples. For within laboratory precision the coefficients of variation were found below or equal to 19.3% for serum samples.

### Accuracy

#### Linearity

The assay demonstrated to be linear from 0.30 to 114.11 pg/mL using serum samples (determined consistent with guidelines in CLSI document EP6-A [6]).

Measurement range (from LoQ to the highest calibrator):

0.37 pg/mL to approximately 100 pg/mL

## LIMITATIONS

- Failure to follow these instructions for use (IFU) may significantly affect results.
  - Failure to blot tubes adequately following decantation may result in poor replication and spurious results.
  - Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.
  - Avoid repeated freezing and thawing of reagents or specimens.
  - For more details, see the data sheet "APPENDIX".
- The possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays. Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.
-

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## APPENDIX

### PERFORMANCE CHARACTERISTICS

**ACTIVE** is a trademark of BECKMAN COULTER Inc. and its subsidiaries.

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary

#### SUMMARY AND EXPLANATION OF THE TEST

**Indications for use:** Radioimmunoassay for the quantitative measurement of free testosterone in human serum. This assay is intended for in vitro diagnostic use. Free testosterone test is used in the diagnosis and treatment of disorders involving the male sex hormones (androgens), including primary and secondary hypogonadism, delayed or precocious puberty, impotence in males and, in female's hirsutism (excessive hair) and virilization (masculinization) due to tumors, polycystic ovaries, and adrenogenital syndromes [7].

Testosterone (17 $\alpha$ -Hydroxy-4-androsten-3-one), a C<sub>19</sub> steroid, is a potent circulating naturally-secreted androgen [8]. In normal postpubertal males, testosterone is secreted primarily by the testes with only a small amount derived from peripheral conversion of 4-Androstene-3,17-dione (ASD) [8]. In adult women, it has been estimated that more than 50% of serum testosterone is derived from peripheral conversion of ASD secreted by the adrenal glands and ovaries, with the remainder derived from direct secretion of testosterone by these glands [9,10]. Testosterone production rates in blood have been estimated at 0.34 mg/day in adult females and more than 20 times this amount in adult males [9,10,11]. Approximately 60% of blood testosterone is normally bound with high affinity to sex hormone-binding globulin (SHBG) [12]; of the remainder, all but 1-2% is loosely bound to albumin.

The percentage bound to SHBG is somewhat lower in males as compared to females [13]. Both the albumin-bound and free fractions may be biologically active, while SHBG effectively inhibits testos-terone action. Approximately half of the blood testosterone is metabolized in the liver to androsterone, etiocholanolone and epiandrosterone; all relatively weak androgens [8,10]. Testosterone is also converted to the potent androgen, dihydrotestosterone (DHT), in certain target tissues. Measurement of the free or unbound fraction of serum testosterone has been proposed as a means of estimating the physiologically bioactive hormone. Free testosterone levels are elevated in women with hyperandro-genism associated with hirsutism in the presence or absence of polycystic ovarian disease [12,13,14].

Traditional methods for measurement of free testosterone involve the addition of [<sup>3</sup>H]-testosterone to serum or plasma, followed by various procedures, e.g. dialysis, filtration and precipitation, to separate the free and bound fractions. The partitioning of the tritiated testosterone is then used to calculate the endogenous bound and free fractions. These indirect methods are time-consuming and labor-intensive and are not practical for routine applications. The DSL4900 ACTIVE Free Testosterone RIA utilizes an [<sup>125</sup>I]-labeled testosterone analog which has low affinity for SHBG and albumin. The analog com-petes with the unbound testosterone in the test sample for binding to specific anti-testosterone polyclonal antibodies, which have been immobilized on the assay tube. This competitive binding format allows direct estimation of unlabeled free testosterone levels in unextracted samples.

#### Interference

Serum samples containing Free Testosterone concentrations of app. 5 and 50 pg/mL were spiked with multiple concentrations of the substances below and assayed using Free Testosterone RIA kit. Values were calculated as described in CLSI EP07 3<sup>rd</sup> Edition [15]. Interference was determined by testing controls (no interfering substance added) and matched test samples

(with interfering substance added). No interference (defined as a shift in dose >10%) was found for interferent up to concentration stated in the table below.

Interferent	Test concentration
Hemoglobin	641 $\mu$ g/mL
Bilirubin (conjugated)	199 $\mu$ g/mL
Bilirubin (unconjugated)	92.7 $\mu$ g/mL
Biotin	1,020 ng/mL
Ascorbic acid	40.5 $\mu$ g/mL
Ibuprofen	122 $\mu$ g/mL
Cholesterol	0.75 mg/mL
Heparin	7,248 ng/mL
Prednisone	384 ng/mL
Prednisolone	1,294 ng/mL
Protein (gamma globulin)	49.5 mg/mL
Rheumatoid factor	3.45 IU/mL
Acetylsalicylic acid	24.2 $\mu$ g/mL
Triglycerides	3.35 mg/mL

#### Specificity

The cross-reactivity of the free testosterone antiserum has been measured against the following compounds according to CLSI recommendations (EP07-3<sup>rd</sup> Edition) [15]. The percent cross-reactivity is expressed as the ratio of measured value minus true value and concentration of cross-reactant.

COMPOUND	Pooled normal serum	
	Crossreactant Conc. (ng/mL)	Cross Reactivity (%)
19-Nortestosterone	8.62	0.068
11-Ketotestosterone	102	0.020
11 $\beta$ -Hydroxytestosterone	14.1	0.029
Androstenediol	14.8	0.010
Methyltestosterone	96.5	0.017
Androstenedione	92.4	0.006
Norethindrone	1,056	0.001
5 $\alpha$ -Androstane-3 $\alpha$ ,17 $\beta$ -diol	1,076	0.002
5 $\beta$ -Androstane-3,17-dione	1,011	ND
5 $\alpha$ -Androstane-3,17-dione	846	0.001
Danazol	921	0.001
Estradiol	946	0.001
Progesterone	956	0.001
Deoxycorticosterone	911	ND
Androsterone	1,002	ND
Epiandrosterone	943	ND
DHEA	896	ND
Estriol	882	ND
Estrone	913	ND
Corticosterone	951	ND
Cortisol	857	ND
Pregnenolone	944	ND
DHT	168	0.020
Androstenediol Glucuronide	1,031	ND
DHEA-S	922	ND
5 $\alpha$ -Androstan-3 $\beta$ ,17 $\beta$ diol	98.0	0.005
rhSHBG	10,407	ND

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Five serum samples were assayed for 20 days, 2 runs per day, triplicates per run. Assays were performed by two lab technicians, by two reagent lots. There were 120 individual measurements per sample with no invalid results.

Serum	Mean (pg/mL)	Repeatability		Within-laboratory precision	
		SD (pg/mL)	C.V. (%)	SD (pg/mL)	C.V. (%)
S1	0.65	0.09	13.23	0.12	19.32
S2	7.89	0.71	8.95	0.90	11.43
S3	13.65	1.44	10.54	1.67	12.20
S4	40.29	2.58	6.42	3.49	8.68
S5	74.29	4.80	6.47	5.78	7.78

#### Expected data for children

Results are sorted according to the age and sex.

Subgroups "Infant, Child and Adolescent" were defined according to Premarket Assessment of Pediatric Medical Devices, Guidance for Industry and Food and Drug Administration Staff, March 2014.

The second tables provide more detailed age discrimination at the onset of puberty.

Free Testosterone (pg/mL)						
Boys	N	Median	Min.	Max.	2.5 <sup>th</sup> percentile	97.5 <sup>th</sup> percentile
Infant	18	< 0.13	< 0.13	0.31	< 0.13	0.28
Child	48	0.56	0.13	12.79	0.16	12.65
Adolescent	30	15.10	1.71	32.09	1.88	28.55

Free Testosterone (pg/mL)						
Boys	N	Median	Min.	Max.	2.5 <sup>th</sup> percentile	97.5 <sup>th</sup> percentile
6 months - 9 years	36	0.20	< 0.13	0.62	< 0.13	0.54
10 - 11 years	20	0.67	0.41	5.11	0.42	5.00
12 - 13 years	20	6.21	0.60	27.21	0.63	23.27
14 - 15 years	20	18.71	5.50	32.09	8.03	28.77

Free Testosterone (pg/mL)						
Girls	N	Median	Min.	Max.	2.5 <sup>th</sup> percentile	97.5 <sup>th</sup> percentile
Infant	34	< 0.13	< 0.13	0.42	< 0.13	0.33
Child	55	0.46	< 0.13	2.53	< 0.13	1.90
Adolescent	30	1.37	0.57	3.27	0.64	3.23

Free Testosterone (pg/mL)						
Girls	N	Median	Min.	Max.	2.5 <sup>th</sup> percentile	97.5 <sup>th</sup> percentile
6 months - 9 years	69	0.24	< 0.13	0.86	< 0.13	0.57
10 - 12 years	20	0.88	0.37	2.53	0.41	2.25
13 - 16 years	20	1.42	0.57	3.27	0.65	3.24

#### Method Comparison:

The current version of Active® Free Testosterone RIA (Method A) was compared to previous version of Active® Free Testosterone RIA (Method B):

n	Range (pg/mL)	Slope	Intercept	r
278	0.26 - 101.21	1.0078	0.0759	0.9956

#### <sup>125</sup>I Characteristics

T<sub>1/2</sub> (<sup>125</sup>I) = 1443 h = 60.14 d

<sup>125</sup> I	E (MeV)	%
Y	0.035	
X	0.027	114
	0.032	25

## Symbols Key

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biyolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktivt / Ραδιενεργό / 放射性 / Radioaktivtyoij medžiaga / Radioaktiv / Radioaktywny / Radioaktivní / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / / Ανιχνευτής / 追踪剂 / Atsekamoji medžiaga / Nyomjelző / Znacznik / Radioindikátor / Indikátor (tracer) / 트레이서 / Tracer'lar / метка / Индикатор / 追蹤劑

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정 물질 / Kalibratör / Калибратор / Калибратор / 校正液

Control / Contrôle / Kontrolle / Contollo / Control / Controllo / Kontrolle / Μάρτυρας / 质控品 / Kontrolinė / Kontroll / Kontrola / Kontrola / Контроль / Контроль / 質控品

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνώματα / 试管 / Mégintüveliai / Csövek / Probówki / Zkumavky / Skúmavky / 튜브 / Tüpler / пробирки / Епруветки / 試管

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用说明

Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferens / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Referenca k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Ürün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnostikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / Ин vitro диагностика / 體外診斷

Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 组成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 制造商 / Gamintojas / Gyártó / Producent / Výrobce / Výrobca / 제조 / Üretici / Изготовлено / Произведено от / 製造商

Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Mercado CE / Marcação CE / CE-márkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Знаčka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Zapoznać się z instrukcją użycia / Postępujcie podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문의 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижете Инструкциите за употреба / 請參閱使用說明

Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 温度范围 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sıcaklık aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Ispejimas / Figyelem / Uwaga / Urozozmení / Urozozomenie / 주의 / Dikkat / Внимание / 注意

Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galijimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirácie / 만료 날짜 / Son Kullanma Tarihi / Срок годности / Срок на годност / 到期日

Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer serii / Číslo šarže / 로트 번호 / Lot Numarası / Номер партии / Номер на партида / 批號

Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Dátum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

21 February 2020

## REFERENCES

1. Yalow R, Berson S: Introduction and general considerations. IN: Odell WD, Daughday WH (eds): J.B. Lippencott Co., Philadelphia, 1971, pp. 1-19.
2. DHHS (NIOSH) Publication No. 78-127, August 1976. Current Intelligence Bulletin 13 - Explosive Azide Hazard. Available <http://www.cdc.gov/niosh>.
3. Tietz Clinical Guide To Laboratory Tests (4th. edition). Alan H.B.Wu. Saunders Elsevier, 2006.
4. Approved Guideline - Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures, EP17-A2. June 2012. Clinical and Laboratory Standards Institute.
5. Approved Guideline – Evaluation of Precision of Quantitative Measurement Procedures, EP05-A3. October 2014. Clinical and Laboratory Standards Institute.
6. Approved Guideline - Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach, EP6-A. April 2003. Clinical and Laboratory Standards Institute.
7. Code of Federal Regulations; Title 21, Volume 8, Section 862.1680, Revised as of April 1, 2019.
8. Dorfman RI, Shipley RA: Androgens. John Wiley and Sons, Inc., New York, 1956, pp. 116-128.
9. Horton R, Tait JF: Androstenedione production and interconversion rates measured in peripheral blood and studies on the possible site of its conversion to testosterone. J Clin Invest 45:301-313, 1966.
10. Pang S, Riddick L: Hirsutism. IN Lifshitz F (ed): Marcel Dekker, Inc., New York, 1990, pp. 259-291.
11. Androgen resistance syndromes: 5 $\alpha$ -reductase deficiency, testicular feminization, and related syndromes. IN Scriver CR, Beaudet AL, Sly WS, Valle D (eds): The Metabolic Basis of Inherited Disease, 6th edition. McGraw-Hill, New York, 1989, pp. 1919-1944.
12. Cumming DC, Wall SR: Non-sex hormone binding globulin-bound testosterone as a marker for hyperandrogenism. J Clin Endocrinol Metab 61:873-876, 1985.
13. Westphal U: Steroid-Protein Interactions. Springer-Verlag, Berlin, 1986
14. Pugeat M, Nicolas MH, Craves JC, Alvarado-Dubost C, Fimbel S, Dechaud, Lejeune H: Androgens in polycystic ovarian syndrome. IN Tolis G, Bringer J, Chrousos GP (eds): Intraovarian regulators and polycystic ovarian syndrome. NY Academy of Sciences, New York, 1993, pp. 124-135.
15. Approved Guideline - Interference Testing in Clinical Chemistry, EP07 3<sup>rd</sup> Edition. April 2018. Clinical and Laboratory Standards Institute.