IRMA α subunit

REF IM1186

TABLE OF CONTENTS

English	2	Polski	22
Français	4	Čeština	25
Deutsch	7	Slovenčina	27
Italiano	10	한국어	29
Español	12	Türkçe	31
Português (Portugal)	15	Русский	33
Svenska	17	APPENDIX	36
Ελληνικά	19		

IRMA α subunit

REF IM1186

IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION OF α SUBUNIT IN HUMAN SERUM AND PLASMA

For *in vitro* diagnostic use.

PRINCIPLE

The α subunit assay is a one-step "sandwich" type assay in which two mouse monoclonal antibodies, directed against two different epitopes of the molecule, are employed.

Samples or calibrators are incubated in tubes, coated with the first monoclonal antibody, in presence of the second ^{125}I -labeled monoclonal antibody. Following incubation, the content of the tubes is aspirated and unbound labeled antibody is eliminated by washing. The amount of bound reactivity measured in a gamma counter is proportional to the α subunit concentration. The unknown values are determined by interpolation from a standard curve.

WARNING AND PRECAUTIONS

General remarks:

- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be included with each assay.
- The correct setting of the shaker is very important for the reproducibility of the assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Materials of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HBsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All serum and plasma samples should be handled as if capable of transmitting hepatitis or AIDS. Waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Wash Solution (20X) DANGER

H360	May damage fertility or the unborn child.
P201	Obtain special instructions before use.
P280	Wear protective gloves, protective clothing and eye/face protection.
P308+P313	IF exposed or concerned: Get medical advice/attention. Boric Acid 0.1-0.3% Sodium Borate Decahydrate 0.1-0.3%

SDS

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in dry tubes or heparinized tubes.
- Separate serum or plasma from cells by centrifugation.
- Serum and plasma samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage keep frozen (at <-18°C) up to 6 months. Samples should be aliquoted to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.

If samples have concentration greater than the highest calibrator, they must be diluted into separately provided Diluent (cat. #B64533). DO NOT USE zero calibrator.

Serum and heparin plasma values for 15 samples (serum values ranging from 0.21 to 1.12 IU/L) were compared using the IRMA α subunit kit. Results are as follows:

$$[\text{heparin-plasma}] = 1.0165 [\text{serum}] + 0,0036$$

$$R = 0.9901$$

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Expiry dates printed on component vial labels apply to the long-term storage by manufacturer only, prior to assembling of the kit. Do not take into account.

Storage conditions for reagents after dilution are indicated in paragraph Procedure.

Anti- α subunit antibody-coated tubes: 2 x 50 tubes (ready-to-use)

^{125}I -labeled monoclonal antibody tracer: one 11 mL vial (ready-to-use)

The vial contains 370 kBq at time of manufacture of ^{125}I -labeled monoclonal antibody with bovine serum albumine, sodium azide (<0.1%) and a dye.

Calibrators: six 1 mL vials + one 2 mL zero calibrator vial (ready-to-use)

The calibrator vials contain between 0 and approximately 10 IU/L of α subunit in buffer containing sodium azide (<0.1%). The exact concentration is indicated on each vial label. The calibrators were calibrated against the international standard WHO 75/569, (1 IU = 1 μg).

Control serum: two 1 mL vials (ready-to-use)

The vials contain α subunit in buffer containing sodium azide (<0.1%). The expected values are in the concentration range indicated on a supplement.

Wash solution (20 X): one 50 mL vial

Concentrated solution has to be diluted before use.

MATERIALS REQUIRED, BUT NOT PROVIDED

Sample diluent: Diluent (US Estradiol + α -subunit), 2 mL (ready to use), ordered separately (cat. n. B64533)

In addition to standard laboratory equipment, the following items are required:

- precision micropipet (100 μL).

- repeating micropipets (100 µL, 2 mL).
- horizontal or orbital shaker.
- aspiration system.
- gamma counter.

RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of α subunit concentrations in samples measured at the same time as the calibrator.

Standard curve

The results in the package insert were calculated using a log-log curve fit ("spline" mode) with B/T (%) or B/B_{max} (%) on vertical axis and the α subunit concentration of the calibrators on the horizontal axis (IU/L). Other data reduction methods may give slightly different results.

Total activity: 136,828 cpm				
Calibrators	α subunit (IU/L)	cpm	B/T (%)	B/Bmax (%)
0	0	83	0.06	-
1	0.38	823	0.60	3.55
2	0.70	1,407	1.03	6.06
3	1.30	2,538	1.85	10.94
4	2.90	6,250	4.57	26.94
5	5.30	12,170	8.89	52.45
6	10.0	23,202	16.96	100

(Example of standard curve, do not use for calculation)

Samples

Locate the B/T (%) or the B/B_{max} (%) for each sample on the vertical axis and read off corresponding α subunit concentration in IU/L on the horizontal axis.

EXPECTED VALUES

We recommend each laboratory to establish its own reference values. The following values obtained from healthy subjects are indicative only.

Men	Premenopausal women	Postmenopausal women
0 - 0.7 IU/L	0 - 0.6 IU/L	0 - 1.3 IU/L

QUALITY CONTROL

Good laboratory practices imply that control samples must be used regularly to ensure the quality of the results obtained. These samples must be processed exactly in the same way as the assay samples, and it is recommended to analyse their results using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PROCEDURE

Preparation of reagents

Let all the reagents come to room temperature.

Preparation of the wash solution

Pour the content of the vial into 950 mL of distilled water and homogenize. The diluted solution can be stored at 2-8°C until the expiry date of the kit.

Procedure

Immunological step	Washing step	Counting
To antibody coated tubes add 100 µL of sample, calibrator, or control and 100 µL of tracer.* Incubate 1 hour at 18-25°C with shaking (>280 rpm).	Aspirate carefully the content of each tube (except "total cpm"), then wash twice with 2 mL of wash solution, immediately aspirate the contents of tubes.	Count bound cpm (B) and total cpm (T) for 1 min.

*Add 100 µL of tracer in 2 additional tubes to obtain "total cpm".

PERFORMANCE CHARACTERISTICS

(For more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Analytical sensitivity: 0.03 IU/L

Functional sensitivity: 0.04 IU/L

Specificity: the kit is specific for the free alpha subunit.

Precision

Intra-assay

Serum samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 4.23%.

Inter-assay

Serum samples were assayed in duplicate in 10 different series. Coefficients of variation were found below or equal to 5.10%.

Accuracy

Dilution test

High-concentration serum samples were serially diluted with diluent. The recovery percentages obtained were between 82.6% and 118%.

Recovery test

Low-concentration serum samples were spiked with known quantities of α subunit. The recovery percentages obtained were between 94.0% and 99.2%.

Measurement range (from analytical sensitivity to highest calibrator):

0.03 to approximately 10 IU/L.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly. Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Do not use hemolyzed, lipemic or icteric samples.

Hook effect

No hook effect was observed until 1,300 IU/L.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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Specificity

Cross-reactivity is expressed in the following table as a percentage calculated as the concentration of cross reacting substance that gives the same binding as 10 IU/L of free α subunit of hCG.

Compound	Cross reactivity
hTSH	< 0.1%
hLH	< 0.1%
hFSH	< 0.1%
hCG	< 0.1%
β (hCG, hLH, hFSH, hTSH)	0

Precision

Intra-assay

Serum	S1	S2	S3
Number of determinations	25	25	25
Mean value, IU/L	0.81	3.56	8.06
C.V., %	4.23	2.22	1.40

Plasma-heparin	P1	P2	P3
Number of determinations	25	25	25
Mean value, IU/L	0.53	1.79	4.35
C.V., %	3.35	1.70	1.67

Inter-assays

Serum	S1	S2	S3
Number of determinations	10	10	10
Mean value, IU/L	0.88	3.67	7.80
C.V., %	5.10	4.56	2.38

Plasma-heparin	P1	P2	P3
Number of determinations	10	10	10
Mean value, IU/L	1.30	5.87	10.60
C.V., %	6.71	9.44	10.65

Accuracy

Dilution test

Three serum/heparin plasma samples were diluted in diluent and assayed according to the assay procedure of the kit.

Serum	Dilution	α subunit (IU/L)		Ratio (%) Measured/ Expected
		Measured	Expected	
S1	-	2.13	-	-
	1:2	1.12	1.07	105.2
	1:4	0.56	0.53	105.2
	1:8	0.25	0.27	93.9
	1:16	0.11	0.13	82.6
	1:32	0.07	0.07	105.2
S2	-	3.61	-	-
	1:2	1.88	1.81	104.2
	1:4	0.92	0.90	101.9
	1:8	0.50	0.45	110.8
	1:16	0.24	0.23	106.4
S3	-	4.91	-	-
	1:2	2.69	2.6	109.6
	1:4	1.31	1.23	106.7
	1:8	0.68	0.61	110.8
	1:16	0.36	0.31	117.3

Heparin plasma	Dilution	α subunit (IU/L)		Ratio (%) Measured/ Expected
		Measured	Expected	
P1	-	3.01	-	-
	1:2	1.58	1.51	105.0
	1:4	0.86	0.75	114.3
	1:8	0.39	0.38	103.7
	1:16	0.18	0.19	95.7
P2	-	4.53	-	-
	1:2	2.45	2.227	108.2
	1:4	1.22	1.13	107.7
	1:8	0.67	0.57	118.3
	1:16	0.29	0.28	102.4
P3	-	6.62	-	-
	1:2	3.42	3.31	103.3
	1:4	1.77	1.66	106.9
	1:8	0.95	0.83	114.8
	1:16	0.46	0.41	111.2
1:32	0.19	0.21	91.8	

Recovery test

To three sera/heparin plasma containing low α subunit concentrations were added known concentrations of α subunit. Samples are assayed according to the kit procedure.

Serum	Endogen. conc. (IU/L)	Added conc. (IU/L)	Expected conc. (IU/L)	Measured conc. (IU/L)	Ratio (%) Measured/ Expected
S1	1.74	0.85	2.59	2.57	99.19
	1.59	1.57	3.16	3.04	96.25
	1.47	2.17	3.64	3.42	94.00
S2	2.37	0.85	3.23	3.12	96.68
	2.18	1.57	3.74	3.55	94.88
	2.01	2.17	4.18	3.95	94.57
S3	3.14	0.85	3.99	3.94	98.72
	2.88	1.57	4.44	4.25	95.68
	2.65	2.17	4.82	4.61	95.58

Heparin plasma	Endogen. conc. (IU/L)	Added conc. (IU/L)	Expected conc. (IU/L)	Measured conc. (IU/L)	Ratio (%) Measured/ Expected
P1	0.57	0.85	1.43	1.50	105.10
	0.53	1.57	2.09	2.11	100.88
	0.48	2.17	2.65	2.48	93.45
P2	0.21	0.85	1.06	1.05	98.72
	0.19	1.57	1.76	1.54	87.58
	0.18	2.17	2.35	2.13	90.79
P3	0.28	0.85	1.14	1.14	100.32
	0.26	1.57	1.83	1.72	94.25
	0.24	2.17	2.41	2.29	95.11

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