Active® Androstanediol Glucuronide RIA

REF DSL9200

TABLE OF CONTENTS

English	2	Ελληνικά (EL)	21
Français (FR)	4	Čeština (CZ)	24
Deutsch (DE)	7	Slovenčina (SK)	26
Italiano (IT)	10	Русский (RU)	28
Español (ES)	13	Србија (SR)	31
Português Portugal (PT-PT)	16	APPENDIX	33
Svenska (SV)	19		

Beckman Coulter do Brasil Com. e Imp. de Prod. de Lab. Ltda
Alameda Rio Negro, 500, 15º andar, Torre B Alphaville Industrial
CEP 06.454-00 Barueri, São Paulo, Brasil
CNPJ: 42.160.812/0001-44 Telefone: 0800-771-8818

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Active® Androstenediol Glucuronide RIA

REF DSL9200

**RADIOIMMUNOASSAY FOR THE QUANTITATIVE
MEASUREMENT OF ANDROSTANEDIOL GLUCURONIDE
(3 α -diol G) IN HUMAN SERUM
THIS ASSAY IS INTENDED FOR IN VITRO DIAGNOSTIC
USE**

PRINCIPLE

The radioimmunoassay of 3 α -androstenediol glucuronide (3 α -diol G, 5 α -Androstane-3 α ,17 β -diol glucuronide), is a competition assay. The procedure follows the basic principle of radioimmunoassay where there is competition between a radioactive and a non-radioactive antigen for a fixed number of antibody binding sites. The amount of [125I]-labeled 3 α -diol G bound to the antibody is inversely proportional to the concentration of unlabeled 3 α -diol G present. The separation of the free and bound antigen is achieved by decanting or aspirating the antibody-coated tubes. A standard curve is constructed and unknown 3 α -diol G values are obtained from the curve by interpolation.

For Summary and Explanation of the Test see APPENDIX.

WARNING AND PRECAUTIONS

General remarks:

- For *in vitro* diagnostic use.
- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- Do not use any component beyond the expiration date shown on its label.
- A standard curve must be established with each assay.
- It is recommended to perform the assay in duplicate.
- Calibrators and controls should be mixed before use by inverting or swirling gently rather than vortexing.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection.

- No eating, drinking, smoking or application of cosmetics should take place in the presence of radioactive materials.
- No pipeting by mouth.
- Avoid all contact with radioactive materials by using gloves and lab coat.
- All manipulation of radioactive substances should be done in an appropriate location, away from corridors and other busy areas.
- Radioactive materials should be stored in the container provided and in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some Reagents contain sodium azide as preservative. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal of liquids, flush with a large volume of water to prevent azide build-up.

Materials of human origin

Patient samples and blood-derived products may be routinely processed with minimum risk using the procedure described. However, handle these

products as potentially infectious according to universal precautions and good clinical laboratory practices, regardless of their origin, treatment, or prior certification. Use an appropriate disinfectant for decontamination. Store and dispose of these materials and their containers in accordance with local regulations and guidelines.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Serum is recommended sample type.
- Allow specimens to clot completely before centrifugation.
- Serum samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage keep frozen at < -20°C, 6 months maximum. It is recommended to prepare aliquots to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.

MATERIALS PROVIDED

All unopened reagents in the kit are stable until the expiration date indicated on the kit label, when stored at 2-8°C. Expiration dates printed on vial labels apply to the long-term storage of components by the manufacturer only.

Storage conditions for opened reagents are indicated in appropriate paragraphs.

Anti-3 α -diol G antiserum coated tubes: 2 x 50 tubes (ready-to-use)

Plastic tubes with rabbit anti-3 α -diol G immunoglobulin immobilized to the inside wall of each tube.

¹²⁵I-labeled 3 α -diol G tracer (YELLOW): one 55 mL vial (ready-to-use)

At the time of manufacture, the vial contains 185 kBq (<5 μ Ci), of ¹²⁵I-labeled 3 α -diol G in buffer with proteins and sodium azide (<0.1%) and a dye.

Calibrators: one 2.0 mL vial labeled 0 and five 1.0 mL vials labeled 1 - 5 (ready-to-use)

The calibrator vials contain from 0 to approximately 100 ng/mL (0 to approximately 213 nmol/L) of 3 α -diol G in serum and sodium azide (< 0.1%). The exact concentration is indicated on each vial label. After opening, store at 2-8°C up to three weeks, or at < -20°C until expiration date of kit. The calibrator values were established using an internal standard.

Controls: two 1.0 mL vials labeled 1, 2 (ready-to-use)

The vials contain 3 α -diol G in serum and sodium azide (<0.1%). The expected values are indicated in the concentration range indicated in the supplement in the kit. Once opened, store at 2-8°C for up to three weeks, or at < -20°C until expiration date of kit.

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- 12 x 75 mm tubes for Total counts.
- Test tube rack for 12 x 75 mm tubes.
- Precision micropipet (100 μ L).
- Semi-automatic pipet (500 μ L).
- Vortex-type mixer.
- Shaker capable of \geq 180 rpm.
- A sponge rack or similar device for decantation or aspiration system.
- Absorbent material for blotting tubes.
- Gamma counter set for 125 iodine.

PROCEDURE

Preparation of reagents

Let all the reagents come to room temperature and mix them thoroughly by gentle inversion before use.

Assay procedure

Bring all reagents to room temperature before pipeting.

Run calibrators, controls and patient samples in duplicate.

Step 1 Additions	Step 2 Incubation	Step 3 Counting
To antibody coated tubes successively add: 100 µL of calibrator, control or sample and immediately add 500 µL of tracer.* Vortex all tubes.	Cover and incubate 3 hours on a shaker (≥180 rpm) at room temp. (18-25°C). Aspirate or decant all tubes, (except «total cpm» tubes), by simultaneous inversion with a sponge rack into a radioactive waste receptacle.	Strike the tubes sharply on absorbent material to facilitate complete drainage and then allow them to drain on absorbent material for a minimum of 2 minutes. Blot the rim of the tubes. Count bound cpm (B) and total cpm (T) for 1 minute.

*Add 500 µL of tracer to 2 additional uncoated tubes to obtain total cpm.

RESULTS

Results are obtained from the standard curve by interpolation. The curve is used for the determination of 3α-diol G concentrations in samples measured at the same time as the calibrators.

Standard curve

The results in the quality control department were calculated using *spline* curve fit with B/T or B/B_0 on the logit vertical axis and analyte concentration of the calibrators on the log horizontal axis (ng/mL).

Other data reduction methods may give slightly different results.

ED50 = 1.53 ng/mL

Total activity: 63,467 cpm				
Calibrators	3α-diol G (ng/mL)	cpm (n=2)	B/T (%)	B/B ₀ (%)
0	0	37,308	58.78	100.0
1	0.50	28,627	45.11	76.73
2	3.20	12,729	20.06	34.12
3	11.0	6,6940	10.93	18.60
4	31.2	3,877	6.11	10.39
5	115.0	1,714	2.70	4.60

(Example of standard curve, do not use for calculation)

Samples

For each sample, locate the ratio B/T or B/B_0 on the vertical axis of the standard curve and read-off the corresponding 3α-diol G concentration of the sample on the horizontal axis in ng/mL.

To convert concentrations from ng/mL to nmol/L multiply results by 2.13.

EXPECTED VALUES

It is suggested that each laboratory establishes its own normal values. The following values obtained with healthy subjects are indicative only.

	Age range	n	Median (ng/mL)	Absolute range (ng/mL)
All Female	20 - 64	50	1.05	0.15 - 7.53
Females	20 - 29	17	0.84	0.20 - 7.53
	32 - 39	14	1.28	0.19 - 7.40
	40 - 49	10	1.50	0.43 - 3.35
	51 - 64	9	0.78	0.15 - 1.62
All Male	21 - 64	50	6.52	1.68 - 35.3
Males	21 - 28	5	13.62	7.12 - 35.3
	30 - 39	23	5.62	1.68 - 18.8
	40 - 48	15	8.13	2.18 - 27.5
	50 - 64	7	4.34	2.82 - 9.64

QUALITY CONTROL

Good laboratory practices require that control samples be used regularly to ensure the quality of the results obtained. These controls must be processed in exactly the same way as the patient samples, and it is recommended that their results be analyzed using appropriate statistical methods.

Failure to obtain the appropriate values for controls may indicate imprecise manipulations, improper sample handling or deterioration of reagents. In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or imunochem@beckman.com.

PERFORMANCE CHARACTERISTICS

(For more details, see APPENDIX)

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Analytical sensitivity: 0.07 ng/mL

Functional sensitivity: 0.34 ng/mL

Specificity

The antibody used in the immunoassay is highly specific for 3α-diol G. Extremely low cross reactivities were obtained with one related molecule (Androsterone).

Precision

Intra-assay

Samples were assayed 25 times in the same run. The coefficients of variation were ≤ 8.1%.

Inter-assay

Samples were assayed in duplicate in 10 different runs. The coefficients of variation were ≤ 11.2%.

Accuracy

Dilution test

The dilution of patient's samples is not recommended.

Recovery test

Five serum samples containing different levels of endogenous 3α-diol G were spiked with known amounts of 3α-diol G and assayed. The recovery percentages obtained ranged from 83.4% to 111%.

Measurement range (from analytical sensitivity to highest calibrator): 0.07 to approximately 100 ng/mL.

LIMITATIONS

- Failure to follow these instructions for use (IFU) may affect results significantly.
- Failure to blot tubes adequately following decantation may result in poor replication and spurious values.
- Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.
- Repeated freezing and thawing of reagents supplied in the kit and of specimens must be avoided.
- Do not use hemolyzed, lipemic or icteric samples.
- The possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays. Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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ACTIVE is a trademark of BECKMAN COULTER Inc. and its subsidiaries.

SUMMARY AND EXPLANATION OF THE TEST

3 α -androstenediol glucuronide (3 α -diol G, 5 α -Androstane-3 α ,17 β -diol glucuronide) is the glucuronide conjugate of 3 α -androstenediol, a major metabolite of dihydrotestosterone (DHT). 3 α -diol and 3 α -diol G are C19 steroids; the presence of a 17-hydroxyl group (retained from testosterone and DHT) determines their androgenic potency. In bioassays, the potency of 3 α -diol is estimated to be as much as 75% that of testosterone. 3 α -diol is the product of intracellular reduction of DHT, which is enzymatically mediated by 3 β -hydroxysteroid dehydrogenase. A significant proportion of serum 3 α -diol may also be derived from dehydroepiandrosterone sulfate and Δ 4 androstenedione. Although the sites of production are not completely defined, it appears that most of the circulating 3 α -diol in women is derived directly or indirectly from adrenal androgen secretion.

Similarly, the site(s) of 3 α -diol glucuronidation has not been conclusively defined; postulated sites include the skin and liver. 3 α -diol G may also be synthesized directly from glucuronidated precursors, including androsterone glucuronide. Glucuronide conjugation of 3 α -diol appears to facilitate inactivation and urinary excretion of this compound. The measurement of 3 α -diol G has been advocated as a means of assessing skin androgen activity. The androgen sensitivity of skin is dependent upon the presence of the 5 α -reductase enzyme, which converts testosterone to DHT. Evidence suggests that it is DHT binding to the androgen receptor which mediates skin androgen bioactivity. Although 3 α -diol and 3 α -diol G do not appear to have significant direct androgenic action, these compounds may reflect testosterone and DHT production.

Furthermore, interconversion of 3 α -diol and DHT can occur, and this may account for the relative biopotency of 3 α -diol G mentioned above. Elevated levels of serum 3 α -diol G have been reported in congenital adrenal hyperplasia and in various forms of female hirsutism, including idiopathic hirsutism, and polycystic ovary disease. In some cases, serum levels of 3 α -diol G are elevated even in the presence of normal levels of other serum androgens. A decline in serum 3 α -diol G levels has been reported with dexamethasone administration in hirsute women, and measurements of 3 α -diol G may be useful in tests of androgen metabolism.

Therefore, measurement of serum 3 α -diol G levels may be useful in the clinical evaluation and management of hirsutism and other conditions with excessive adrenal androgen production. Previous assays for 3 α -diol G require prior sample extraction and preparation; furthermore, 3 α -diol G levels reported in early literature may have been overestimated due to technical inadequacies. Unlike other 3 α -diol G immunoassays, the DSL method does not require prior sample extraction. The DSL radioimmunoassay for 3 α -diol G uses a specific rabbit polyclonal antiserum with relatively low cross-reactivity to related biological compounds.

Sensitivity

The analytical sensitivity, or minimum detection limit, calculated as 2 standard deviations from the mean of 20 replicates of the 0 ng/mL Androstenediol Glucuronide calibrator, is 0.07 ng/mL.

Specificity

Data on cross-reactivity with several compounds are presented in the following table:

COMPOUND	% CROSS- REACTIVITY
Androstenediol Glucuronide	100.0
Androsterone	0.51
Epiandrosteron Glucuronide	ND
Testosterone	ND
17- β -Estradiol 3-Glucuronide	ND
Androstenedione	ND
Estrone	ND
Estrilol	ND
Progesterone	ND
β -Estradiol	ND
Trans-Dehydroandrosterone	ND
Corticosterone	ND
Etiocolanolone	ND
Cortisone	ND
Dehydroepiandrosterone (DHEA)	ND
5 α -Dihydrotestosterone	ND
DHEA-S	ND
Cortisol	ND

ND = Non-Detectable (<0.1%)

Precision

Intra-assay

The intra-assay precision was determined from the mean of 25 replicates for each of three samples.

Sample	N	Mean (ng/mL)	Coefficient of variation (%)
I	25	1.33	6.23
II	25	20.8	8.05
III	25	56.0	5.83

Inter-assay

The inter-assay precision was determined from the mean of average duplicates for three samples in 10 separate runs.

Sample	N	Mean (ng/mL)	Coefficient of variation (%)
I	10	2.35	5.96
II	10	9.87	6.74
III	10	86.32	11.20

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Recovery test

Five serum samples containing different levels of endogenous 3 α -diol G were spiked with known amounts of 3 α -diol G and assayed.

Sample	Endogen. conc. (ng/mL)	Added conc. (ng/mL)	Measured conc. (ng/mL)	Expected conc. (ng/mL)	Ratio (%) Measured/Expected
S1	1.40	4.42	5.83	5.09	87.35
	1.38	6.51	7.89	6.58	83.43
	1.33	10.45	11.78	11.47	87.35
S2	12.60	4.42	17.02	17.44	102.5
	12.36	6.51	18.87	20.07	106.4
	11.91	10.45	22.36	23.27	104.1
S3	4.35	4.42	8.77	7.46	85.07
	4.26	6.51	10.77	10.96	101.7
	4.11	10.45	14.56	14.73	101.1
S4	9.23	4.42	13.65	14.79	108.3
	9.06	6.51	15.57	17.33	111.3
	8.73	10.45	19.18	19.81	103.3
S5	4.88	4.42	9.31	9.10	97.77
	4.79	6.51	11.30	10.59	93.70
	4.62	10.45	15.07	12.98	86.12

¹²⁵I Characteristics

T1/2 (125I) = 1443 h = 60.14 d

¹²⁵ I Characteristics	E (MeV)	%
γ	0.035	
X	0.027	114
	0.032	25

Symbols Key

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biyolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktiv / Ραδιενεργό / 放射性 / Radioaktivtyoiji medžiaga / Radioaktiv / Radioaktywny / Radioaktivni / Rádioaktívny / 방사성 / Radioaktiv / Радиоактивный / Радиоактивен / 具放射性

Tracer / Tracur / Tracer / Marcato / Trazador / Marcador / Tracer / Ανιχνευτής / 追踪剂 / Ateksamoji medžiaga / Nyomjelző / Znacznik / Radioindikator / Indikator (tracer) / 트레이서 / Tracer /ar / метка / Индикатор / 追蹤劑

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정 물질 / Kalibratör / Калибратор / Калибратор / 校正液

Control / Contrôle / Kontrolle / Controlla / Control / Controla / Kontrolle / Μάρτυρας / 质控品 / Kontrolinė / Kontroll / Kontrola / Kontrola / 컨트롤러 / Kontrol / Контроль / Контролна / 质控品

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνάρια / 试管 / Megintuvėliai / Csövek / Probówki / Zkumavky / Skúmavky / 튜브 / Tüpler / пробирки / Bulgarian / 試管

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用说明

Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferens / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Referenca k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Ürün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználasra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Диагностика in vitro / Ин vitro диагностика / 體外診斷

Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 组成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 制造商 / Gamintojas / Gyártó / Producent / Výrobce / Výrobca / 제조 / Üretici / Изготовлено / Произведено от / 製造商

Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Räckert till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida testi için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Mercado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Знаčka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Sakerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문외 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижете Инструкциите за употреба / 請參閱使用说明

Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 温度范围 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sicaklik aralikları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Ispėjimas / Figyelem / Uwaga / Urozozomení / Urozozomenie / 주의 / Dikkat / Внимание / 注意

Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galijimo data / Lejárati idő / Data ważności / Datum expirace / Dátum expirácie / 만료 날짜 / Son Kullanma Tarihi / Срок годности / Срок на годност / 到期日

Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer serii / Číslo šarže / 로트 번호 / Lot Numarası / Номер партии / Номер на партида / 批號

Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Dátum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期