ACTIVE Androstenedione RIA

REF DSL3800

TABLE OF CONTENTS

English	2	Polski	30
Français	4	Čeština	33
Deutsch	7	Slovenčina	36
Italiano	10	Türkçe	39
Español	13	Русский	42
Português (Portugal)	16	中文 ZH-TW	45
Dansk	19	Србија	47
Ελληνικά	22	APPENDIX	50
中文 ZH-CH	25	REFERENCES	53
Magyar	27		

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ACTIVE Androstenedione RIA

REF DSL3800

RADIOIMMUNOASSAY KIT FOR THE QUANTITATIVE MEASUREMENT OF ANDROSTENEDIONE IN HUMAN SERUM OR PLASMA THIS ASSAY IS INTENDED FOR IN VITRO DIAGNOSTIC USE

PRINCIPLE

The radioimmunoassay of androstenedione (ASD, 4-Androstene-3,17-dione) is a competitive assay. The procedure follows the basic principle of radioimmunoassay where there is a competition between a radioactive and a non-radioactive antigen for a fixed number of antibody binding sites. The amount of [¹²⁵I]-labeled androstenedione bound to the antibody is inversely proportional to the concentration of unlabeled androstenedione present [1]. The separation of the free and bound antigen is achieved by decanting or aspirating the antibody-coated tubes. A standard curve is constructed and unknown androstenedione values are obtained from the curve by interpolation.

For Summary and Explanation of the Test see APPENDIX.

WARNING AND PRECAUTIONS

General remarks:

- For *in vitro* diagnostic use.
- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- Do not use any component beyond the expiration date shown on its label.
- A standard curve must be established with each assay.
- It is recommended to perform the assay in duplicate.
- Calibrators and controls should be mixed before use by inverting or swirling gently rather than vortexing.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material are subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should take place in the presence of radioactive materials.
- No pipeting by mouth.
- Avoid all contact with radioactive materials by using gloves and lab coat.
- All manipulation of radioactive substances should be done in an appropriate location, away from corridors and other busy areas.
- Radioactive materials should be stored in the container provided and in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Material of human origin

Patient samples and blood-derived products may be routinely processed with minimum risk using the procedure described. However, handle these products as potentially infectious according to universal precautions and good clinical laboratory practices, regardless of their origin, treatment, or prior certification. Use an appropriate disinfectant for decontamination. Store and dispose of these materials and their containers in accordance with local regulations and guidelines.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

SDS

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Serum and plasma-EDTA are the recommended sample types.
- Allow serum samples to clot completely before centrifugation.
- Samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage, keep frozen at <-18°C for 1 year maximum, after aliquoting so as to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- Frozen samples should be thawed and mixed thoroughly by gentle swirling or inversion prior to use.
- Any sample reading above the concentration of the highest calibrator should be diluted twice with the 0 ng/mL calibrator and reassayed. Further dilution give cause erroneous results.

Serum and EDTA-plasma values for 20 samples (serum values ranging from 0.50 to 2.07 ng/mL) were compared using the DSL3800 Androstenedione RIA kit. Results are as follows:

[EDTA-plasma] = 1.03 [serum] - 0.03

R = 0.9758

MATERIALS PROVIDED

All unopened reagents in the kit are stable until the expiration date indicated on the kit label, when stored at 2-8°C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Storage conditions for opened reagents are indicated in appropriate paragraphs.

Anti-Androstenedione-Coated Tubes: 2 x 50 tubes (ready-to-use)

Plastic tubes with rabbit anti-androstenedione immunoglobulin immobilized to the inside wall of each tube.

¹²⁵I-labeled Androstenedione Tracer: one 55 mL vial (ready-to-use)

At the time of manufacture, the vial contains 185 kBq (<5 µCi) of ¹²⁵I-labeled androstenedione in buffer with proteins and sodium azide (<0.1%).

Calibrators: one vial labeled 0 and five vials labeled 1-5 (lyophilized)

The calibrator vials contain from 0 to approximately 10.0 ng/mL (0 to approximately 34.9 nmol/L) of androstenedione in human serum with sodium azide (<0.1%). The exact concentration is indicated on each vial label. The volume for reconstitution is indicated on the vial label. After reconstitution, store at 2-8°C for up to 3 weeks or at <-20°C until expiration date of kit.

The calibrator values were established using an internal standard.

Controls: two vials labeled 1, 2 (lyophilized)

The vials contain androstenedione in human serum with sodium azide (<0.1%). The expected values are indicated in a supplement found in the kit. The volume for reconstitution is indicated on the vial label. After reconstitution, store at 2-8°C for up to 3 weeks or at <-20°C until expiration date of kit.

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- 12 x 75 mm plastic or glass test tubes for Total Counts.
- Test tube rack for 12 x 75 mm tubes.
- Precision micropipet (50 µL).
- Deionized water.
- Semi-automatic pipet (500 µL).
- Waterbath, 37°C ± 2°C.
- A sponge rack for decantation or similar device.
- Absorbent material for blotting tubes.
- Gamma counter set for 125 iodine.

RESULTS

Results are obtained from the standard curve by interpolation. The curve is used for the determination of androstenedione concentrations in samples measured at the same time as the calibrators.

Standard curve

The results in the quality control department were calculated using *weighted cubic regression* curve fit with B/T or B/B_0 on the logit vertical axis and analyte concentration of the calibrators on the logarithmic horizontal axis (ng/mL).

$ED_{50} = 1.23$ ng/mL.

Other data reduction methods may give slightly different results.

Total activity: 62,710 cpm				
Calibrators	Andro. (ng/mL)	cpm (n=3)	B/T (%)	B/B ₀ (%)
0	0	23,317	37.2	100
1	0.13	20,425	32.6	87.6
2	0.31	17,375	27.7	74.5
3	1.00	12,542	20.0	53.8
4	3.00	8,284	13.2	35.5
5	10.5	4,498	7.2	19.3

(Example of standard curve, do not use for calculation)

Samples

For each sample locate the ratio B/T or B/B_0 on the vertical axis of the standard curve and read off the corresponding androstenedione concentration of the sample on the horizontal axis in ng/mL. To convert concentrations from ng/mL to nmol/L, multiply results by 3.49.

EXPECTED VALUES

Each laboratory should establish its own reference ranges. Results of a normal range study conducted by an independent laboratory are reported below and are provided for reference only.

Population	N	Median	Min-Max	2.5th - 97.5th percentile
				(ng/mL)
Males	132	1.15	0.62 - 3.12	0.64 - 2.97
Females	99	1.09	0.24 - 3.44	0.35 - 2.78
Postmenopausal Females	50	0.86	0.22 - 2.24	0.30 - 2.07

(For more details, see APPENDIX)

QUALITY CONTROL

Good laboratory practices require that control samples be used regularly to ensure the quality of the results obtained. These controls must be processed in exactly the same way as the patient samples, and it is recommended that their results be analyzed using appropriate statistical methods. Failure to obtain the appropriate values for controls may indicate imprecise manipulations, improper sample handling or deterioration of reagents.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PROCEDURE

Preparation of reagents

Let all the reagents come to room temperature and mix them thoroughly by gentle inversion before use.

Reconstitution of calibrators and control sample

The content of the vials is reconstituted with the volume of deionized water indicated on the vial label. Wait for 10 min following reconstitution and mix gently to avoid foaming before dispensing.

Assay procedure

Bring all reagents to room temperature before pipeting.

Step 1 Additions*	Step 2 Incubation	Step 3 Counting
To antibody coated tubes successively add: 50 µL of calibrator, control or sample, Immediately add 500 µL of tracer. Mix rack vigorously by hand.	Incubate 1 hour at 37°C in water bath.	Aspirate or decant tubes, (except «total cpm» tubes), by simultaneous inversion with a sponge rack into a radioactive waste receptacle. Strike tubes sharply and drain on absorbent material for >2 minutes and gently blot the tubes. Count bound cpm (B) and total cpm (T) for 1 minute.

*Add 500 µL of tracer to 2 additional tubes to obtain «total cpm».

PERFORMANCE CHARACTERISTICS

(For more details, see APPENDIX)

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Analytical sensitivity: 0.05 ng/mL

Functional sensitivity: 0.09 ng/mL

Specificity

The antibody used in the immunoassay is highly specific for androstenedione. Low (<1%) cross reactivities were obtained with several related molecules (Androsterone, 17 OH progesterone, cortisone etc).

Precision

Intra-assay

Serum samples were assayed 25 times in the same run. The coefficients of variation were $\leq 7.5\%$.

Inter-assay

Serum samples were assayed in duplicate in 10 different runs. The coefficients of variation were $\leq 11.3\%$.

Accuracy

Dilution test

High-concentration serum samples were serially diluted with zero calibrator. The recovery percentages ranged from 81.6% to 99.4%.

Recovery test

Low-concentration serum samples were spiked with known quantities of androstenedione. The recovery percentages ranged from 93.0% to 111%.

Measurement range (from analytical sensitivity to highest calibrator): 0.05 to approximately 10.0 ng/mL.

LIMITATIONS

- Failure to follow these instructions for use (IFU) may significantly affect results.
- Failure to blot tubes adequately following decantation may result in poor replication and spurious results.
- Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.
- Avoid repeated freezing and thawing of reagents or specimens.
- Do not use hemolyzed, lipemic or icteric samples.
- The possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays. Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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APPENDIX

PERFORMANCE CHARACTERISTICS

ACTIVE is a trademark of BECKMAN COULTER Inc. and its subsidiaries.

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Summary and explanation of the test

Androstenedione (ASD, 4-Androstene-3,17-dione), a C19 steroid, is produced in the adrenal gland and gonads. ASD is an immediate precursor to both testosterone and estrone, both of which may be subsequently converted to estradiol. Due to the presence of a 17-oxo (rather than hydroxyl) group, ASD has relatively weak androgenic activity, estimated at <20% of testosterone [2]. Although it is a weak androgen, serum ASD levels may exceed testosterone in both normal and disease states, ASD secretion and production rates exceed those of testosterone in women, and significant extra-adrenal conversion of ASD to testosterone occurs. Furthermore, the affinity of sex hormone-binding globulin for ASD is much less than for testosterone or estradiol [2, 3, 4].

The physiologic role of ASD is not well defined. Serum ASD levels are high in fetal and neonatal serum, decrease during childhood, and increase during puberty. In normal pubertal and adult men, the major portion of ASD is derived from the testis, either directly or from conversion of testosterone, while in normal adult women essentially equivalent amounts of ASD are produced by the adrenal gland and ovary [3, 4]. Increased ASD levels may play a role in the development of secondary sexual hair during adrenarche. Serum ASD levels show significant diurnal variation dependent on the secretion of ACTH. Ovarian ASD production is stimulated by luteinizing hormone, and serum ASD levels vary with the menstrual cycle [4]. Adrenal ASD production gradually declines with advanced age in both men and women. In addition, ovarian ASD production decreases after menopause [4].

Measurement of serum ASD provides a useful marker of androgen biosynthesis. Elevated ASD levels have been demonstrated in virilizing congenital adrenal hyperplasia; additionally ASD levels may have advantages over 17-hydroxy-progesterone levels in monitoring treatment of this condition, e.g. less marked diurnal variation and less suppression after brief glucocorticoid exposure [5]. Serum ASD levels are also increased in polycystic ovary syndrome, ovarian stromal hyperthecosis, 3 β -hydroxysteroid dehydrogenase deficiency, and other causes of hirsutism in women [4,6,7]. By definition, ASD levels are normal in idiopathic hirsutism. Elevated serum ASD levels may also occur in adrenal and ovarian virilizing tumors [4].

In a prospective study of over 1000 men, a dose-response relationship of androstenedione and prostate cancer risk was demonstrated [8]; additional studies will be needed to confirm these findings. Assays for ASD include gas-liquid chromatography, mass spectrometry, and immunoassay. This ASD Assay Kit uses a specific and sensitive rabbit anti-human ASD polyclonal antiserum [6, 7, 9].

Sensitivity

The analytical sensitivity, or minimum detection limit, calculated by the interpolation of the mean minus two standard deviations of 10 replicates of the 0 ng/mL Androstenedione calibrator, is 0.05 ng/mL.

Specificity

The cross-reactivity of the androstenedione antiserum has been measured against various compounds. The percent cross-reactivity is expressed as the ratio of the androstenedione concentration to the concentration of the reacting compound at 50% binding of the 0 ng/mL calibrator.

COMPOUND	% CROSS-REACTIVITY	COMPOUND	% CROSS-REACTIVITY
Androstenedione	100	Corticosterone	0.04
Androsterone	0.33	Cortisol	0.04
17OH-Progesterone	0.25	Estrone	0.03
Cortisone	0.16	Pregnenolone	0.02
Isoandrosterone	0.10	Cholesterol	ND
Deoxycortisone	0.09	DHEA-sulfate	ND
Etiocholanolone	0.08	Estradiol	ND
5 α -Dihydrotestosterone	0.08	Estriol	ND
DHEA	0.07	17OH-Pregnenolone	ND
Progesterone	0.06		

ND – not detectable

Precision

Intra-assay

The intra-assay precision was determined from the mean of 25 replicates each with three serum and EDTA plasma samples.

Sample	Serum			EDTA plasma		
	S1	S2	S3	P1	P2	P3
Number of determinations	25	25	25	25	25	25
Mean value, ng/mL	0.51	0.94	8.10	0.55	0.93	7.89
C.V., %	7.46	5.42	5.06	6.14	6.57	5.94

Inter-assay

The inter-assay precision was determined from the mean of average duplicates for 10 separate runs with three serum and EDTA plasma samples.

Sample	Serum			EDTA plasma		
	S1	S2	S3	P1	P2	P3
Number of determinations	10	10	10	10	10	10
Mean value, ng/mL	0.27	0.66	6.71	0.38	1.15	9.18
C.V., %	11.32	6.31	4.12	12.44	7.19	6.70

Accuracy

Dilution test

Five human serum/EDTA-plasma samples were diluted twice with 0 ng/mL androstenedione calibrator and assayed.

Serum	Dilution factor	Measured (ng/mL)	Expected (ng/mL)	Ratio (%) Measured/Expected
S1	-	5.18	-	-
	1:1	2.33	2.59	89.96
S2	-	5.07	-	-
	1:1	2.52	2.54	99.41
S3	-	9.23	-	-
	1:1	3.93	4.62	85.16
S4	-	7.82	-	-
	1:1	3.19	3.91	81.59
S5	-	5.51	-	-
	1:1	2.25	2.76	81.67

EDTA plasma	Dilution factor	Measured (ng/mL)	Expected (ng/mL)	Ratio (%) Measured/Expected
P1	-	4.55	-	-
	1:1	1.92	2.28	84.40
P2	-	4.70	-	-
	1:1	1.98	2.35	84.26
P3	-	9.17	-	-
	1:1	4.18	4.59	91.17
P4	-	6.50	-	-
	1:1	2.85	3.25	87.69
P5	-	6.00	-	-
	1:1	2.53	3.00	84.33

Recovery test

Three serum/EDTA-plasma samples containing different levels of endogenous androstenedione were spiked with known amounts of androstenedione and assayed.

Serum	Endogen. (ng/mL)	Added (ng/mL)	Expected (ng/mL)	Measured (ng/mL)	Ratio (%) Measured/Expected
S1	1.20	0.35	1.55	1.72	110.6
	1.16	0.69	1.85	1.86	100.5
	1.13	1.00	2.13	2.30	108.1
S2	0.94	0.35	1.29	1.30	100.5
	0.91	0.69	1.60	1.59	99.57
	0.88	1.00	1.88	1.75	93.00
S3	0.30	0.34	0.64	0.64	100.2
	0.29	0.66	0.95	0.97	102.4
	0.28	0.95	1.24	1.31	106.0

EDTA plasma	Endogen. (ng/mL)	Added (ng/mL)	Expected (ng/mL)	Measured (ng/mL)	Ratio (%) Measured/Expected
P1	0.21	0.34	0.55	0.60	108.8
	0.21	0.66	0.86	0.86	99.71
	0.20	0.95	1.15	1.32	114.3
P2	0.31	0.34	0.65	0.69	106.4
	0.30	0.66	0.96	1.03	107.7
	0.29	0.95	1.25	1.41	113.2
P3	0.21	0.34	0.55	0.58	105.1
	0.21	0.66	0.86	0.90	104.3
	0.20	0.95	1.15	1.25	108.3

Expected values

Results of normal range studies are reported below and are provided for reference only.

Population	Age range	N	Median	Min - Max	2.5th - 97.5th percentile
(ng/mL)					
Males	20 - 67	132	1.15	0.62 - 3.12	0.64 - 2.97
	20 - 30	31	1.51	0.24 - 2.29	0.65 - 2.20
	31 - 40	38	1.29	0.66 - 2.93	0.67 - 2.56
	41 - 50	37	1.33	0.62 - 2.74	0.74 - 2.61
	51 - 67	26	1.15	0.62 - 3.12	0.64 - 2.97
Females	19 - 62	99	1.09	0.24 - 3.44	0.35 - 2.78
	19 - 30	25	1.47	0.61 - 3.44	0.67 - 3.05
	31 - 40	25	1.12	0.47 - 2.76	0.48 - 2.55
	41 - 49	25	1.17	0.71 - 2.97	0.72 - 2.28
	51 - 62	24	0.72	0.24 - 1.44	0.26 - 1.31
Postmenopausal female		50	0.86	0.22 - 2.24	0.30 - 2.07

Symbols Key

REF Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferens / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Úrün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

IVD In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / Ин витро диагностика / 體外診斷

CONTENTS Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 組成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

 Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 制造商 / Gamintojas / Gyártó: / Producent / Výrobce / 제조 / Üretici / Изготовлено / Произведено от / 製造商

 Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Räckert till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

CE CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Mercado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

SDS Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

 Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문의 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明

 Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度範圍 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Стаклик аралıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

 Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Upozornění / Upozornenie / 주의 / Dikkat / Внимание / 注意

 Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirácie / 만료 날짜 / Son Kullanma Tarihi / Срок годности / Срок на годност / 到期日

LOT Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer seri / Číslo sarže / 로트 번호 / Lot Numarasi / Номер партии / Номер на партида / 批號

 Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Datum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Blyojolik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktiv / Ραδιενεργό / 放射性 / Radioaktivioj medžiaga / Radioaktiv / Radioaktywny / Radioaktivní / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / Tracer / Ανιχνευτής / 追踪剂 / Atsekamoji medžiaga / Nyomjelző / Znacznik / Radioindikátor / Indikátor (tracer) / 트래이서 / Tracer'lar / метка / Индикатор / 追蹤劑

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정 물질 / Kalibratör / Калибратор / Калибратор / 校正液

Control / Contrôle / Kontrolle / Controllo / Control / Controllo / Kontrolle / Μέτρηση / 质控品 / Kontrolliné / Kontroll / Kontrola / Kontrola / Kontrola / 정도관리 / Kontrol / Контроль / Контролна / 質控品

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνάρια / 试管 / Mégintüveliai / Csövek / Probówki / Zkumavky / Skúmvavky / 튜브 / Tüpler / пробирки / Bulgarian / 試管

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用說明

27 September 2017

REFERENCES

1. Yalow R, Berson S: Introduction and general considerations. IN: Odell WD and Daughaday WH, eds. Principles of Competitive Protein Binding Assays. J.B. Lippincott Co., Philadelphia, 1971, pp. 1-19.
2. Dorfman RI, Shipley RA: Androgens. John Wiley and Sons, New York, pp. 116-128, 1956.
3. Horton R, Tait J: Androstenedione production and interconversion rates measured in peripheral blood and studies on the possible site of its conversion to testosterone. J Endocrinol Invest 45:301-313, 1966.
4. Pang S, Riddick L: Hirsutism. IN Lifshitz F (ed): Marcel Dekker, Inc., New York, 1990, pp. 259-291.
5. Cavallo A, Corn C, Bryan GT, Meyer WJ III: The use of plasma androstenedione in monitoring therapy of patients with congenital adrenal hyperplasia. J Pediatr 95:33-37, 1979. Bull NY Acad Med 53, 347, 1977.
6. Rittmaster RS, Thompson DL: Effects of leuprolide and dexamethasone on hair growth and hormone levels in hirsute women: The relative importance of the ovary and adrenal in the pathogenesis of hirsutism. J Clin Endocrinol Metab 70:1096-1102, 1990.
7. Zwicker H, Rittmaster RS: Androsterone sulfate: Physiology and significance in hirsute women. J Clin Endocrinol Metab 76:112-116, 1993.
8. Barrett-Connor E, Garland C, McPhillips JB, Khaw KT, Wingard DL: A prospective population-based study of androstenedione, estrogens and prostatic cancer. Canc Res 50:159-173, 1990.
9. Rosen M, Nouri N, Alexander L: Evaluation of a direct assay for the measurement of androstenedione. Clin Chem 33:891, 1987.