RRA Anti-R TSH

REF A15728, A54464

TABLE OF CONTENTS

English	2	Čeština (cs-cz)	24
Français (France) (fr-fr)	4	Slovák (sk-sk)	27
Deutsche (de-de)	7	한국어 (ko-kr)	30
Italiano (it-it)	10	Türkçe (tr-tr)	32
Español (es-es)	13	Русский (ru-ru)	34
Ελληνικά (el-gr)	16	Srpski (sr-rs)	37
Lietuvių (lt-lt)	19	APPENDIX	39
Polski (pl-pl)	21		

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RRA Anti-R TSH

REF A15728, A54464

RADIORECEPTOR ASSAY FOR THE IN VITRO DETERMINATION OF ANTI-TSH RECEPTOR AUTOANTIBODIES IN HUMAN SERUM

PRINCIPLE

The radioreceptor assay of anti-TSH receptor autoantibodies (TRAb) is a competitive binding assay. Samples, controls and calibrators are incubated in TSH receptor-coated tubes. TRAb in patients' sera are detected by their ability to inhibit the binding of ¹²⁵I-labeled TSH to the receptor coated tubes. A calibration curve is established and unknown values are determined by interpolation from the curve. Alternatively, results can be expressed as an inhibition of TSH binding index.

WARNING AND PRECAUTIONS

General remarks

- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Bring all reagents to room temperature before pipeting.
- Do not mix the reagents from kits of different lots.
- A standard curve must be included with each assay.
- The correct setting of the shaker is very important for the reproducibility of the assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Sodium azide disposal must be in accordance with appropriate local regulations.

Material of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HbsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All patient specimens should be handled as potentially infectious and waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Wash solution (10x)	WARNING	
	H316	Causes mild skin irritation.
	EUH208	May produce an allergic reaction.
	P332+P313	If skin irritation occurs: Get medical advice/attention.
		Tris(hydroxymethyl)–aminomethane 1 - 3%
		2-Chloroacetamide < 0.1%

Safety Data Sheet is available at beckmancoulter.com/techdocs

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in dry tubes.
- Separate serum from cells by centrifugation.
- Serum samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage keep frozen (< -18°C) after aliquoting so as to avoid repeated freezing and thawing.
- If samples have concentrations greater than the highest calibrator, they must be diluted with the negative control (C). Note that due to the heterogeneity of anti-R TSH autoantibodies, the results of dilution test may show non-linear dependence.

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Expiry dates printed on component vial labels apply to the long-term storage by manufacturer only, prior to assembling of the kit. Do not take into account.

Storage conditions for reagents after dilution are indicated in paragraph Procedure.

Kit for the determination of anti-TSH receptor autoantibodies: 60 tubes (REF. #A15728)

TSH receptor-coated tubes: 3 x 20 tubes (ready-to-use)

¹²⁵I-labeled TSH: one 6.5 mL vial (ready-to-use)

The vial contains 100 kBq, at the date of manufacture, of ¹²⁵I-labeled TSH in buffer with BSA, sodium azide (<0.1%), and a dye.

Calibrators: four 0.7 mL vials (ready-to-use)

The calibrator vials contain from 1 to 40 IU/L of anti-TSH receptor autoantibodies in human serum and sodium azide (<0.1%). The exact concentration is indicated on each vial label. The calibrators were calibrated against the reference preparation NIBSC 08/204.

Start buffer: one 10 mL vial (ready to use)

This buffer contains bovine serum albumin and sodium azide (<0.1%).

Control sera: two 0.7 mL vials (ready-to-use)

The vials contain human serum and sodium azide (<0.1%). The positive control contains anti-R TSH antibodies. The expected values are in the concentration range indicated on the vial label.

Wash solution (10x): one 50 mL vial (ready-to-use)

Concentrated solution has to be diluted before use.

Self-seal polythene bag and dessicant for storing opened bags of tubes

Kit for determination of anti-TSH receptor autoantibodies: 100 tubes (REF. # A54464)

TSH receptor-coated tubes: 5 x 20 tubes (ready-to-use)

¹²⁵I-labeled TSH: one 11 mL vial - 180 kBq (ready-to-use)

Calibrators: four 0.7 mL vials (ready-to-use)

Start buffer: one 10 mL vial (ready to use)

Control sera: two 0.7 mL vials (ready-to-use)

Wash solution (10x): one 50 mL vial (ready-to-use)

Self-seal polythene bag and dessicant for storing opened bags of tubes

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- Precision micropipet (100 µL).
- Semi-automatic pipets (100 µL, 50 µL).
- Vortex type mixer.
- Horizontal or orbital shaker (200 rpm).
- Aspiration system.
- Gamma counter set for 125 iodine.

PROCEDURE

Preparation of reagents

Let all the reagents come to room temperature.

Coated-tubes handling

TSH receptor coated tubes must be allowed to reach room temperature (18-25°C for 30 min) before removing from their foil bag(s). After opening a bag of tubes, keep any unused tubes in the foil packet (reseal with adhesive tape), seal in the self-seal plastic bag provided (the bag contains a sachet of desiccant) and store at 2-8°C for up to 1 week. However we recommend that tubes are used on the same day the bag is opened. Tubes which have been stored after opening should not be used in the same assay run as tubes from a freshly opened bag.

Preparation of the wash solution

Pour the content of the vial into 450 mL of distilled water and homogenize. The diluted solution may be stored at 2-8°C until the expiry date of the kit.

Assay procedure

Step 1 Addition 1	Step 2 Addition 2'	Step 3 Counting
To coated tubes, add successively: 50 µL of start buffer and 100 µL of calibrator, control or sample. Mix, cover the tubes. Incubate 120 minutes at 18-25°C with shaking (200 rpm).	Wash twice with 1 mL of wash solution. Aspirate carefully the content of each tubes. To each tube, add: 100 µL of tracer. Mix, cover the tubes. Incubate 60 minutes at 18-25°C with shaking (200 rpm).	Wash twice with 1 mL of wash solution. Aspirate carefully the content of tubes (except the 2 tubes «total cpm».) Count bound cpm (B) and total cpm (T) for 1 min.

* Add 100 µL of tracer to 2 additional tubes to obtain total cpm.

RESULTS

Results are obtained from the standard curve by interpolation. The standard curve serves for the determination of Anti-TSH receptor autoantibodies concentrations in samples measured at the same time as the calibrators.

Standard curve

The results below were calculated using a log-linear curve fit with B/T (%) on vertical axis and the anti-TSH receptor concentration of the calibrators on the horizontal axis (IU/L). Other data reduction methods may give slightly different results.

Total activity : 83,978 cpm				
Calibrators	TRAb (IU/L)	cpm (n=3)	B/T (%)	Inhibition (%)
Control -	0	13,856	16.50	0
1	1	11,665	13.89	16
2	2	10,279	12.24	26
3	8	5,316	6.33	62
4	40	1,822	2.17	87

(Example of standard curve, do not use for calculation. Use the concentration of calibrators indicated on each vial label. The concentrations are lot specific, check carefully.)

The negative control of the kit ("Control -") serves as zero calibrator for plotting the standard curve.

Samples

For each sample, locate the B/T (%), on the vertical axis and read off the corresponding anti-TSH receptor autoantibodies concentration on the horizontal axis. Alternatively, the results can be expressed as inhibition of TSH binding using the formula: Inhibition of TSH binding = $100 \times (1 - B/B_0)$.

EXPECTED VALUES

It is suggested that each laboratory establishes its own normal values. The following values obtained with healthy subjects are indicative only.

	IU/L
Negative	≤ 1
Borderline Positive	$> 1-1.5$
Positive	> 1.5

242 samples from healthy blood donors were assayed and all of them were found negative for TSH receptor autoantibodies.

QUALITY CONTROL

Good laboratory practices imply that control samples be used regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended that their results be analyzed using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PERFORMANCE CHARACTERISTICS

Analytical sensitivity: 0.33 IU/L

Specificity

Analysis of sera from patients with autoimmune diseases other than Graves' disease indicated no interference from autoantibodies to: thyroglobulin, thyroid peroxidase, glutamic acid decarboxylase, 21-hydroxylase, acetylcholine receptor, rheumatoid factor or dsDNA. In addition, human LH up to 10 U/mL, hCG up to 160 U/mL, human FSH up to 15 U/mL and human TSH up to 0.30 U/L had no effect on TSH binding.

Precision

Inter-assay

Samples were assayed in duplicate in 20 different series. The coefficients of variation were found below or equal to 17%.

Intra-assay

Samples were assayed in duplicate in 25 different series. The coefficients of variation were found below or equal to 5.30%.

Measurement range (from analytical sensitivity to the highest calibrator): 0.33 – 40 IU/L.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly. Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Tubes which have been stored after opening should not be used in the same assay run as tubes from a freshly opened bag.

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