# IRMA GH

**REF** IM1397

## TABLE OF CONTENTS

English . . . . .	2	Čeština . . . . .	20
Français . . . . .	5	Slovenčina . . . . .	23
Deutsch . . . . .	8	한국어 . . . . .	26
Italiano . . . . .	11	Русский . . . . .	28
Español . . . . .	14	APPENDIX . . . . .	31
Ελληνικά . . . . .	17		

---

# IRMA GH

**REF** IM1397

## IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION OF GROWTH HORMONE IN HUMAN SERUM AND PLASMA

For *in vitro* diagnostic use.

### PRINCIPLE

The immunoradiometric assay of growth hormone (GH) is a sandwich-type assay. The kit utilizes mouse monoclonal antibodies directed against two different epitopes of the molecule. The antibodies recognize the 22 kDa monomer, the dimer and GH bound to its binding protein. Samples and calibrators are incubated in a tube, coated with the first monoclonal antibody, with a second <sup>125</sup>I-labeled monoclonal antibody. The liquid contents of the tubes are aspirated or decanted after incubation and bound radioactivity is measured. Values are calculated by interpolation from the standard curve. The radioactivity bound is directly proportional to the concentration of GH in the sample.

### WARNING AND PRECAUTIONS

#### General remarks

- The vials with calibrators should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be included with each assay.
- The correct setting of the shaker is very important for the reproducibility of the assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

#### Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

#### Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

#### Material of human origin

All serum and plasma samples should be handled as if capable of transmitting hepatitis or AIDS. Waste should be discarded according to the country rules.

### GHS HAZARD CLASSIFICATION

Calibrator / Control

DANGER

H360	May damage fertility or the unborn child.
P201	Obtain special instructions before use.
P280	Wear protective gloves, protective clothing and eye/face protection.
P308+P313	IF exposed or concerned: Get medical advice/attention. Boric Acid 0.1 - 0.3%

Wash solution (20x)

DANGER

H360	May damage fertility or the unborn child.
P201	Obtain special instructions before use.
P280	Wear protective gloves, protective clothing and eye/face protection.
P308+P313	IF exposed or concerned: Get medical advice/attention. Boric Acid 0.1-0.3% Sodium Borate Decahydrate 0.1-0.3%

**SDS**

Safety Data Sheet is available at [techdocs.beckmancoulter.com](http://techdocs.beckmancoulter.com)

### SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in dry tubes or in tubes with EDTA.
- Separate serum or plasma from cells by centrifugation.
- Serum and plasma samples may be stored at 2-8 °C, if the assay is to be performed within 24 hours. For longer storage keep frozen (at <-20 °C) for up to 3 months after aliquoting so as to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- If samples have concentrations greater than the highest calibrator, they must be diluted in the zero calibrator.

Serum and EDTA plasma values for 15 samples (serum values ranging from 2.46 to 81.38 mIU/L) were compared using the IM1397 IRMA GH kit. Results are as follows:

$$[\text{EDTA-plasma}] = 1.0781[\text{serum}] - 1.055$$

$$R = 0.9933$$

### MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8 °C. Expiry dates printed on component vial labels apply to the long-term storage by manufacturer only, prior to assembling of the kit. Do not take into account.

Storage conditions for reagents after reconstitution or dilution are indicated in paragraph Procedure.

**Anti-GH antibody-coated tubes: 2 x 50 tubes** (ready-to-use)

**<sup>125</sup>I-labeled monoclonal anti-GH antibody: one 11 mL vial** (ready-to-use)

The vial contains 370 kBq, at the date of manufacture, of <sup>125</sup>I-labeled immunoglobulins in liquid form containing bovine serum albumin, sodium azide (<0.1%), and a dye.

**Calibrators: five vials** (lyophilized) + **one 5 mL vial zero calibrator** (ready-to-use)

The calibrator vials contain from 0 to approximately 100 mIU/L of GH, in bovine serum with sodium azide (<0.1% ) for liquid zero calibrator, and in buffer with bovine serum albumine and sodium azide (<0.1%) for the other lyophilized calibrators. The exact concentration is indicated on each vial. The calibrators are calibrated against the international standard WHO 2<sup>nd</sup> IS 98/574 in human serum. 1 mIU correspond to 0.33 µg.

Note: Do not use for recovery tests.

**Control: one vial** (lyophilized)

The vial contains GH in buffer with bovine serum albumin and sodium azide (<0.1%). The expected values are in the concentration range indicated on a supplement.

**Wash solution (20x): one 50 mL vial**

Concentrated solution has to be diluted before use.

## MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- precision micropipet (50 µL).
- semi-automatic pipets (100 µL; 500 µL; 2 mL).
- vortex type mixer.
- horizontal or orbital shaker.
- aspiration system.
- gamma counter set for 125 iodine.

## RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of GH concentrations in samples measured at the same time as the calibrator.

### Standard curve

The results in the package insert were calculated using a log-log curve fit (spline mode) with B/T (%) or B/B<sub>max</sub> (%) on vertical axis and the GH concentration of the calibrators on the horizontal axis (mIU/L). Other data reduction methods may give slightly different results.

Total activity: 147,619 cpm				
Calibrators	GH (mIU/L)	cpm (n=3)	B/T (%)	B/Bmax (%)
0	0	93	0.06	-
1	0.51	829	0.56	0.83
2	2.00	2,689	1.82	2.96
3	10.3	12,677	8.59	14.3
4	51.0	54,005	36.6	61.4
5	100	87,843	59.5	100

(Example of standard curve, do not use for calculation)

### Samples

Locate the ratio B/T or B/B<sub>max</sub> on the vertical axis of the standard curve, read off the GH concentration of the sample on the horizontal axis in mIU/L. In order to convert concentration from mIU/L to ng/mL, multiply results by 0.33.

## EXPECTED VALUES

It is suggested that each laboratory establishes its own normal values.

It is difficult to interpret the basal serum level of GH in view of its circadian variation. However, in the normal adult, a GH concentration of below 20 mIU/L is found usually.

For children, in a stimulation test the lower threshold value indicative of a GH deficiency has been fixed at 20 mIU/L.

## QUALITY CONTROL

Good laboratory practices imply that control samples be used regularly to ensure the quality of the results obtained. These samples must be processed exactly in the same way as the assay samples, and it is recommended that their results be analyzed using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

## PROCEDURE

### Preparation of reagents

Let all the reagents come to room temperature.

### Reconstitution of calibrators and control sample

Calibrators, with the exception of the liquid zero calibrator, and control sample must be brought to the temperature of the laboratory before reconstitution with distilled water, volume is indicated on each vial. Wait at least one-half hour after solubilization before dispensing. **Mix gently to avoid foaming.**

**Do not use a vortex system.** Reconstituted calibrators and control sample may be kept for one week at 2-8 °C or at <-18 °C for longer storage.

### Preparation of the wash solution

Pour the content of the vial into 950 mL of distilled water and homogenize. The diluted solution can be stored at 2-8 °C until the expiry date of the kit.

### Assay procedure

Bring all reagents to room temperature before pipeting.

Step 1 Additions*	Step 2 Incubation	Step 3 Counting
To antibody coated tubes, add successively:  50 µL of calibrator, control or sample and 100 µL of tracer	Incubate 1 hour at 18-25 °C with shaking (>280 rpm).	Aspirate carefully the content of tubes (except the 2 tubes «total cpm») Wash twice with 2 mL of wash solution Count bound cpm (B) and total cpm (T) for 1 minute
Mix		

\*Add 100 µL of tracer to 2 additional tubes to obtain total cpm.

## PERFORMANCE CHARACTERISTICS

(for more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

### Sensitivity

**Analytical sensitivity:** 0.10 mIU/L

**Functional sensitivity:** 0.24 mIU/L

### Specificity

The antibody used in the immunoassay is highly specific for GH. Extremely low cross reactivities were obtained against several related molecules (hCG, HPL, hpGH, prolactin).

### Precision

#### Intra-assay

Serum samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 2.7 %.

#### Inter-assay

Serum samples were assayed in duplicate in 10 different series. The coefficients of variation were found below or equal to 7.1 %.

### Accuracy

#### Dilution test

High-concentration serum samples were serially diluted with the zero calibrator. The recovery percentages obtained were between 81.2 % and 113 %.

#### Recovery test

Low-concentration serum samples were spiked with known quantities of GH. The recovery percentages obtained were between 96.2 % and 107 %.

**Measurement range** (from analytical sensitivity to highest calibrator): 0.10 to approximately 100 mIU/L.

## LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly. Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Do not use hemolyzed, lipemic or icteric samples.

There is no hook effect up to 3000 mIU/L.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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#### Specificity

##### Prolactin

The absence of any interference with the assay by prolactin was verified by the assay of a sample, containing a high concentration of prolactin, to which purified GH was added; 97 % of the GH was recovered.

##### hCG, HPL and placental GH.

The absence of any interference by human chorionic gonadotrophin, human placental lactogen and placenta GH was verified by diluting a human serum with a high level of GH (52.1 mIU/L) into two sera obtained from pregnant women with low GH levels.

	Serum from pregnant women		Serum from pregnant women	
	1 <sup>st</sup> trimester (GH: 1.75 mIU/L)		3 <sup>rd</sup> trimester (GH: 1.76 mIU/L)	
	GH (mIU/L)			
	Expected	Measured	Expected	Measured
Undiluted	-	52.1	-	52.1
1/2	26.9	25.2	26.9	26.3
1/4	14.3	12.5	14.3	12.9
1/8	8.0	7.0	8.0	7.1

#### Precision

##### Intra-assay

Serum	S1	S2	S3
Number of determinations	25	25	25
Mean value (mIU/L)	5.95	24.7	41.2
CV%	2.65	2.62	2.71

Plasma - EDTA	P1	P2	P3
Number of determinations	25	25	25
Mean value (mIU/L)	2.57	11.19	25.8
CV%	4.75	2.46	2.90

##### Inter-assay

Serum	S1	S2	S3	S4	S5
Number of determinations	10	10	10	10	10
Mean value (mIU/L)	5.97	6.57	19.60	32.61	59.80
CV%	5.58	7.12	2.80	6.54	6.56

Plasma - EDTA	P1	P2	P3	P4	P5
Number of determinations	10	10	10	10	10
Mean value (mIU/L)	0.43	18.37	19.66	32.82	64.99
CV%	16.61	3.24	2.98	3.84	4.46

#### Accuracy

##### Dilution test

Serum / EDTA plasma samples were diluted in the zero calibrator and assayed according to the procedure of the kit.

Serum	Dilution	GH (mIU/L)		Ratio (%) Measured/Expected
		Measured	Expected	
S1	-	18.53	-	-
	1:2	9.43	9.27	101.8
	1:4	4.92	4.63	106.2
	1:8	2.57	2.32	111.0
	1:16	1.15	1.16	99.3
	1:32	0.47	0.58	81.2
S2	-	29.61	-	-
	1:2	14.84	14.81	100.2
	1:4	7.83	7.40	105.8
	1:8	4.18	3.70	112.9
	1:16	2.09	1.85	112.9
	1:32	0.88	0.93	95.1
S3	-	55.23	-	-
	1:2	27.09	27.62	98.1
	1:4	12.93	13.81	93.6
	1:8	6.70	6.90	97.0
	1:16	3.35	3.45	97.0
	1:32	1.59	1.73	92.1

EDTA plasma	Dilution	GH (mIU/L)		Ratio (%) Measured/Expected
		Measured	Expected	
P1	-	18.71	-	-
	1:2	10.00	9.36	106.9
	1:4	4.79	5.00	95.8
	1:8	2.51	2.50	100.4
	1:16	1.11	1.25	88.8
	1:32	0.59	0.63	94.4
P2	-	31.59	-	-
	1:2	16.54	15.80	104.7
	1:4	8.35	7.90	105.7
	1:8	4.30	3.95	108.9
	1:16	2.23	1.97	112.9
	1:32	1.07	0.99	108.4
P3	-	62.99	-	-
	1:2	30.06	31.50	95.4
	1:4	14.02	15.75	89.0
	1:8	7.25	7.87	92.1
	1:16	3.87	3.94	98.3
	1:32	1.90	1.97	96.5

#### Recovery test

Known amount of GH was added to serum/EDTA plasma samples and assayed according to the procedure of the kit.

Serum	GH (mIU/L)				Ratio (%) Measured/Expected
	Endogen. conc.	Added conc.	Expected conc.	Measured conc.	
S1	0.75	0.56	1.31	1.37	104.8
	10.73	1.10	1.83	1.93	105.3
	0.71	1.79	2.50	2.68	107.0
S2	5.12	1.28	6.40	6.57	102.7
	5.25	2.81	8.06	8.15	101.1
	5.09	7.26	12.34	12.36	100.1
S3	18.70	5.51	24.21	23.77	98.17
	18.70	11.02	29.73	28.60	96.22
	18.23	17.91	36.15	35.24	97.49

EDTA plasma	GH (mIU/L)				Ratio (%) Measured/Expected
	Endogen. conc.	Added conc.	Expected conc.	Measured conc.	
P1	1.84	0.56	2.41	2.45	101.9
	1.81	1.10	2.91	3.07	105.5
	1.76	1.79	3.55	3.68	103.6
P2	7.36	2.81	10.17	9.96	97.91
	7.22	5.51	12.73	12.61	99.04
	7.04	8.96	16.00	14.62	91.39
P3	25.46	9.25	34.71	34.27	98.74
	24.67	17.91	42.58	43.30	101.7
	23.92	26.05	49.97	48.64	97.33

## <sup>125</sup>I Characteristics

$$T_{1/2} (^{125}\text{I}) = 1443 \text{ h} = 60.14 \text{ d}$$

<sup>125</sup> I	E (MeV)	%
Y	0.035	
X	0.027	114
	0.032	25

## Symbols Key

- REF** Reference Number / Numéro de référence / Referenznummer / Numero di riferimento / Número de referencia / Número de Referència / referensnummer / Αριθμός αναφοράς / 参考编号 / katalogo numeris / referenciaszám / Numer referencyjny / Referenční číslo / Referenčné číslo / 참조 번호 / Referans Numarası / Кат. № / Референтен номер / 参考編號
- IVD** In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Диагностика in vitro / Ин vitro диагностика / 體外診斷
- CONTENTS** Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 组成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄
-  Manufacturer / Fabricant / Hersteller / Produttore / Fabricante / Κατασκευαστής / 制造商 / Gamintojas / Gyártó / Producent / Výrobce / Výrobca / 제조자 / Üretici / Производитель / Производител / 製造業者
-  Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Rácker till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 内容物足夠執行 <n> 次测试
- CE** CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Marcado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識
- SDS** Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adattal / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表
-  Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문의 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижете Инструкциите за употреба / 請參閱使用說明
-  Temperature Limitation / Limites de température / Temperaturbegrenzung / Limiti di temperatura / Limite de temperatura / Limites de temperatura / Temperaturbegränsning / Περιουριοι θερμοκρασίας / 溫度限制 / Temperatūros ribojimai / Hőmérséklet-korlátozás / Ograniczenie temperatury / Teplotní rozmezí od do / Teplotné obmedzenie / 온도 제한 / Sıcaklık sınırlaması / Температурное ограничение / Температурно ограничение / 溫度限制
-  Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Upozornění / Urozornenie / 주의 / Dikkat / Внимание / 注意
-  Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirácie / 만료 날짜 / Son Kullanma Tarihi / Срок годности / Срок на годност / 到期日
- LOT** Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer serii / Číslo šarže / 로트 번호 / Lot Numarası / Номер партии / Номер на партида / 批號
-  Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Dátum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktivt / Ραδιενεργό / 放射性 / Radioaktyvioji medžiaga / Radioaktiv / Radioaktywny / Radioaktivní / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

**DANGER**

DANGER / DANGER / GEFAHR / PERICOLO / PELIGRO / PERIGO / FARA / ΚΙΝΔΥΝΟΣ / 危險 / PAVOJUS / VESZÉLY / NIEBEZPIECZENSTWO / NEBEZPEČÍ / NEBEZPEČENSTVO / 위험 / ТЕΗΛΙΚΕ / ОПАСНО / ОПАСНОСТ / 危險

**Ag** <sup>125I</sup>

**Ab** <sup>125I</sup>

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / Tracer / Ανχνευτής / 追踪剂 / Atsekamoji medžiaga / Nyomjelző / Znacznik / Radioindikátor / Indikátor (tracer) / 트레이서 / Tracer'ar / метка / Индикатор / 追蹤劑

**CAL**

**CAL** **0**

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정물 질 / Kalibratör / Калибратор / Калибратор / 校正液

**CTRL**

Control / Contrôle / Kontrolle / Controllo / Control / Controllo / Kontrolle / Μάρτυρας / 质控品 / Kontroliné / Kontroll / Kontrola / Kontrola / Kontrola / 정도관리 / Kontrol / Контроль / Контролна / 質控品

**TUBE**

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνάρια / 试管 / Mégintüveliai / Csövek / Probówki / Zkumavky / Skúmavky / 튜브 / Túpler / пробирки / Bulgarian / 試管

**IFU**

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用说明

**SOLN** **WASH** **20x**

Wash Solution Concentrate 20X / Solution de lavage concentrée 20X / Waschlösungskonzentrat 20X / Concentrato di soluzione di lavaggio 20X / Solución de lavado concentrada 20X / Concentrado de solução de lavagem 20X / Tvättlösningskoncentrat 20X / Συμπυκνωμένο διάλυμα πλύσης 20X / 浓缩清洗液 20X / Plovimo tirpalo koncentratas 20X / 20X mosóoldat-koncentrátum / Koncentrat 20X roztworu płuczącego / Koncentrát mycího roztoku 20X / Koncentrát premývacieho roztoku 20X / 농축 세척액(20배) / Yıkama Çözeltisi Konsantresi 20X / Концентрат промывочного раствора 20X / Концентрат на разтвор за промиване 20X / 清洗液濃縮 20X