## EIA Histamine

**REF** IM2562

### TABLE OF CONTENTS

English . . . . .	2	Ελληνικά (EL) . . . . .	26
Français (FR) . . . . .	6	Čeština (CZ) . . . . .	31
Deutsch (DE) . . . . .	10	Русский (RU) . . . . .	35
Italiano (IT) . . . . .	14	Србија (SR) . . . . .	40
Español (ES) . . . . .	18	APPENDIX . . . . .	44
Dansk (DA) . . . . .	22		

---

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# EIA Histamine

**REF** IM2562

## ENZYME IMMUNOASSAY FOR THE IN VITRO DETERMINATION OF HISTAMINE IN BIOLOGICAL SAMPLES

For *in vitro* diagnostic use.

### PRINCIPLE

The enzyme immunoassay of histamine is a competition assay. Histamine of the samples, controls or calibrators is chemically modified into acylated histamine during a first step. Acylated histamine is incubated in monoclonal antibody coated wells in presence of alkaline phosphatase acylated histamine conjugate. After incubation, the wells are rinsed in order to remove non-bound components. The bound enzymatic activity is then measured after the addition of a chromogenic substrate.

### WARNING AND PRECAUTIONS

#### General remarks:

- Do not mix the reagents from kits of different lots.
- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- A standard curve must be established with each assay.
- It is recommended to perform the assay in duplicate.
- Bring all reagents to room temperature before pipeting.
- Contact of substrate tablet with skin may adversely affect the product. Handle tablets with forceps.
- Gloves must be worn to avoid contamination by extraneous histamine, e. g. from perspiration.
- Histamine adsorbs on glass, use only plastic pipets, tubes etc.

#### Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Sodium azide disposal must be in accordance with appropriate local regulations.

#### Materials of human origin

All blood samples should be handled as if capable of transmitting hepatitis or AIDS and waste should be discarded according to the country rules.

### GHS HAZARD CLASSIFICATION

Wash Solution (20x) DANGER

H360 May damage fertility or the unborn child.  
 P201 Obtain special instructions before use.  
 P280 Wear protective gloves, protective clothing and eye/face protection.  
 P308+P313 IF exposed or concerned: Get medical advice/attention. Boric Acid 0.1 - 0.3% Sodium Borate Decahydrate 0.1 - 0.3%

Stop Solution (NaOH) DANGER

H314 Causes severe skin burns and eye damage.  
 P280 Wear protective gloves, protective clothing and eye/face protection.  
 P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

Substrate buffer

DANGER

H316 Causes mild skin irritation.  
 H318 Causes serious eye damage.  
 H373 May cause damage to organs through prolonged or repeated exposure.

P280 Wear protective gloves, protective clothing and eye/face protection.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor/physician.

P314 Get medical advice/attention if you feel unwell.

P332+P313 If skin irritation occurs: Get medical advice/attention. Diethanolamine 7 - 10%

Conjugate diluent

DANGER

H360 May damage fertility or the unborn child.  
 P201 Obtain special instructions before use.  
 P280 Wear protective gloves, protective clothing and eye/face protection.  
 P308+P313 IF exposed or concerned: Get medical advice/attention. Boric Acid 0.1 - 0.5%

Conjugate

DANGER

H360 May damage fertility or the unborn child.  
 P201 Obtain special instructions before use.  
 P280 Wear protective gloves, protective clothing and eye/face protection.  
 P308+P313 IF exposed or concerned: Get medical advice/attention. Boric Acid 0.5 - 1%

DMSO

WARNING

H227 Combustible Liquid  
 H315 Causes skin irritation.

	H319	Causes serious eye irritation.
	H335	May cause respiratory irritation.
	P261	Avoid breathing vapours.
	P280	Wear protective gloves, protective clothing and eye/face protection.
	P304+P340	IF INHALED: Remove person to fresh air and keep at rest in a position comfortable for breathing.
	P312	Call a POISON CENTER or doctor/physician if you feel unwell.
	P337+P313	If eye irritation persists: Get medical advice/attention.
	P403+P233	Store in a well-ventilated place. Keep container tightly closed.
		Dimethyl Sulfoxide > 95%
Acylation Buffer	DANGER	
		
	H360	May damage fertility or the unborn child.
	P201	Obtain special instructions before use.
	P280	Wear protective gloves, protective clothing and eye/face protection.
	P308+P313	IF exposed or concerned: Get medical advice/attention. Boric Acid 1 - 5%

 Safety Data Sheet is available at [beckmancoulter.com/techdocs](http://beckmancoulter.com/techdocs)

## SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

### Plasma

- collect blood into a chilled tube containing EDTA only and chill immediately on ice
- centrifuge 10 minutes at 900 g at 4°C within 20 minutes of sample collection
- aspirate the upper 2/3 of plasma
- plasma may be diluted in PBS

### Urine

- Collect urine specimen in plastic container. Either freeze sample or add bacteriostatic agent.
- dilute samples 1:25 in PBS

### Histamine release in whole blood

Do not vortex cell suspension.

Blood sampling: collect 1-2 mL of blood into heparine tubes only, do not use EDTA. Keep samples at 18-25°C for a maximum of 24 hours if the assay cannot be done immediately.

Total histamine after cell lysis: Add 50 µL of undiluted blood to 950 µL of distilled water (dilution 1:20). Freeze and thaw twice.

#### Cell challenge:

- **challenge with histamine release agent:** dilute whole blood 1:7 in histamine release buffer (e. g. 500 µL of whole blood + 3 mL of buffer), add 50 µL of the different dilutions of releasing agent to the wells of a plate, then add 100 µL of the 1:7 dilution of whole blood. Incubate plate covered with lid, for 30 minutes at 37°C (all volumes may be doubled if cells are challenged in tubes).
- **spontaneous histamine release:** add 50 µL of histamine release buffer to wells of a plate, then add 100 µL of the 1:7 dilution of whole blood. Incubate plate covered with lid, for 30 minutes at 37°C.

- **recovery of released histamine:** centrifuge 5 minutes at 900 g at 4°C and keep the plate on crushed ice after centrifugation, carefully collect 100 µL of supernatant fluid without disturbing cell pellet.

### Liquid samples

Collect samples in plastic tubes.

### Solid samples

Extract histamine with 10 µL per mg of tissue of 0.2N HClO<sub>4</sub>, homogenize by sonication or crush and centrifuge at 10 000 g for 5 minutes at < 4°C. Collect supernatant and clarify by filtration if necessary. Neutralize (pH 6-8) supernatant fluid by addition of an equal volume of 1M of potassium borate. Centrifuge at 10 000 g, 1 minute at <4°C. Collect supernatant fluid in plastic tubes.

### Sample storage

Except for cell challenge whole blood, store all type of sample at < -18°C in plastic tubes if the assay cannot be done immediately.

## MATERIALS PROVIDED

Before opening, all reagents of the kit are stable until the expiry date indicated on the kit labels, if stored at 2-8°C.

Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

**Anti-histamine antibody-coated plate: 12 x 8 wells** (ready-to-use)

Unused strips have to be stored at 2-8°C in the self-lock bag provided.

**Histamine-alkaline phosphatase conjugate: two vials** (lyophilised)

The vial contains histamine-alkaline phosphatase conjugate lyophilized in presence of bovine serum albumine. Reconstitute the conjugate with the volume of conjugate diluent stated on the vial label and pour the content of both vials with reconstituted conjugate into Vial for conjugate. After reconstitution the conjugate is stable one week at 2-8°C or at < -18°C until the expiration date of the kit.

**Diluent: one 25 mL vial** (ready to use)

The diluent is used to reconstitute the conjugate. The diluent is stable at 2-8°C until expiration date of the kit.

**Calibrators: six 1 mL vials** (ready-to-use)

The vials contain from 0 to approximately 100 nM of non-acylated histamine, in buffer with bovine serum albumine and sodium azide (<0.1%). After opening, store the vial one week at 2-8°C or at < -18°C until the expiration date of the kit. Calibrators are verified to an internal reference standard.

**Control sample: one 1 mL vial** (ready-to-use)

The vial contains non-acylated histamine, in buffer with bovine serum albumine and sodium azide (<0.1%). The expected values are in the concentration range indicated on the vial label. After opening, store the vial one week at 2-8°C or at < -18°C until the expiration date of the kit.

**Acylation reagent: one vial** (powder)

The content of the vial must be dissolved just before use in the volume of DMSO stated on the label. The dissolved reagent is stable at < -18°C until the expiration date of the kit. Avoid repeated freezing and thawing.

**Acylation buffer: one 5 mL vial** (ready-to-use)

After opening, store the vial at 2-8°C until the expiration date of the kit.

**Histamine release buffer: one vial** (lyophilized)

The buffer contains bovine serum albumine. Reconstitute with 25 mL of distilled water. After opening, store the vial at < -18°C until the expiration date of the kit.

**DMSO: one 3 mL vial** (ready-to-use)

Let the solvent equilibrate at 18-25°C until completely melted.

**Wash solution (20x): one 50 mL vial**

Concentrated solution has to be diluted before use. Stability after dilution: one month at 2-8°C or at < -18°C until the expiration date of the kit.

**Substrate buffer: one 30 mL vial** (ready-to-use)

The substrate buffer is a diethanolamine-HCl solution.

**Substrate: two tablets**

The substrate paranitrophenyl phosphate working solution is prepared at least 30 minutes before use, by dissolving one tablet in 15 mL of substrate

buffer. The solution is stable 24 hours at 2-8°C or stored at < -18°C until the expiration date of the kit.

**Stop solution: one 6 mL vial** (ready-to-use)

This solution is a 1N NaOH solution. The solution is stable at 2-8°C until the expiration date of the kit.

**Vial for conjugate: one vial**

Empty vial for fusion of reconstituted conjugates.

## MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

### For the immunoassay:

- plastic tubes (3 or 5 mL)
- precision micropipets (20-200, 200-1,000 µL)
- repeating micropipets (25, 200 µL)
- "vortex" type mixer
- microtiter plate shaker
- microtiter plate washer (optional)
- microtiter plate reader (405-414 nm)

For histamine release:

- spherical or conical wells microtiter plates (or plastic tubes)
- 37°C water bath
- refrigerated centrifuge for plates (or tubes)

For solid sample processing:

- Perchloric acid (HClO<sub>4</sub>).
- Potassium borate (K<sub>2</sub>B<sub>4</sub>O<sub>7</sub>).

## PROCEDURE

### Preparation of reagents

- Let the reagents come to room temperature.
- Reconstitute the histamine release buffer with 25 mL of distilled water.
- Dissolve just before use acylation reagent with the volume of DMSO stated on the label. If several vials are needed, pool after dissolution.
- Reconstitute the conjugates with the volume of diluent stated on the vial labels. Wait 20 minutes before mixing. Then pour the content of both vials with reconstituted conjugate into Vial for conjugate.
- Pour the content of the wash solution concentrated vial into 950 mL of distilled water and homogenize.
- Prepare the substrate in dissolving one tablet in 15 mL of substrate buffer at least 30 minutes before use.

### Preparation of samples: acylation

The assay procedure includes an acylation step for samples, calibrators or controls. The acylated samples are quite stable and may be stored 3 days at 2-8°C.

### Procedure

Acylation step	Immunological step	Enzymatic step	Reading
To clean plastic tubes, add 25µL of acylation buffer, 100 µL of calibrator, control or sample, 25 µL of acylation solution and vortex immediately.	To antibody coated wells add 50 µL of acylated sample, control or calibrator and 200 µL of conjugate. Incubate 2 hours at 2-8°C with shaking (350 rpm).*	Wash 3 times with 350 µL of diluted wash solution.** Add 200 µL of substrate. Incubate 30 minutes at 18-25°C with shaking (350 rpm)* in the dark.	Add 50 µL of stop solution.  Read absorbance at 405-414 nm.

\* Alternatively incubate 18 hours without shaking at 2-8°C. Then, enzymatic step needs to be reduced to 20 minutes.

\*\* Plate washing: This step is essential to obtain the expected kit performance. After washing, wells must not dry prior to the addition of the next reagent.

Using a microtiter plate washer

Select a three cycle plate washing program that meets the following cycle criteria:

- Fluid in wells must be completely aspirated
- Wells must be filled to the rim with the wash solution
- Wash solution must be injected rapidly (typically one second to fill one well)
- After the three cycles, wash solution must be completely aspirated

Manual procedure

Repeat at least three cycles as follows:

- Turn plate upside-down and shake vigorously over the sink
- Fill wells with wash solution, the solution may run over the rim of the wells
- Turn plate upside-down and shake vigorously over the sink, and firmly tap the inverted microtiter plate onto a clean absorbent paper

## RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of histamine concentration in samples measured at the same time as the calibrators.

### Standard curve

The results in the quality control department were calculated using *spline* curve fit with *B/T* or *B/B<sub>0</sub>* on the logit vertical axis and analyte concentration of the calibrators on the log horizontal axis (nM).

Other data reduction methods may give slightly different results.

Calibrators	Histamine (nM)	Absorbance	B/B <sub>0</sub> (%)
0	0	1.523	100
1	0.90	1.428	93.2
2	2.60	1.269	81.8
3	8.20	0.969	60.1
4	24.0	0.623	35.2
5	83.0	0.328	13.9

(Example of standard curve, do not use for calculation.)

### Samples

- Plasma concentrations correspond to values read off the standard curve.
- Values obtained for released histamine must be multiplied by 10.5.
- Values obtained for total histamine must be multiplied by 20.
- Other samples results must be corrected for the proper dilution factor.

## EXPECTED VALUES

We recommend each laboratory to establish its own reference values. The following values obtained from healthy subjects are indicative only.

- Plasma: < 10 nM (1 ng/mL)
- Whole blood: 200-2,000 nM (20-200 ng/mL)
- Positive cell challenge: >10% of total histamine
- Spontaneous histamine release: < 5% of total histamine
- Urine: 10-35 µg/g of creatinin

## QUALITY CONTROL

Good laboratory practices imply that control samples be used regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended to analyze their results using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

## PERFORMANCE CHARACTERISTICS

(For more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

**Sensitivity**

**Analytical sensitivity:** 0.52 nM

**Specificity**

The immunoassay is specific of histamine.

**Precision****Intra-assay**

EDTA plasma samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 11.0%.

**Inter-assay**

EDTA plasma samples were assayed in duplicate in 10 different series. Coefficients of variation were found below or equal to 13.6%.

**Accuracy****Dilution test**

High-concentration EDTA plasma samples were serially diluted with the zero calibrator. The recovery percentages obtained were between 95.1% and 119%.

**Recovery test**

Low-concentration EDTA plasma samples were spiked with know quantities of histamine. The recovery percentages obtained were between 92.2% and 119%.

**Measurement range** (from analytical sensitivity to highest calibrator):

0.52 to approximately 100 nM.

**LIMITATIONS**

The non-respect of the instructions in this package insert may affect results significantly. Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

**Interference**

Do not use hemolyzed, lipemic or icteric samples.

Pregnant women blood contains high levels of diamine oxydase and cannot be used.

Material containing amine other than histamine at a concentration in excess of 10 mM interferes with the acylation of histamine.

Phenol inhibits histamine release.

Histamine may be a component of allergen such as hymenopteran venom.

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#### Specificity

Data on cross-reactivity with several related molecules are presented in the following table.

Cross-reactivity (%) = related molecule / acylated histamine at 50% of the binding of the zero calibrator.

Related molecule	Cross-reactivities (%)
acylated histamine	100
acylated 3 methylhistamine	0.027
acylated 1 methylhistamine	$2 \times 10^{-3}$
N-acetylimidazole	$<4 \times 10^{-5}$
acylated serotonin	$<4 \times 10^{-5}$
acylated histidine	$<4 \times 10^{-5}$

#### Accuracy

##### Dilution test

5 EDTA plasma samples were after acylation serially diluted by zero calibrator and assayed according to kit procedure.

EDTA plasma	Dilution factor	Expected conc. (nM)	Measured conc. (nM)	Ratio (%) Measured/Expected
1	-	33.89	-	-
	1:2	16.91	16.95	99.79
	1:4	8.48	8.47	100.1
	1:8	4.18	4.24	98.67
	1:16	2.35	2.12	110.9
2	-	55.57	-	-
	1:2	29.53	27.79	106.3
	1:4	15.55	13.89	111.9
	1:8	7.30	6.95	105.1
	1:16	3.85	3.47	110.9
3	-	16.82	-	-
	1:2	8.74	8.41	103.9
	1:4	4.21	4.21	100.1
	1:8	2.35	2.10	111.8
	1:16	1.00	1.05	95.12
4	-	42.06	-	-
	1:2	20.42	21.03	97.10
	1:4	11.00	10.52	104.6
	1:8	6.24	5.26	118.7
	1:16	2.79	2.63	106.1
5	-	26.54	-	-
	1:2	13.37	13.27	100.8
	1:4	7.06	6.64	106.4
	1:8	3.94	3.32	118.8
	1:16	1.90	1.66	114.5

#### Recovery test

5 EDTA plasma samples were spiked by calibrator with higher concentration, acylated and assayed according to kit procedure.

EDTA plasma	Endogen. conc. (nM)	Added conc. (nM)	Measured conc. (nM)	Expected conc. (nM)	Ratio (%) Measured/Expected
P1	4.54	4.20	8.74	9.01	103.1
	4.48	5.19	9.67	9.03	93.39
	4.35	7.55	11.89	12.66	106.5
P2	3.55	1.91	5.47	5.16	94.40
	3.66	4.20	7.87	8.42	107.0
	3.51	7.55	11.05	10.73	97.06
P3	1.81	0.92	2.73	2.79	102.2
	1.75	1.64	3.39	3.21	94.65
	1.78	4.20	5.99	6.29	105.1
P4	2.07	1.22	3.28	3.07	93.46
	2.11	2.68	4.79	5.69	118.9
	2.02	6.15	8.17	7.98	97.71
P5	2.43	1.50	3.93	3.94	100.3
	2.49	3.19	5.68	5.24	92.21
	2.38	6.62	9.00	9.11	101.2

#### Precision

##### Intra-assay

EDTA plasma	P1	P2	P3
Number of determinations	25	25	25
Mean concentration (nM)	1.56	10.51	35.38
C.V., (%)	10.74	11.00	5.28

##### Inter-assay

EDTA plasma	P1	P2	P3
Number of determinations	10	10	10
Mean concentration (nM)	1.81	1.10	17.79
C.V., (%)	12.13	13.62	7.12

## Symbols Key

**REF** Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferens / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Úrün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

**IVD** In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / За ин витро диагностика / 體外診斷

**CONTENTS** Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 組成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

 Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 制造商 / Gamintojas / Gyártó: / Producent / Výrobce / Výrobca / 제조 / Üretici / Изготовлено / Произведено от / 製造商

 Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Räckert till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

**CE** CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Marcado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

**SDS** Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

 Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문의 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明

 Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度範圍 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sıcaklık aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

 Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Upozornění / Upozornenie / 주의 / Dikkat / Внимание / 注意

 Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirácie / 만료 날짜 / Son Kullanma Tarihi / Срок годности / Срок на годност / 到期日

**LOT** Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer serii / Číslo sarže / 로트 번호 / Lot Numarasi / Номер партии / Номер на партида / 批號

 Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Datum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biyolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

**DANGER** Danger / Danger / Gefahr / Pericolo / Peligro / Perigo / Fara / Κίνδυνος / 危險 / Pavojus / Veszély! / Niebezpieczeństwo / Nebezpečí / Nebezpečnost / 위험 / Tehlike / Опасно! / Опасност / 危險

**CAL** Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정 물질 / Kalibrátor / Калибратор / Калибратор / 校正液

**CTRL** Control / Contrôle / Kontrolle / Controllo / Control / Control / Kontrolle / Μέτρηση / 质控品 / Kontrollinė / Kontroll / Kontroll / Kontrola / Kontrola / 컨트롤 / 质控品 / Kontroll / Контроль / Контролна / 质控品

**PLATE** Plate / Plaque / Platte / Piastra / Placa / Platta / Πλάκα / プレート / Plokselē / Lemez / Płytk / Deska / Doštička / 판 / Plaka / Планшет / Чашка / 滴盤

**IFU** Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanna Talimati / Инструкции / Инструкции за употреба / 使用说明

**REAG ACYL** Acylation Reagent, Réactif Acylation, Acylierungsreagenz, Reattivo di acilazione, Reactivo de acilación, Reagente de Acilação, Acylierungsreagens, Αντιδραστήριο ακυλίωσης, 酰化剂, Acilnimo reagentas, Acilezési reagens, Odczynnik do acylacji, Reagencie na acylaci, Acylačné činidlo, 아실화 시약, Açilleme Reaktifi, Реагент для реакции ацилирования, Реактив за ацилиране, 酰化劑

**BUF ACYL** Acylation Buffer, Tampon Acylation, Acylierungspuffer, Tampone di acilazione, Tampón de acilación, Tampão de Acilação, Acylierungspuffer, Ρυθμιστικό διάλυμα ακυλίωσης, 酰化缓冲液, Acilnimo buferinis tirpalas, Acilezési puffer, Bufor do acylacji, Pufr na acylaci, Acilačný tlmivý roztok, 아실화 완충액, Açilleme Tamponu, Буфер для реакции ацилирования, Буфер за ацилиране, 酰化緩衝劑

**RELEASE BUFFER** Release Buffer / Tampon de libération / Freisetzungspuffer / Tampone di rilascio / Tampón de liberación / Tampão de libertação / Frigöringsbuffert / Ρυθμιστικό διάλυμα αποδέσμευσης / 解除緩衝劑 / Buferinio tirpalo išleidimas / Felszabadítási puffer / Bufor uwalniania / Uvoľňovací pufr / Uvoľňovací tlmivý roztok / 해제 버퍼 / Salim Tamponu / Буфер высвобождения / Пуснат на пазара буфер / 釋放緩衝劑

**BUF SUBS** Substrate Buffer / Tampon substrat / Substrat / Tampone per substrato / Tampón Sustrato / Tampão de Substrato / Substrate Buffer / Ρυθμιστικό διάλυμα υποστρώματος / 基質/ツッファー / Substrato buferinis tirpalas / Szubsztrátpuffer / Bufor substratu / Substrátový pufr / Substrátový tlmivý roztok / 기질 완충액 / Substrat Tamponu / Буфер для субстрата / Субстратен буфер / 基質緩衝劑

**DMSO** Dimethyl Sulfoxide, Diméthyle sulfoxyde, Dimethylsulfoxid, Dimetilsolfossido, Dimetilsulfóxido, Dimetilsulfóxido, Dimetilsulfoxid, Διμεθυλοσουλφοξείδιο, 二甲亚砜, Dimetilsulfoksid, Dimetil-sulfoxid, Dimetylosulfoltenek, Dimethylsulfoxid, Dimetyl sulfoxid, 디메틸 술폭사이드, Dimetil Sulfoksit, Диметил сульфоксид, Диметилсулфоксид, 二甲基亞砜

**Diluent** Diluent / Diluant / Diluent / Diluente / Diluyente / Diluente / Späddemedel / Αραιωτικό / 希釈液 / Skiediklis / Diluens / Rozcieńczalnik / Redidlo / Riediaci roztok / 희석액 / Dilüent / Дилуент / Дилуент / 稀釋劑

**CONJ** Conjugate / Conjugué / Konjugat / Coniugato / Conjugado / Konjugat / Συζευγμα / コンジュゲート / konjugatas / Konjugátum / Koniugat / Konjugát / 포합체 / Konjugat / Конъюгат / Конюгат / 結合物

**SOLN WASH 20x** Wash Solution (20x) / Solution de lavage (20x) / Waschlösung (20x) / Soluzione di lavaggio (20x) / Solución de lavado (20 unidades) / Solução de lavagem (20x) / Tvättlösning (20x) / Διάλυμα πλύσης (20x) / 洗浄液 (20x) / Ploviklis (20x) / 20-szoros töménységű mosóoldat / Roztwór płuczący (20x) / Promyvací roztok (20x) / 세척액(20x 농도) / Yıkama Çözeltisi (20x) / Промывочный раствор (концентрированный в 20 раз) / Разтвор за промиване (20x) / 清洗溶液 (20x)

**SOLN STOP** Stopping Solution / Solution d'arrêt / Stopperlösung / Soluzione di arresto / Solución de Detención / Solução de paragem / Διάλυμα αναστολής / 终止液 / Stabdymo tirpalas / Leállító oldat / Roztwór zatrzymujący / Roztok pro zastavení / Zastavovací roztok / 정지 용액 / Durdurma Çözeltisi / Завершающий раствор / Разтвор за спиране / 终止溶液

**SUBS pNPP** Substrate / Substrat / Substrat / Substrato / Sustrato / Substrato / Substrate / Υπόστρωμα / 基質 / 基質 / Substratas / Szubsztrát / substrat / Substrát / 기질 / Substrat / Субстрат / 基材

**VIAL CONJ** Vial for conjugates / Flacon pour les conjugués / Fläschchen für Konjugate / Fiala per coniugati / Vial para conjugados / Frasco para conjugados / Flaska för konjugat / Φιαλίδιο για συζεύγματα / 複合体用バイアル / Buteliukas konjugatams / Konjugátumnak való üveg / Fiolka na koniugaty / Lahvička na konjugát / Flaštička na konjugáty / 집합체용 바이알 / Konjugat flakonu / Флакон для конъюгатов / Флакон за конюгати / 結合物小瓶