INHIBIN A ELISA

REF DSL-10-28100T-1, DSL-10-28100T-4

TABLE OF CONTENTS

English	2	Čeština (CZ)	54
Français (FR)	7	Türkçe (TR)	59
Deutsch (DE)	13	Русский (RU)	64
Italiano (IT)	19	Română (RO)	70
Español (ES)	25	Србија (SR)	75
Português Portugal (PT-PT)	31	APPENDIX	80
Dansk (DA)	37	REFERENCES	84
Ελληνικά (EL)	42		
Magyar (HU)	48		

Beckman Coulter do Brasil Com. e Imp. de Prod. de Lab. Ltda
Alameda Rio Negro, 500, 15º andar, Torre B Alphaville Industrial
CEP 06.454-00 Barueri, São Paulo, Brasil
CNPJ: 42.160.812/0001-44 Telephone: 0800-771-8818

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REF DSL-10-28100T-1, DSL-10-28100T-4

FOR PROFESSIONAL USE ONLY,

The Inhibin A enzyme linked immunosorbent assay (ELISA) kit provides materials for the quantitative measurement of dimeric inhibin A in human serum or plasma. It is intended strictly for *in vitro* use as an aid in the diagnosis and monitoring of various hormonal reproductive disorders. Down syndrome (Trisomy 21) screening in the 2nd trimester using combined biochemical and ultrasound markers may be assessed with appropriate algorithms that can be run on commercially available risk calculation software.^{24, 25} It is strongly recommended to use validated (CE marked) software designated specifically for evaluating the risk of Trisomy 21, e.g. software alpha.* The Inhibin A kit is intended for use in combination with hAFP + hCG + uE3.

SUMMARY

Inhibins are heterodimeric protein hormones secreted by granulosa cells of the ovary in the female and Sertoli cells of the testis in the male. They selectively suppress the secretion of pituitary follicle stimulating hormone (FSH) and also have local paracrine actions in the gonads.^{1,2}

The fully processed form of the inhibin molecule has a molecular weight of approximately 32 kD and consists of the two distinct chains (α and β), linked by disulfide bridges. Higher molecular weight forms, with precursor forms of the α -subunit, also occur in follicular fluid and serum. In addition, free α -subunit forms, unassociated with a β -subunit, and lacking inhibin bioactivity, are also present.^{3,4,5,6}

Inhibin A consists of an α -subunit and a β_A -subunit. Measurements of inhibin A are shown to be useful in studying its role in the human reproductive physiology.^{7,8,9} Several published reports indicate the utility of inhibin A measurement as an endocrine marker for monitoring ovarian function.^{10,11,12,13,14,15,16,17} Until recently, it was not possible to distinguish between circulating functional dimeric inhibin and free α -subunit in the normal human menstrual cycle.

However, by using the highly characterized pair of antibodies, the two-site sandwich ELISA is found to be able to specifically measure only the dimeric inhibin A.¹⁸

PRINCIPLE

The Inhibin A ELISA is an enzymatically amplified two-step "sandwich" assay. In the assay, calibrators, controls and samples are incubated in microtitration wells which have been coated with anti-inhibin β_A subunit antibody. After incubation and washing, anti-inhibin alpha subunit detection antibody labeled with horseradish peroxidase (HRP) is added to each well. After a second incubation and washing step, the substrate tetramethylbenzidine (TMB) is added to the wells.

Lastly, an acidic stopping solution is added. The degree of enzymatic turnover of the substrate is determined by dual wavelength absorbance measurement at 450 nm and between 600 and 630 nm. The absorbance measured is directly proportional to the concentration of inhibin A in the samples. A set of Inhibin A calibrators is used to plot a standard curve of absorbance versus inhibin A concentration. The inhibin A concentrations in the samples can then be calculated from this standard curve.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use.

Use good laboratory practices.¹⁹

Samples and blood-derived products may be routinely processed with minimum risk using the procedure described. However, handle these products as potentially infectious according to universal precautions and good clinical laboratory practices, regardless of their origin, treatment or prior certification.²⁰ Use an appropriate disinfectant for decontamination. Store and dispose of these materials and their containers in accordance with local regulations and guidelines.

GHS HAZARD CLASSIFICATION

Calibrators / Controls WARNING

H317	May cause an allergic skin reaction.
H412	Harmful to aquatic life with long lasting effects.
P273	Avoid release to the environment.
P280	Wear protective gloves, protective clothing and eye/face protection.
P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
P362+P364	Take off contaminated clothing and wash it before use. reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) <0.05%

Antibody Enzyme Conjugate Concentrate

WARNING

H316	Causes mild skin irritation.
H317	May cause an allergic skin reaction.
H412	Harmful to aquatic life with long lasting effects.
P273	Avoid release to the environment.
P280	Wear protective gloves, protective clothing and eye/face protection.
P332+P313	If skin irritation occurs: Get medical advice/attention.
P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
P362+P364	Take off contaminated clothing and wash it before use. 3-(N-Morpholino)-2-hydroxy Propane Sulfonic Acid, Sodium Salt 1 - 2% reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) <0.05%

Sample Buffer A

DANGER

H316	Causes mild skin irritation.
H318	Causes serious eye damage.
P280	Wear protective gloves, protective clothing and eye/face protection.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER or doctor/physician.
P332+P313	If skin irritation occurs: Get medical advice/attention. Sodium Lauryl Sulfate 1 - < 3% Tris(hydroxymethyl)-aminomethane 1 - 2%

Each laboratory should determine the acceptability of its own blood collection tubes and serum separation products. Variations in these products may exist between manufacturers and, at times, from lot-to-lot.

Avoid repeated freezing and thawing of samples.

Avoid assaying lipemic, icteric or hemolyzed samples.

MATERIALS PROVIDED

Kit for the determination of Inhibin A: 96 wells (Cat. #DSL-10-28100T-1)

Anti-Inhibin A Antibody Coated Microtitration strips: One strip holder, containing 96 microtitration wells

Polystyrene microtitration wells with anti-inhibin β_A immobilized to the inside wall of each well.

Store at 2 to 8°C until expiration date in the resealable pouch with a desiccant to protect from moisture.

Inhibin A Antibody-Enzyme Conjugate Concentrate: One 0.6 mL vial

Mouse monoclonal anti-inhibin α -subunit antibody-HRP conjugate, MOPSO, BSA and < 1.0% ProClin** 300.

Store at 2 to 8°C until expiration date.

Dilute prior to use in Conjugate Diluent.

Calibrators: Six 1 mL vials and one 2 mL vial of “zero” calibrator (ready-to-use)

Concentrations of approximately 0, 10, 30, 100, 250, 500 and 1000 pg/mL recombinant dimeric inhibin A in bovine serum with < 0.5% ProClin 300.

Refer to vial labels for exact concentrations.

Store unopened at 2 to 8°C until the expiration date.

Stable at 2 to 8°C for 28 days after initial use. For longer periods store in a -20°C freezer (-15°C to -24°C) until expiration date.

The measurand (analyte) in the Inhibin A calibrators is traceable to the World Health Organization First International Standard for Inhibin-A (code 91/624) containing recombinant 32 kDa human inhibin A (1 pg/mL or 0.037 IU/mL). The assigned values were established using representative samples from this lot of standard and are specific to the assay methodologies of the reagents. Values assigned by other methodologies may be different. Such differences, if present, may be caused by inter-method bias.

Controls: Two 1.0 mL vials (ready-to-use)

Low and high concentrations of recombinant dimeric inhibin A in bovine serum with < 0.5% ProClin 300.

The expected values are in the concentration range indicated on a supplement.

Store unopened at 2 to 8°C until the expiration date.

Stable at 2 to 8°C for 28 days after initial use. For longer periods store in a -20°C freezer (-15°C to -24°C) until expiration date.

Sample Buffer A: One 10.0 mL bottle

Buffer with bovine serum albumin (BSA), animal serum (goat, mouse), surfactant, and sodium azide.

Store at 2 to 8°C until expiration date.

Sample Buffer B: One 10.0 mL bottle

Buffer with < 20% urea hydrogen peroxide, < 0.5% ProClin 300, and sodium azide.

Store at 2 to 8°C until expiration date.

TMB Chromogen Solution: One 11.0 mL bottle

Tetramethylbenzidine (TMB) in citrate buffer with hydrogen peroxide.

Store at 2 to 8°C until expiration date.

Stopping Solution A: One 11.0 mL bottle

0.2 M sulfuric acid.

Store at 2 to 8°C or room temperature (18-25°C) until expiration date.

Conjugate Diluent: One 15.0 mL bottle

Buffer with BSA, animal serum (goat, mouse), surfactant, and < 1.0% ProClin 300.

Store at 2 to 8°C until expiration date.

Wash solution U (20X): One 50 mL vial

Concentrated solution has to be diluted before use.

Store at 2 to 8°C or room temperature (18-25°C) until expiration date.

Kit for the determination of Inhibin A: 384 wells (Cat. #DSL-10-28100T-4)

Anti-Inhibin A Antibody Coated Microtitration strips: Four strip holders, each containing 96 microtitration wells.

Inhibin A Antibody-Enzyme Conjugate Concentrate: Four 0.6 mL vials

Calibrators: Twelve 1 mL vials and two 2 mL vials of “zero” calibrator (ready-to-use)

Controls: Four 1.0 mL vials (ready-to-use)

Sample Buffer A: Four 10.0 mL bottles

Sample Buffer B: Three 10.0 mL bottles

TMB Chromogen Solution: One 50.0 mL bottle

Stopping Solution A: Four 11.0 mL bottles

Conjugate Diluent: Four 15.0 mL bottles

Wash solution (20X): two 50 mL vials

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- Microtitration plate reader capable of absorbance measurement at 450 nm and preferentially capable of dual wavelength correction between 600 and 630 nm
- Deionized Water
- Precision pipette to deliver 50-100 μ L
- Microtitration plate shaker capable of 500–700 orbital revolutions per minute (rpm)
- Microtitration plate washer
- Vortex mixer
- Absorbent materials for blotting the strips
- Graph paper for manual data reduction

PROCEDURE

Procedural notes

- A thorough understanding of this package insert is necessary for successful use of the Inhibin A ELISA.
- It is the responsibility of the customer to validate the assay for their use.
- Reliable results will only be obtained by using precise laboratory techniques and accurately following the package insert.
- A standard curve must be included with each assay.
- Bring all kit reagents to room temperature (18-25°C) before use.
- Thoroughly mix the reagents before use by gentle inversion.
- Do not mix various lots of any kit component within an individual assay.
- Do not use any component beyond the expiration date shown on its label.
- Incomplete washing will adversely affect the outcome and assay precision.
- To minimize potential assay drift due to variation in the substrate incubation time, care should be taken to add the stopping solution into the wells in the same order and speed used to add the TMB chromogen solution.
- Handle all reagents carefully to avoid introducing microbial contaminants which can damage the reagents, especially the Conjugate Diluent and Sample Buffer A.
- Avoid contamination of the TMB chromogen solution with the conjugates.
- Use a clean disposable pipette tip for each reagent, calibrator, control or sample.
- For dispensing sulfuric acid and TMB chromogen solution, avoid pipettes with metal parts.
- The enzyme used as the label is inactivated by oxygen, and is highly sensitive to microbial contamination, sodium azide, hypochlorous acid and aromatic chlorohydrocarbons often found in laboratory water supplies.

- Use deionized water.
- Avoid exposure of the reagents to excessive heat or direct sunlight during storage and incubation.

Preparation of reagents

1. **Wash Solution:** Pour the contents of the vial into 950 mL of distilled water and homogenize. The diluted solution can be stored at 18-25°C one month or at 2-8°C until the expiry date of the kit.
2. **Inhibin A Antibody-Enzyme Conjugate Solution:** The Inhibin A Antibody-Enzyme Conjugate Concentrate should be diluted at a ratio of 1 part into 50 parts of the Inhibin A conjugate diluent, according to the number of wells used. For an entire plate, pipet exactly 220 µL of the Inhibin A antibody-enzyme conjugate concentrate into 11 mL of the conjugate diluent.
NOTE: The Inhibin A antibody-enzyme conjugate concentrate should be freshly diluted 10–15 minutes prior to use.
3. **Microtitration Wells:** Select the number of coated wells required for the assay. The remaining unused wells should be placed in the resealable pouch with a desiccant. The pouch must be resealed to protect from moisture.

Assay procedure

Allow all samples and reagents to reach room temperature (18-25°C). Mix reagents thoroughly by gentle inversion before use. After reconstitution of reagents, mix thoroughly, avoiding foam. Calibrators, controls and samples should be assayed in duplicate.

1. Mark the microtitration strips to be used.
2. Pipet 50 µL of the calibrators, controls and samples to the appropriate wells.
3. Add 50 µL of the Sample Buffer A to each well using a precision pipette.
4. Add 50 µL of the Sample Buffer B to each well using a precision pipette.
5. Incubate the wells, shaking at 500–700 rpm on an orbital microplate shaker, for three hours at room temperature (18-25°C).
6. Prepare the antibody-enzyme conjugate solution by diluting the antibodyenzyme conjugate concentrate in the Inhibin A conjugate diluent as described under the “Preparation of Reagents” section of this package insert.
7. Aspirate and wash each well six times with the wash solution using an automatic microplate washer or manually using a precision pipette. Blot and dry by inverting plate on absorbent material.

NOTE: Use of an automatic microplate washer is strongly recommended. Incomplete washing will adversely affect assay precision. If a microplate washer is not available, follow these steps to wash the plate manually:

- (a) Completely aspirate the liquid from each well
 - (b) Dispense 350 µL of the wash solution into each well using a precision pipette
 - (c) Aspirate the liquid again
 - (d) Repeat steps (b) and (c) five times
8. Add 100 µL of the antibody-enzyme conjugate solution to each well using a precision pipette.
 9. Incubate the wells, shaking at 500–700 rpm on an orbital microplate shaker, for one hour at room temperature (18-25°C).
 10. Aspirate and wash each well six times with the wash solution using an automatic microplate washer. Blot dry by inverting plate on absorbent material.
 11. Add 100 µL of the TMB chromogen solution to each well using a precision pipette.

Avoid exposure to direct sunlight.

12. Incubate the wells, shaking at 500–700 rpm on an orbital microplate shaker, for 15 minutes at room temperature (18-25°C).

NOTE: Be aware that the color may develop more quickly or more slowly than the recommended incubation time depending on the localized room temperature. Visually monitor the color development to optimize the incubation time.

13. Add 100 µL of the stopping solution to each well using a precision pipette.
14. Read the absorbance of the solution in the wells within 30 minutes, using a microplate reader set to 450 nm.

NOTE:

1) While reading the absorbance of the microtitration well, it is necessary to program the zero calibrator as a “Blank”.

2) If wavelength correction is available, set the instrument to dual wavelength measurement at 450 nm with background wavelength correction set between 600 and 630 nm.

RESULTS

1. Calculate the mean absorbance for each calibrator, control or sample.
2. Plot the log of the mean absorbance readings for each of the calibrators along the y-axis versus log of the inhibin A concentrations in pg/mL along the x-axis, using a linear curve-fit. Alternatively, the data can be plotted linear vs. linear and a smoothed spline curve-fit can be used.
3. Draw the best fitting curve through the mean of the duplicate points.
4. Determine the Inhibin A concentrations of the controls and samples from the standard curve by matching their mean absorbance readings with the corresponding Inhibin A concentrations.
5. Any sample reading higher than the highest calibrator should be appropriately diluted using Inhibin A Calibrator 0 and reassayed.
6. Any sample reading lower than the analytical sensitivity should be reported as such.
7. Multiply the value by a dilution factor, if required.

NOTE: If the absorbance readings exceed the limitations of the plate reader, a second reading at 405 nm is needed (reference filter between 600 and 630 nm if available). In this case, proceed to construct a second standard curve as above with the absorbance readings of all calibrators at 405 nm. The concentration of the off-scale samples at 450 nm is then read from the new standard curve. The readings at 405 nm should not replace the on-scale readings at 450 nm.

Standard curve

Calibrators	Conc. (pg/mL)	A	B/B _{max} (%)
0	0.00	0.022 (Blank)	-
1	9.20	0.028	1.021
2	24.3	0.082	2.991
3	90.0	0.305	11.11
4	214	0.733	26.75
5	415	1.302	47.48
6	833	2.741	100.0

(Example of standard curve, do not use for calculation).

Samples

To convert from IU/mL, use the following equation:

$$1 \text{ IU/mL (WHO 91/624)} = 26.7 \text{ pg/mL}$$

EXPECTED VALUES

Each laboratory should establish its own reference ranges to assure proper representation of specific populations. The following values presented were obtained with the Inhibin A ELISA using serum samples from apparently normal healthy adults. All values are reported in pg/mL. For normally cycling females, numbers beneath each phase represent days before or after LH surge.

Population	n	Mean	Median	95% Confidence range
Normally Cycling Females				
Early Follicular Phase (-14 to -10)	136	13.48±0.93	10.54	†5.46-28.16
Mid Follicular (-9 to -4)	228	18.63±0.59	17.06	†7.87-34.54
Late Follicular (-3 to -1)	130	56.10±2.30	52.19	19.49-102.3
Mid Cycle (Day 0)	42	98.87±30.99	98.23	49.92-155.5
Early Luteal (1 to 3)	115	71.59±3.39	67.24	35.93-132.7
Mid Luteal (4 to 11)	268	75.22±2.63	72.51	13.15-159.6
Late Luteal (12 to 14)	82	27.87±3.20	17.44	†7.28-89.95
IVF Peak levels	43	792.4±70.63	705.91	354.2-1,690
PCOS - Ovulatory	26	†9.32±2.61		†5.65-15.99
Postmenopausal	23	†1.15±0.24		†<1-3.88
Normal Males	40	†1.33±0.42		†<1-3.58

† Inhibin A values below Calibrator 1 are extrapolated.

The values represented in the following table are maternal serum inhibin A in the second trimester.

Completed week	No. of samples	Median Inhibin A (pg/mL)	LOG ₁₀ SD
15	30	157.55	0.27
16	124	153.29	0.27
17	74	151.20	0.27
18	45	155.14	0.27
19	20	165.18	0.27
20	16	185.76	0.27

Table adapted from values provided by a US screening center (2005).

QUALITY CONTROL

- Inhibin A ELISA controls or other commercial controls should fall within established confidence limits.
- The confidence limits for Inhibin A ELISA controls are printed on the supplement.
- A full standard curve, plus low and high level controls, should be included in each assay.
- The TMB chromogen solution should be colorless to very light yellow. Development of a blue color may indicate reagent contamination or instability.
- Quality control materials simulate the characteristics of patient samples and are essential for monitoring the system performance of immunochemical assays. Include QC or other commercially available quality control materials that cover at least two levels of analyte. More frequent use of controls or the use of additional controls is left to the discretion of the user based on good laboratory practices or laboratory accreditation requirements and applicable laws. Follow manufacturer's instructions for reconstitution and storage. Each laboratory should establish mean values and acceptable ranges to assure proper performance. Quality control results that do not fall within acceptable ranges may indicate invalid test results.

PERFORMANCE CHARACTERISTICS

(For more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Limit of Detection (LoD): 2.61 pg/mL

Specificity

The antibody used in the immunoassay is highly specific for Inhibin A.

Precision

Intra-assay

Samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 4.6% for serum samples.

Inter-assay

Samples were assayed in duplicate in 10 different series. The coefficients of variation were found below or equal to 6.9% for serum samples.

Accuracy

Dilution test

High-concentration samples were serially diluted with the zero calibrator. The recovery percentages obtained were between 97.8% and 117%.

Recovery test

Low-concentration samples were spiked with known quantities of Inhibin A. The recovery percentages obtained were between 83.5% and 99.5%.

Measurement range (from Limit of detection to highest calibrator): 2.61 to approximately 1,000 pg/mL.

LIMITATIONS

- The reagents supplied in this kit are optimized to measure inhibin A levels in serum or plasma.
- For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the sample. Samples from individuals which have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays. Additionally, other heterophile antibodies such as human antgoat antibodies may be present in samples.^{22,23} Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.
- If there is evidence of microbial contamination or excessive turbidity in a reagent, discard the vial.
- The Inhibin A ELISA results should be interpreted in light of the total clinical presentation of the patient, including: symptoms, clinical history, data from additional tests, and other appropriate information.

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APPENDIX

PERFORMANCE CHARACTERISTICS

Summary and explanation of the test

Sensitivity

Limit of Detection (LoD): 2.61 pg/mL

Specificity

The Inhibin A ELISA assay is highly specific for inhibin A. The following substances did not exhibit cross-reactivity in the Inhibin A ELISA.

HORMONE	CONCENTRATION TESTED
Inhibin B	1 µg/mL
Activin-A	1 µg/mL
Activin-B	1 µg/mL
Pro-Alpha C	0.27 µg/mL

Interference

The following substances in the amounts listed did not significantly interfere with the measurement of inhibin A in the Inhibin A ELISA.

Substance	Concentration Tested	Substance	Concentration Tested
Bilirubin (conjugated)	5 mg/dL	Aminophylline	100 µg/mL
Triglycerides	500 mg/dL	Acetaminophen	100 µg/mL
Hemoglobin	100 mg/dL	Caffeine	100 µg/mL
Follistatin	1 µg/mL	Ibuprofen	100 µg/mL
Alpha 2 Macroglobulin	10 µg/mL	Indomethacin	100 µg/mL
HSA	6 g/dL	Acetylsalicylic Acid	100 µg/mL
Catalase	1 units/mL	Ascorbic Acid	100 µg/mL

Precision

Intra-assay

Serum samples	S1	S2	S3
Number of determinations	25	25	25
Mean value (pg/mL)	47.26	206.4	741.2
C.V., %	4.41	3.20	4.58

EDTA plasma samples	P1	P2	P3
Number of determinations	25	25	25
Mean value (pg/mL)	32.56	305.73	710.76
C.V., %	4.37	4.11	5.24

Inter-assay

Serum samples	S1	S2	S3
Number of determinations	10	10	10
Mean value (pg/mL)	37.36	282.0	374.9
C.V. (%)	6.86	4.87	3.18

EDTA plasma samples	P1	P2	P3
Number of determinations	10	10	10
Mean value (pg/mL)	26.13	77.44	575.38
C.V., %	5.11	7.86	3.59

Accuracy

Dilution test

Three serum and EDTA plasma samples were diluted with 0 pg/mL Inhibin A calibrator and assayed.

Serum samples	Dilution factor	Inhibin A (pg/mL)		Ratio (%) Measured/Expected
		Measured	Expected	
S1	-	720.1	-	-
	1:2	367.6	360.3	102.0
	1:4	201.7	180.1	112.0
	1:8	100.7	90.07	111.9
S2	1:16	47.12	45.04	104.6
	-	847.2	-	-
	1:2	414.1	423.6	97.76
	1:4	230.3	211.8	108.8
S3	1:8	112.2	105.9	106.0
	1:16	57.06	52.95	107.8
	-	644.6	-	-
	1:2	348.4	322.3	108.1
S3	1:4	188.1	161.2	116.7
	1:8	87.05	80.58	108.0
	1:16	42.13	40.29	104.6

EDTA plasma samples	Dilution factor	Inhibin A (pg/mL)		Ratio (%) Measured/Expected
		Measured	Expected	
P1	-	704.0	-	-
	1:2	362.4	352.0	103.0
	1:4	187.5	176.0	106.5
	1:8	93.89	88.00	106.7
P2	1:16	45.20	44.00	102.7
	-	678.6	-	-
	1:2	304.2	339.3	89.66
	1:4	154.6	169.6	91.11
P3	1:8	76.73	84.82	90.46
	1:16	35.21	42.41	83.00
	-	484.8	-	-
	1:2	239.5	242.4	98.79
P3	1:4	107.3	121.2	88.54
	1:8	53.30	60.60	87.96
	1:16	25.09	30.30	82.81

Recovery test

Known amount of Inhibin A was added to 5 serum or plasma samples and assayed according to the procedure of the kit.


Serum samples	Endog. conc. (pg/mL)	Added conc. (pg/mL)	Expected conc. (pg/mL)	Measured conc. (pg/mL)	Ratio (%) Measured/Expected
S1	14.77	6.18	20.96	20.86	99.54
	14.03	17.51	31.54	28.54	90.50
	14.64	36.29	50.92	46.44	91.19
S2	29.84	9.74	39.57	38.52	97.34
	30.24	29.62	59.86	55.34	92.44
	29.39	57.57	86.95	81.24	93.43
S3	32.01	12.97	44.98	42.20	93.82
	32.59	39.59	72.18	66.45	92.07
	31.37	76.20	107.57	94.09	87.47
S4	41.45	16.81	58.26	55.13	94.63
	42.43	51.65	94.08	87.66	93.18
	40.64	92.45	133.09	116.59	87.60
S5	15.48	5.96	21.44	17.90	83.48
	14.71	17.11	31.82	27.70	87.04
	15.33	35.81	51.14	43.18	84.43

EDTA plasma	Endog. conc. (pg/mL)	Added conc. (pg/mL)	Expected conc. (pg/mL)	Measured conc. (pg/mL)	Ratio (%) Measured/Expected
P1	579.6	9.42	589.0	553.2	93.92
	587.2	28.66	615.9	590.0	95.80
	571.1	55.75	626.9	618.6	98.69
P2	25.31	8.03	33.34	35.34	106.0
	25.59	24.32	49.91	50.38	100.9
	25.00	47.51	72.50	72.49	99.99
P3	60.37	17.01	77.38	75.24	97.23
	61.83	52.10	113.9	111.6	97.94
	59.25	92.45	151.7	146.5	96.57
P4	48.29	16.81	65.10	71.01	109.1
	49.44	51.65	101.1	92.69	91.70
	47.35	92.45	139.8	135.6	96.97
P5	351.8	5.96	357.8	320.6	89.61
	334.4	75.33	409.7	383.8	93.67
	348.4	35.81	384.2	427.3	111.2

Expected values

Figure 1: Physiological Profiles of Inhibin A, Estradiol and Progesterone During the Normal Ovulatory Cycle

The mean ± SE concentrations of inhibin A, estradiol, and progesterone in 35 subjects (9 subjects with two cycles each, total 44 normal cycles) are shown, aligned relative to the day of the mid cycle LH peak. The inhibin A levels remained low (13.6 pg/mL) during the early follicular phase and rose in the late follicular phase to reach a peak concentration (98.9 pg/mL) on the day of the LH peak. After the LH peak, inhibin A levels fell briefly, then rose again during early to mid luteal phase, and fell again in concert with the drop in progesterone during the late luteal phase (28.3 pg/mL). Estradiol and inhibin A were highly correlated during the follicular phase (day -14 to -2) as evident by the following linear regression data:


Inhibin A = 0.18 [Estradiol] + 3.3	Inhibin A range†: 3.0 - 82.8 pg/mL
n = 433	Estradiol range: 5.0 - 359 pg/mL
r = 0.80	95% Confidence interval for slope = 0.17 - 0.19
p < 0.001	95% Confidence interval for intercept = 1.8 - 4.8

†Inhibin A values below Calibrator 1 are extrapolated.

Figure 2: Inhibin A and Estradiol Profiles in Ovulation Induction Cycles

The mean ± SE concentrations of inhibin A and estradiol in 20 subjects during ovarian stimulation for IVF are shown, aligned to the day of hCG treatment.

Serum and Plasma Comparisons

Figures 3 and 4 compare inhibin A levels in paired serum vs. EDTA plasma and serum vs. Heparin plasma samples, respectively.



Figure 4 - Serum vs. Heparin Plasma

Comparison of Inhibin A and Estradiol

Two hundred forty-six samples were collected from subjects scheduled for IVF at two different centers. Inhibin A was measured by ELISA and estradiol was measured by two separate commercially available kits. Regression analysis yielded the following:

Inhibin A = 0.247 [Estradiol] + 67.6	Inhibin A range† = 2.2 - 1,395 pg/mL
n = 246	Estradiol range = 10 - 5,960 pg/mL
r = 0.85	95% Confidence interval for slope = 0.227 - 0.267
p < 0.0001	95% Confidence interval for intercept = 44.1 - 91.2

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Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

*alpha is a trademark of Logical Medical Systems Ltd.

**ProClin™ is a trademark of The Dow Chemical Company ("Dow") or an affiliated company of Dow.

Symbols Key

REF Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferens / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Úrűn Referans / Ссылка на продукт / Референция за продукт / 產品參考

IVD In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnostikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / За ин витро диагностика / 體外診斷

CONTENTS Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 組成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

 Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 製造商 / Gamintojas / Gyártó: / Producent / Výrobce / Výrobca / 제조 / Üretici / Изготовлено / Произведено от / 製造商

 Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Räckert till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

CE CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Mercado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

SDS Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

 Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 請參閱使用說明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문의 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明

 Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo(i) di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度範圍 / Temperaturüros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Szaklık aralliklat / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

 Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Upozornění / Upozornenie / 주의 / Dikkat / Внимание / 注意

 Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirácie / 만료 날짜 / Son Kulanma Tarihi / Срок годности / Срок на годност / 到期日

LOT Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer seri / Číslo sarže / 로트 번호 / Lot Numarasi / Номер партии / Номер на партида / 批號

 Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Datum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fare/ Βιολογικός κίνδυνος / Veszélyes biológiai anyag / Biologické riziko/ Βιολογική τεχληία / Биологическая опасность / Pericol biologic/ biološki rizik

CONJ CONC Conjugate Concentrate/ Conjugué concentré/ Concentrato di coniugato/ Concentrado de conjugado/ Concentrado de conjugado/ Konjugatkoncentration/ Συμπύκνωμα συζε ύματος/ Konjugátum koncentrátum/ Koncentrát konjugátu / Konjugat Konsantre/ Концентрат на конюгат/ Conjugat concentrat/ Koncentrat konjugata

CONJ DIL Conjugate Diluent/ Diluant pour le conjugué/ Diluente del coniugato/ Diluyente de conjugado/ Diluente de conjugado / Konjugatfortynder/ Αραιωτικό συζευγματος/ Konjugátum-oldószér/ Diluent pro konjugát/ Konjugat Dilüent/ Разбавитель для конъюгата/ Solvent conjugat/ Razredivač konjugata

CAL 0 Calibrator/ Calibrateur/ Calibratore/ Calibrador/ Calibrador / Kalibrator/ Βαθμονομητής/ Kalibrátor/ Kalibrátor/ Kalibratör / Калибратор/ Calibrator/ Kalibrator

CAL

CTRL Control/ Contrôle/ Controllo/ Control/ Control/ Kontrol / Δεδομένα/ Kontroll/ Kontrolni vzorek/ Kontrol/ Контроль / Ser de control/ Kontrola

SOLN STOP Stopping Solution/ Solution d'arrêt/ Soluzione di arresto/ Solución de Detención/ Solução de parada/ Stopopløsning / Διάλυμα αναστολής/ Leállító oldat/ Stop roztok/ Durdurma Çözeltisi/ Завершающий раствор/ Soluție de oprire/ Rastvor za zaustavljanje

SOLN TMB TMB Chromogen Solution/ Solution chromogène TMB/ Soluzione cromogena TMB/ Solución de cromógeno TMB/ Solução de cromogénio TMB/ TMB-kromogenopløsning/ Διάλυμα χρωμογόνου TMB/ TMB kromogén oldat/ Roztok chromogenu TMB / TMB Kromojen Çözeltisi/ Раствор хромогена TMB/ Soluție cromogenă TMB/ Rastvor TMB hromogena

SOLN WASH 20x Wash Solution/ Solution de lavage/ Soluzione di lavaggio / Rastvor za pranje/ Solução de lavagem/ Skyløpløsning / Διάλυμα Πλύσης/ Mosooldat/ Promyvaci roztok/ Yıkama Çözeltisi/ Промывочный раствор/ Soluție de spălare/ Rastvor za pranje

SAMPLE BUF A Sample Buffer/ Tampon pour échantillon/ Tampone per campione/ Tampón de muestra/ Tampão de amostra/ Prøvebuffer/ Ρυθμιστικό διάλυμα δειγμάτων/ Minta puffer/ Vzorkový pufr/ Numune Tamponu/ Буфер для проб/ Eşantion soluție tampon/ Uzorak pufera

SAMPLE BUF B

PLATE Plate/ Plaque/ Piastra/ Placa/ Placa/ Plade/ Πλάκα/ Lemez / Deska/ Plaka/ Планшет/ Placă/ Pločica

IFU Instruction For Use/ Mode d'emploi/ Istruzioni per l'uso / Instrucciones de uso/ Instruções de utilização/ Bruksanvisning / Οδηγίες χρήσης/ Használati utasítás/ Návod k použití / Kullanma Talimatı/ Инструкции по применению / Instrucțiuni de utilizare/ Uputstvo za upotrebu

REFERENCES

1. Vale WW, Hseuh A, Rivier C, Yu J. 1990. The inhibin/activin family of hormones and growth factors. In: Sporn MA and Robert AB, eds. Peptide growth factors and their receptors: handbook of experimental physiology. Vol 95. Berlin: Springer-Verlag; 211-248.
2. Burger HG. 1992. Inhibin. *Reprod Med Rev.* 1:1-20.
3. Knight PG, Beard AJ, Wrathall HM, Castillo RJ. 1989. Evidence that the bovine ovary secretes large amounts of monomeric inhibin α -subunit and its isolation from bovine follicular fluid. *J Mol Endocrinol* 2:189-200.
4. Schneyer AL, Sluss PM, Whitcomb RW, Martin KA, Sprengel R, Crowley WF Jr. 1991. Precursors of α inhibin modulate FSH receptor binding and biological activity. *Endocrinology* 129:1987-1999.
5. Robertson DM, Sullivan J, Watson M, Cahir N. 1995. Inhibin forms in human plasma. *J Endocrinol* 144:261-269.
6. Robertson D, Burger HG, Sullivan J, Cahir N, Groome N, Poncelet E, Franchimont P, Woodruff T, Mather JP. 1996. Biological and immunological characterization of inhibin forms in human plasma. *J Clin Endocrinol & Metab.* 81:669-676.
7. Lockwood GM, Muttukrishna S, Ledger WL. 1998. Inhibins and activins in human ovulation, conception and pregnancy. *Human Reproduction Update* 4:284-295.
8. Groome NP, Illingworth PJ, O'Brien M, Cooke I, Ganesan TS, Baird DT, and McNeilly AS. 1994. Detection of dimeric inhibin throughout the human menstrual cycle by two-site enzyme immunoassay. *Clin Endocrinol* 40:717- 723.
9. Lambert-Messerlian, G, Hall JE, Sluss PM, Taylor AE, Martin KA, Groome NP, Crowley WF Jr, Schneyer A. 1994. Relatively low levels of dimeric inhibin circulate in men and women. *J Clin Endocrinol Metab* 79:45-50.
10. Muttukrishna S., Fowler PA, Groome NP, Mitchell GC, Robertson WR and Knight PG. 1994. Serum concentrations of dimeric inhibin during the spontaneous human menstrual cycle and after treatment with exogenous gonadotrophin. *Hum Reprod* 9:1634-1642.
11. Lockwood GM, Muttukrishna S, Groome NP, Knight PG, Ledger L. 1996. Circulating inhibins and activin-A during GnRH down-regulation and ovarian hyperstimulation with recombinant FSH for in-vitro fertilization embryo transfer. *Clinical Endocrinol* 45:741-748.
12. Lindheim SR, Chang PL, Vidali A, Ferin M, Sauer MV. 1998. The utility of progesterone and inhibin A for monitoring natural-cycle IVF-ET. *J Assist Reprod Genet* 15:538-541.
13. Rombauts L, Verhoeven G, Meuleman C, Koninckx PR, Poncelet E, Franchimont P. 1996. Dimeric inhibin A and alpha-subunit immunoreactive material in maternal serum during spontaneous and in vitro fertilization pregnancies. *J Clin Endocrinol Metabol* 81:985-989.
14. Hall JE, Welt CK, Cramer DW. 1999. Inhibin A and inhibin B reflect ovarian function in assisted reproduction but are less useful at predicting outcome. *Hum Reprod* 14:409-415.
15. Ehret W, Heil W, Schmitt Y, Töpfer G, Wisser H, Zawta B, et al. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. *WHO/DIL/LAB/99.1 Rev.2:22,46pp.*
16. Crofton PM, Illingworth PJ, Groome NP, Stirling HF, Swanston I, Gow S, Wu FCW, McNeilly A, Kelnar CJH. 1997. Changes in dimeric inhibin A during normal early puberty in boys and girls. *Clin Endocrin* 46:109-114.
17. Lockwood GM, Ledger WL, Barlow DH, Groome NP, Muttukrishna S. 1997. Measurement of inhibin and activin in early human pregnancy: demonstration of fetoplacental origin and role in prediction of early pregnancy outcome. *Biol Reprod* 57:1490-1494.
18. Groome NP, O'Brien M. 1993. Immunoassays for inhibin and its subunits. Further applications of the synthetic peptide approach. *J Immunol Methods* 165:167-176.
19. Clinical Laboratory Safety; Approved Guideline Second Edition, GP17-A2. 2004. Clinical and Laboratory Standards Institute.
20. HHS Publication, 5th ed., December 2009. Biosafety in Microbiological and Biomedical Laboratories. Available <http://.cdc.gov/biosafety/publications/bmb15>
21. Approved Guideline – Procedures for the Handling and Processing of Blood Specimens, H18-A3. 2004. Clinical and Laboratory Standards Institute.
22. Kricka, L. 2000. Interferences in immunoassays - still a threat. *Clin Chem* 46:1037.
23. Bjerner J, Nustad K, Norum LF, Olsen KH, Børner OP. 2002. Immunometric assay interference: incidence and prevention. *Clin Chem* 48:613-621.
24. Maymon R and Shulman A. 2004. Integrated first- and second-trimester Down syndrome screening test among unaffected IVF pregnancies. *Prenat. Diagn.* 24: 125-129.
25. Wald NJ, Bestwick JP, Morris JK. 2006. Cross-trimester marker ratios in prenatal screening for Down syndrome. *Prenat. Diagn.* 26: 514-523.