LH IRMA KIT

REF IM1381, IM3302

TABLE OF CONTENTS

English	2	Čeština	29
Français	5	Slovenčina	32
Deutsch	8	한국어	35
Italiano	11	Türkçe	37
Español	14	Русский	40
Ελληνικά	17	中文 ZH-TW	43
Lietuviškai	20	Србија	45
Magyar	23	APPENDIX	48
Polski	26		

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LH IRMA KIT

REF IM1381, IM3302

IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION OF LUTEINIZING HORMONE IN HUMAN SERUM AND PLASMA For *in vitro* diagnostic use.

PRINCIPLE

The immunoradiometric assay of luteinizing hormone (LH) is a sandwich-type assay. The kit utilizes mouse monoclonal antibodies directed against two different epitopes of LH and hence not competing. Samples or calibrators are incubated in tubes, coated with the first monoclonal antibody in the presence of the second ¹²⁵I labeled monoclonal antibody. The content of tubes is aspirated and rinsed after incubation and bound radioactivity is measured. Values are calculated by interpolation from the standard curve. The radioactivity bound is directly proportional to the concentration of LH in the sample.

WARNING AND PRECAUTIONS

General remarks:

- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be established with each assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipeting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- Radioactive materials should be stored in the container provided in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Every case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Materials of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HBsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All serum and plasma samples should be handled as if capable of transmitting hepatitis or AIDS and waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Wash Solution (20X) DANGER

H360	May damage fertility or the unborn child.
P201	Obtain special instructions before use.
P280	Wear protective gloves, protective clothing and eye/face protection.
P308+P313	IF exposed or concerned: Get medical advice/attention. Boric Acid 0.1-0.3% Sodium Borate Decahydrate 0.1-0.3%

SDS Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in tubes containing EDTA or no additives.
- Separate serum or plasma from cells by centrifugation.
- Serum and plasma samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage keep frozen (< -20°C, 1 year maximum) after aliquoting so as to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- If samples have concentrations greater than the highest calibrator, they must be diluted in the zero calibrator.

Serum and EDTA plasma values for 16 samples (serum values ranging from 2.26 to 40.70 IU/L) were compared using the IM1381 LH IRMA kit. Results are as follows:

$$[\text{EDTA-plasma}] = 1.0846 [\text{serum}] - 0.107$$

$$R = 0.9993$$

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Storage conditions for reagents after reconstitution or dilution are indicated in paragraph Procedure.

Kit for determination of LH: 100 tubes (Cat. #IM1381)

Anti-LH monoclonal antibody coated tubes: 2 x 50 tubes (ready-to-use)

¹²⁵I-labeled monoclonal anti-LH antibody: one 5.5 mL vial (ready-to-use)

The vial contains 370 kBq, at the date of manufacture, of ¹²⁵I-labeled immunoglobulins in buffer with bovine serum albumin, sodium azide (<0.1%), and a dye.

Calibrators: six vials (lyophilized)

The calibrator vials contain from 0 to approximately 180 IU/L of LH in horse serum with sodium azide (<0.1%). The exact concentration is indicated on each vial. The calibrators were calibrated against the international standard, WHO 2nd IS 80/552.

Control serum: one vial (lyophilized)

The vial contains human LH lyophilized in human serum. The expected values are in the concentration range indicated on a supplement.

Wash solution (20x): one 50 mL vial

Concentrated solution has to be diluted before use.

Kit for determination of LH: 400 tubes (Cat. # IM3302)**Anti-LH monoclonal antibody coated tubes: 8 x 50 tubes** (ready-to-use)**¹²⁵I-labeled monoclonal anti-LH antibody: four 5.5 mL vials** (ready-to-use)**Calibrators: six vials** (lyophilized)**Control serum: one vial** (lyophilized)**Wash solution (20X): two 50 mL vials****MATERIALS REQUIRED, BUT NOT PROVIDED**

In addition to standard laboratory equipment, the following items are required:

- Precision micropipets (100 µL).
- Semi-automatic pipets (50 µL; 2 mL).
- Vortex type mixer.
- Horizontal or orbital shaker.
- Aspiration system.
- Gamma counter set for 125 iodine.

RESULTS

Results are obtained from the standard curve by interpolation. The standard curve serves for the determination of LH concentrations in samples measured at the same time as the calibrator.

Standard curveThe results in the quality control department were calculated using *spline* curve fit with determined radioactivity ($cpm_{cal} - cpm_{cal0}$) on the log vertical axis and analyte concentration of the calibrators on the log horizontal axis (IU/L).

Other data reduction methods may give slightly different results.

Calibrators	LH (IU/L)	cpm (n=3)	B/T (%)	$cpm_{cal} - cpm_{cal0}$
0	0	101	0.09	-
1	0.47	282	0.24	181
2	2.40	1,170	1.00	1,069
3	14.7	6,310	5.39	6,209
4	60.0	22,674	19.4	22,573
5	170	49,663	42.4	49,562

*(Example of standard curve, do not use for calculation)***Samples**

For each sample, locate the cpm or B/T value on the vertical axis and read off the corresponding analyte concentration on the horizontal axis.

The concentrations of diluted samples must be corrected by the dilution factor.

EXPECTED VALUES

We recommend each laboratory to establish its own reference values. LH concentrations measured in serum samples from healthy subjects in the following table are indicative only.

	N	Min.	Max.	Median	2.5 th percentile	97.5 th percentile
		IU/L				
Men	98	0.81	9.05	3.32	1.04	8.02
Women						
Follicular phase	351	0.77	80.3	6.65	2.81	20.4
Preovulatory peak	70	11.9	112	40.7	13.9	96.2
Luteal phase	309	0.39	38.5	7.17	1.04	24.9
Postmenopausal	48	9.03	70.6	29.5	10.5	58.8

(For more details, see APPENDIX)

Children	LH (IU/L)			
	Girls		Boys	
	Min.	Max.	Min.	Max.
Pre-pubertal		<0.19		<0.6
Stage of puberty PII		<0.6		<2.4
PIII, PIV, PV	0.32	2.7	0.24	5.9

QUALITY CONTROL

Good laboratory practices imply that control samples must be used regularly to ensure the quality of the results obtained. These samples

must be processed exactly the same way as the assay samples, and it is recommended to analyze their results using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PROCEDURE**Preparation of reagents**

Let all the reagents come to room temperature.

Reconstitution of calibrators and control serum

The content of the vials is reconstituted with the volume of distilled water indicated on the label. Wait for 30 min following reconstitution and mix gently to avoid foaming before dispensing. Store the reconstituted calibrators aliquoted at < -18°C until the expiry date of the kit. Store the reconstituted control samples at 2-8°C for one week or aliquoted at < -18°C for a longer time, until the expiry date of the kit.

Preparation of the wash solution

Pour the content of the vial into 950 mL of distilled water and homogenize. The diluted solution may be stored at 2-8°C until the expiry date of the kit.

Procedure**Bring all reagents to room temperature before pipeting.**

Step 1 Additions*	Step 2 Incubation	Step 3 Counting
To antibody coated tubes, add successively: 100 µL of calibrator, control or sample and 50 µL of tracer. Mix.	Incubate 90 minutes at 18-25°C with shaking (>350 rpm).	Aspirate carefully the content of tubes (except the 2 tubes «total cpm»). Wash twice with 2 mL of wash solution. Count bound cpm (B) and total cpm (T) for 1 minute.

*Add 50 µL of tracer to 2 additional tubes to obtain total cpm.

PERFORMANCE CHARACTERISTICS*(For more details, see the data sheet "APPENDIX")*

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity**Analytical sensitivity:** 0.16 IU/L**Functional sensitivity:** 0.48 IU/L**Specificity**

The antibody used in the immunoassay is highly specific for LH. Extremely low cross reactivities were obtained against several related molecules (hTSH, hFSH, hCG).

Precision**Intra-assay**

Serum samples were assayed 25 times in the same series. The coefficients of variation were found below or equal to 7.33%.

Inter-assay

Serum samples were assayed in duplicate in 10 different series. The coefficients of variation were found below or equal to 8.42%.

Accuracy**Dilution test**

High-concentration serum samples were serially diluted with the zero calibrator. The recovery percentages obtained were between 80.7% and 102%.

Recovery test

Low-concentration serum samples were spiked with known quantities of LH. The recovery percentages obtained were between 103% and 119%.

Measurement range (from analytical sensitivity to highest calibrator):

0.16 to approximately 180 IU/L.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly.

Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Do not use hemolyzed, lipemic or icteric samples.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been

regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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