OSTASE IRMA KIT

REF A44176

TABLE OF CONTENTS

English	2	Čeština	23
Français	5	Slovenčina	26
Deutsch	8	Türkçe	29
Italiano	11	Русский	32
Español	14	APPENDIX	35
Português (Portugal)	17		
Ελληνικά	20		

OSTASE IRMA KIT

REF A44176

IMMUNORADIOMETRIC ASSAY FOR THE IN VITRO DETERMINATION OF BONE ALKALINE PHOSPHATASE (BAP) IN HUMAN SERUM For *in vitro* diagnostic use.

PRINCIPLE

The immunoradiometric assay of bone alkaline phosphatase (BAP) is a "sandwich" type assay. In the kit, mouse monoclonal antibodies directed against two different epitopes of BAP and hence not competing are used.

Serum samples, controls and calibrators are incubated in tubes coated with the first monoclonal antibody in the presence of the second monoclonal antibody which is labeled with iodine 125. After incubation, the content of tubes is rinsed so as to remove unbound 125I-labeled antibody. The bound radioactivity is then determined in a gamma counter. The BAP concentrations in the samples are obtained by interpolation from the standard curve. The concentration of BAP in the samples is directly proportional to the radioactivity.

WARNING AND PRECAUTIONS

General remarks:

- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be included with each assay.
- The size of the assay should be limited so as all the samples can be pipetted within 20 minutes.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use.

Adherence to the basic rules of radiation safety should provide adequate protection:

- No eating, drinking, smoking or application of cosmetics should be carried out in the presence of radioactive materials.
- No pipetting of radioactive solutions by mouth.
- Avoid all contact with radioactive materials by using gloves and laboratory overalls.
- All manipulation of radioactive substances should be done in an appropriate place, distant from corridors and other busy places.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as a preservative. Sodium azide can react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Materials of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HbsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All serum samples should be handled as if capable of transmitting hepatitis or AIDS and waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Wash solution (20x) DANGER

H360	May damage fertility or the unborn child.
P201	Obtain special instructions before use.
P280	Wear protective gloves, protective clothing and eye/face protection.
P308+P313	IF exposed or concerned: Get medical advice/attention. Boric Acid 0.1-0.3% Sodium Borate Decahydrate 0.1-0.3%

Safety Data Sheet is available at techdocs.beckmancoulter.com

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in dry tubes.
- Separate serum from cells by centrifugation.
- Serum samples may be stored at 2-8 °C, if the assay is to be performed within 24 hours. For longer storage keep frozen (<-20 °C, maximum 2 months) after aliquoting so as to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- If samples have concentrations greater than the highest calibrator, they must be diluted in the zero calibrator.

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8 °C. Expiry dates printed on vial labels apply to the long-term storage of components by the manufacturer only, prior to assembly of the kit. Do not take into account.

Storage conditions for reagents after reconstitution or dilution are indicated in paragraph Procedure.

Anti- BAP monoclonal antibody coated tubes: 2 x 50 tubes (ready-to-use)

¹²⁵I-labeled monoclonal anti- BAP antibody: one 11 mL vial (ready-to-use)

The vial contains 650 kBq, at the date of manufacture, of 125I-labeled immunoglobulins in buffer with bovine and horse proteins and sodium azide (<0.1 %) and a dye.

Note: Occasional presence of clotted particles in the tracer does not affect assay performance.

Calibrators: five 1 mL vials and one 6 mL "zero calibrator" vial (ready-to-use)

The calibrator vials contain from 0 to approximately 120 µg/L of BAP in buffer with bovine serum albumin and sodium azide (<0.1 %). The exact concentration is indicated on each vial label. The calibrators are verified to an internal reference standard.

Control sera: two vials (ready-to-use)

The vials contain BAP in buffer with bovine serum albumin and sodium azide (<0.1 %). The expected values are in the concentration range indicated on a supplement.

Wash solution (20x): one 50 mL vial

Concentrated solution has to be diluted before use.

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- Precision micropipette (100 µL).
- Repeating pipettes (100 µL and 2 mL).
- Refrigerator or cold room (2-8 °C).
- Vortex-type mixer.
- Aspiration system.

- Gamma counter set for ¹²⁵I.

RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of BAP concentrations in samples assayed at the same time as the calibrators.

Standard curve

The results in the package insert were calculated using a log-log curve fit (weighted quadratic regression) with determined radioactivity ($\text{cpm}_{\text{cal}} - \text{cpm}_{\text{cal0}}$) on vertical axis and the BAP concentration of the calibrators on the horizontal axis ($\mu\text{g/L}$). Other data reduction methods may give slightly different results.

Total activity: 225,220 cpm				
Calibrators	BAP ($\mu\text{g/L}$)	cpm (n=3)	B/T (%)	$\text{cpm}_{\text{cal}} - \text{cpm}_{\text{cal0}}$
0	0	196	0.09	-
1	7.5	2,481	1.10	2,285
2	15	5,122	2.27	4,926
3	30	9,448	4.20	9,252
4	60	18,208	8.08	18,012
5	120	32,125	14.3	31,929

(Example of standard curve, do not use for calculation)

Samples

For each sample or control, locate the ($\text{cpm}_{\text{sample}} - \text{cpm}_{\text{cal0}}$) value on the vertical axis and read off the corresponding BAP concentration on the horizontal axis. The concentrations of diluted samples must be corrected by the dilution factor.

Enzymatic activity 100 U/L of BAP corresponds approximately to 38.4 $\mu\text{g/L}$.

EXPECTED VALUES

We recommend each laboratory to establish its own reference values. The below presented values obtained with healthy subjects are indicative only.

A multicenter clinical trial was conducted to test the effectiveness of BAP as an aid in the management of postmenopausal osteoporosis and Paget's disease. Although the BAP results in this trial were generated with the Tandem®-R OSTASE assay, the Immunotech OSTASE IRMA kit assay has been developed using the same monoclonal antibodies employed in the Tandem-R OSTASE assay and has been standardized to provide the same clinical performance.

Assays measuring BAP using other monoclonal antibodies or other methods may not provide similar clinical performance.

Study results are provided in the following table.

	n	Mean ($\mu\text{g/L}$)	SD ($\mu\text{g/L}$)	Median ($\mu\text{g/L}$)	95 th Percentile ($\mu\text{g/L}$)
Males	217	12.3	4.3	11.6	20.1
Premenopausal Females	228	8.7	2.9	8.5	14.3
Postmenopausal Females	529	13.2	4.7	12.5	22.4

QUALITY CONTROL

Good laboratory practices imply that control samples be used regularly to ensure the quality of the results obtained. These samples must be processed exactly in the same way as the assay samples, and it is recommended that their results be analyzed using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PROCEDURE

Preparation and storage of reagents

Let all the reagents come to room temperature.

Preparation of the wash solution

Pour the content of the vial into 950 mL of distilled water and homogenize. The diluted solution may be stored at 2-8 °C until the expiry date of the kit.

Assay procedure (see table below)

Bring all reagents to room temperature before pipetting.

Step 1 Additions	Step 2 Incubation	Step 3 Washing and Counting
To antibody coated tubes, add successively: 100 μL of calibrator, control or sample and 100 μL of tracer.* Shake gently test tube rack by hand for 15 seconds.	Put the rack into pre-cooled water bath placed in refrigerator. Incubate 19 ± 1 hours at 2-8 °C without shaking. The temperature must not exceed 8 °C.	Add to all tubes 2 mL of wash solution and aspirate** (with exception of the 2 tubes «total cpm»). Add once more 2 mL of wash solution and aspirate.** Determine activity (cpm) for 1 min.

*Add 100 μL of tracer to 2 additional tubes to obtain total cpm.

**Aspirate each rack of tubes separately to decrease the time of washing solution presence in the tube to minimum. After finalization of second washing step, wait approximately 60 seconds and aspirate once more.

PERFORMANCE CHARACTERISTICS

(for more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Analytical sensitivity: 1.75 $\mu\text{g/L}$

Functional sensitivity: 5.00 $\mu\text{g/L}$

Specificity

Intestinal BAP: 100 U/L gives a result of 2.4 $\mu\text{g/L}$ in this assay.

Placental BAP: 100 U/L did not produce detectable result in this assay.

Liver BAP: 100 U/L gives a result of 6.9 $\mu\text{g/L}$ in this assay.

Precision

Intra-assay

Samples were assayed in 18 replicates in the same series. The coefficients of variation were found below or equal to 8.4 %.

Inter-assay

Samples were assayed in duplicate in 10 different series. Coefficients of variation were found below or equal to 13.6 %.

Accuracy

Dilution test

High-concentration serum samples were serially diluted in zero calibrator. The recovery percentages were obtained between 88.5 % and 115 %.

Recovery test

Serum samples were spiked with known quantities of BAP. The recovery percentages were obtained between 86.7 % and 106 %.

Measurement range (from analytical sensitivity to highest calibrator): 1.75 to approximately 120 $\mu\text{g/L}$.

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly.

Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

It is not recommended to use haemolysed, icteric or lipaemic samples with the OSTASE assay.

Hook effect: "Hook effect" of this IRMA kit may appear in samples containing BAP concentration equal or higher than 1000 $\mu\text{g/L}$. The samples with clinical suspicion to very high concentrations should be diluted with zero calibrator. Values obtained must be corrected by the dilution factor.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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Accuracy

Dilution test

Serum samples were diluted with zero calibrator and assayed according to the assay procedure of the kit.

Sample	Dilution	BAP (µg/L)		Recovery measured/expected (%)
		expected	measured	
K	Undiluted	-	103	-
	1/2	51.4	47.1	91.7
	1/4	25.7	23.5	91.5
	1/8	12.8	11.7	91.1
	1/16	6.42	6.24	97.2
L	Undiluted	-	82.7	-
	1/2	41.4	36.6	88.5
	1/4	20.7	18.9	91.4
	1/8	10.3	9.75	94.3
	1/16	5.17	5.90	114
M	Undiluted	-	56.3	-
	1/2	28.2	27.4	97.3
	1/4	14.1	13.8	98.0
N	Undiluted	-	50.2	-
	1/2	25.1	24.9	99.2
	1/4	12.6	13.3	106
	1/8	6.28	7.20	115

Recovery test

Serum samples containing known concentrations of BAP were added to three sera and assayed according to the assay procedure of the kit.

Endogenous conc.	BAP (µg/L)		Recovery (%) Measured/Expected
	Added	Obtained	
59.4	40.2	101.9	106.0
	58.7	119.8	102.9
	29.7	87.8	95.7
26.6	40.2	68.9	105.2
	20.0	46.1	97.7
	10.5	35.7	86.7
19.0	58.7	76.1	97.3
	29.7	48.4	99.1
	15.2	33.0	92.0

Precision

Intra-assay

Sample	A	B	C	D	E	F
No. of determinations	19	18	18	18	18	19
Mean value (µg/L)	12.9	17.8	23.2	40.6	53.4	99.9
CV (%)	8.39	4.61	3.55	6.18	2.93	2.83

Inter-assay

Sample	1	2	3	4	5	6
No. of determinations	10	10	10	10	10	10
Mean value (µg/L)	6.77	8.40	11.8	25.0	45.5	74.1
CV (%)	10.9	13.6	9.35	6.53	7.06	9.30

Interference by drugs

Analogous IRMA sandwich system Tandem-R OSTASE, employing the same monoclonal antibodies, was tested for interference by drugs. The drugs and the highest concentrations tested were:

Acetaminophen (350 µg/mL), alendronate (20 µg/mL), aspirin (350 µg/mL), calcitonin-human (80 µg/mL), calcitonin-salmon (60 IU/mL), calcium (500 µg/mL), estrogen (100 µg/mL), etidronate (350 µg/mL), ibuprofen (150 µg/mL), norethidronel/ethinyl estradiol mixture (oral contraceptive) (3.0 mg/mL), pamidronate (180 µg/mL), vitamin D (400 IU/mL).

These drugs did not interfere in the Tandem-R OSTASE assay. Due to analogous design of the Immunotech OSTASE IRMA assay, no interference is expected as well.

¹²⁵I Characteristics

T_{1/2} (¹²⁵I) = 1443 h = 60.14 d

¹²⁵ I	E (MeV)	%
Gamma	0.035	
X	0.027	114
	0.032	25

Symbols Key

REF Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferens / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Úrün Referansı / Ссылка на продукт / Референция за продукт / 產品參考

IVD In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / Ин витро диагностика / 體外診斷

CONTENTS Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 組成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄

 Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 製造商 / Gamintojas / Gyártó: / Producent / Výrobce / Výrobca / 제조 / Üretici / Изготовлено / Произведено от / 製造商

 Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Räcker till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 內容物足夠執行 <n> 次測試

CE CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Mercado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識

SDS Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表

 Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 請參閱使用說明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문의 / Kullanma Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明

 Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度範圍 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sıcaklık aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍

 Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Upozornění / Upozornenie / 주의 / Dikkat / Внимание / 注意

 Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirație / 만료 날짜 / Son Kullanma Tarihi / Срок годности / Срок на годност / 到期日

LOT Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer serii / Číslo sarže / 로트 번호 / Lot Numarası / Номер партии / Номер на партида / 批號

 Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Datum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biyolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktiv / Ραδιενεργό / 放射性 / Radioaktyvioji medžiaga / Radioaktiv / Radioaktywny / Radioaktivní / Rádioaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

DANGER

DANGER / DANGER / GEFAHR / PERICOLO / PELIGRO / PERIGO / FARA / ΚΙΝΔΥΝΟΣ / 危險 / PAVOJUS / VESZÉLY / NIEBEZPIECZENSTWO / NEBEZPEČÍ / NEBEZPEČENSTVO / 위험 / TEHLÍKE / ОПАСНО / ОПАСНОСТ / 危險

Ag^{125I}

Ab^{125I}

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / Tracer / Ανιχνευτής / 追踪剂 / Atsekamoji medžiaga / Nyomjelző / Znacznik / Radioindikator / Indikator (tracer) / 트레이서 / Tracer'lar / метка / Индикатор / 追蹤劑

CAL

CAL 0

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrator / Kalibrátor / 보정 물질 / Kalibratör / Калибратор / Калибратор / 校正液

CTRL

Control / Contrôle / Kontrolle / Controllo / Control / Controllo / Kontrolle / Μάρτυρας / 质控品 / Kontrolliné / Kontroll / Kontrola / Kontrola / Kontrola / 정도관리 / Kontrol / Контроль / Контролна / 質控品

TUBE

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνάκια / 试管 / Mégintüveliai / Csövek / Probówki / Zkumavky / Skúmavky / 튜브 / Tüpler / пробирки / Bulgarian / 試管

IFU

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用說明

SOLN | **WASH** | **20x**

Wash Solution Concentrate 20X / Solution de lavage concentrée 20X / Waschlösungskonzentrat 20X / Concentrato di soluzione di lavaggio 20X / Solución de lavado concentrada 20X / Concentrado de solução de lavagem 20X / Tvättlösningkoncentrat 20X / Συμπυκνωμένο διάλυμα πλύσης 20X / 浓缩清洗液 20X / Plovimo tirpalu koncentratas 20X / 20X mosóoldat-koncentrátum / Koncentrat 20X roztworu płuczającego / Koncentrát mycího roztoku 20X / Концентрат промывочного раствора 20X / Концентрат на разтвор за промиване 20X / 清洗溶液濃縮 20X

Beckman Coulter do Brasil Com. e Imp. de Prod. de Lab. Ltda,
Calçada Aldebarã, 39, Centro de Apoio 2 - Alphaville,
CEP 06541-055 - Santana de Parnaíba, SP, Brasil
Telefone: (11) 4154-8818 - CNPJ: 42.160.812/0001-44

IMMUNOTECH s.r.o., Radiova 1, 102 27 Prague 10, Czech Republic