RIA PROGESTERONE

REF IM1188**TABLE OF CONTENTS**

English	2	Čeština (CZ)	23
Français (FR)	4	Slovenčina (SK)	25
Deutsch (DE)	6	한국어 (KO)	27
Italiano (IT)	8	Türkçe (TR)	29
Español (ES)	10	Русский (RU)	31
Português Portugal (PT-PT)	12	Србија (SR)	34
Ελληνικά (EL)	14	Português Brasil (PT-BR)	36
Lietuviškai (LT)	17	APPENDIX	38
Magyar (HU)	19	REFERENCES	42
Polski (PL)	21		

Beckman Coulter do Brasil Com. e Imp. de Prod. de Lab. Ltda
Alameda Rio Negro, 500, 15º andar, Torre B Alphaville Industrial
CEP 06.454-00 Barueri, São Paulo, Brasil
CNPJ: 42.160.812/0001-44 Telephone: 0800-771-8818

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RIA PROGESTERONE

REF IM1188

RADIOIMMUNOASSAY FOR THE IN VITRO DETERMINATION OF PROGESTERONE IN HUMAN SERUM AND PLASMA

For *in vitro* diagnostic use.

PRINCIPLE

The radioimmunoassay of progesterone is a competition assay. Samples and calibrators are incubated with ¹²⁵I-labeled progesterone, as tracer, in antibody-coated tubes. After incubation the content of tubes is aspirated and bound radioactivity is measured. A calibration curve is established and unknown values are determined by interpolation from the curve.

WARNING AND PRECAUTIONS

General remarks:

- The vials with calibrators and controls should be opened as shortly as possible to avoid excessive evaporation.
- Do not mix the reagents from kits of different lots.
- A standard curve must be included with each assay.
- The correct setting of the shaker is very important for the reproducibility of the assay.
- It is recommended to perform the assay in duplicate.
- Each tube must be used only once.

Basic rules of radiation safety

The purchase, possession, utilization, and transfer of radioactive material is subject to the regulations of the country of use. Adherence to the basic rules of radiation safety should provide adequate protection.

- No eating, drinking, smoking or application of cosmetics should take place in the presence of radioactive materials.
- No pipeting by mouth.
- Avoid all contact with radioactive materials by using gloves and lab coat.
- All manipulation of radioactive substances should be done in an appropriate location, away from corridors and other busy areas.
- Radioactive materials should be stored in the container provided and in a designated area.
- A record of receipt and storage of all radioactive products should be kept up to date.
- Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes.
- Each case of radioactive contamination or loss of radioactive material should be resolved according to established procedures.
- Radioactive waste should be handled according to the rules established in the country of use.

Sodium azide

Some reagents contain sodium azide as preservative. Sodium azide may react with lead, copper or brass to form explosive metal azides. Dispose of the reagents by flushing with large amounts of water through the plumbing system.

Materials of human origin

The materials of human origin, contained in this kit, were found negative for the presence of antibodies to HIV 1 and HIV 2, antibodies to HCV, as well as of Hepatitis B surface antigen (HBsAg). However, they should be handled as if capable of transmitting disease. No known test method can offer total assurance that no virus is present. Handle this kit with all necessary precautions.

All serum and plasma samples should be handled as if capable of transmitting hepatitis or AIDS. Waste should be discarded according to the country rules.

GHS HAZARD CLASSIFICATION

Not classified as hazardous

SDS

Safety Data Sheet is available at beckmancoulter.com/techdocs

SPECIMEN COLLECTION, PROCESSING, STORAGE AND DILUTION

- Collect blood in dry tubes or tubes containing EDTA.
- Separate serum or plasma from cells by centrifugation.
- Serum and plasma samples may be stored at 2-8°C, if the assay is to be performed within 24 hours. For longer storage keep frozen (< -20°C, for up to 1 year) after aliquoting so as to avoid repeated freezing and thawing. Thawing of sample should be performed at room temperature.
- If samples have concentrations greater than the highest calibrator, they must be diluted into a serum with a low concentration of progesterone (e.g. male serum with concentration ≤ 0.1 ng/mL).

Serum and EDTA plasma values for 30 samples (serum samples ranging from 0.09 to 3.93 ng/mL) were compared using the IM1188 Progesterone RIA kit. Results are as follows:

$$[\text{EDTA-plasma}] = 0.8823[\text{serum}] + 0.0379$$

$$R = 0.9973$$

MATERIALS PROVIDED

All reagents of the kit are stable until the expiry date indicated on the kit label, if stored at 2-8°C. Expiry dates printed on component vial labels apply to the long-term storage by manufacturer only, prior to assembling of the kit. Do not take into account.

Anti-progesterone antibody-coated tubes: 2 x 50 tubes (ready-to-use)

¹²⁵I-labeled progesterone: one 55 mL vial (ready-to-use)

The vial contains 200 kBq, at the date of manufacture, of ¹²⁵I-labeled progesterone in buffer with proteins and sodium azide (<0.1%).

Note: Occasional presence of clotted particles in the tracer does not affect assay performance.

Calibrators: six 0.5 mL vials (ready-to-use)

The calibrator vials contain concentration of progesterone for a standard range from 0 to approximately 50 ng/mL in human serum with sodium azide (<0.1%). The exact concentration is indicated on each vial. The calibrators were calibrated against the reference preparation ERM®-DA347.

Control samples: two 0.5 mL vials (ready-to-use)

The vials contain progesterone in human serum with sodium azide (<0.1%). The expected values are in the concentration range indicated on a supplement.

MATERIALS REQUIRED, BUT NOT PROVIDED

In addition to standard laboratory equipment, the following items are required:

- precision micropipet (50 µL).
- semi-automatic pipet (500 µL).
- vortex type mixer.
- horizontal or orbital shaker.
- aspiration system.
- gamma counter set for 125 iodine.

PROCEDURE

Bring all reagents to room temperature before pipeting.

Step 1 Additions	Step 2 Incubation	Step 3 Counting
To coated tubes add sequentially: 50 µL of calibrator, control or sample and 500 µL of tracer* Mix	Incubate 1 hour at 18-25°C with shaking (≥280 rpm).	Aspirate carefully the contents of tubes (except the 2 tubes «total cpm») Count bound cpm (B) and total cpm (T) for 1 min.

*Add 500 µL of tracer to 2 additional tubes to obtain total cpm.

RESULTS

Results are obtained from the standard curve by interpolation. The curve serves for the determination of progesterone concentrations in samples measured at the same time as the calibrator.

Standard curve

The results in the quality control department were calculated using *spline* curve fit with B/T or B/B_0 on the logit vertical axis and analyte concentration of the calibrators on the log horizontal axis (ng/mL).

Other data reduction methods may give slightly different results.

Total activity: 73,748 cpm				
Calibrators	Progesterone (ng/mL)	cpm (n=3)	B/T (%)	B/B ₀ (%)
0	0	27,619	37.45	100
1	0.09	22,020	29.86	79.73
2	0.43	17,861	24.22	64.67
3	2.00	9,958	13.50	36.05
4	12.5	4,279	5.80	15.49
5	65.0	1,676	2.27	6.07

(Example of standard curve, do not use for calculation)

Samples

For each sample, locate the B/T or B/B_0 on the vertical axis and read off the corresponding progesterone concentration on the horizontal axis.

To convert ng/mL into nmol/L (nM), multiply results by 3.18.

EXPECTED VALUES

It is suggested that each laboratory establishes its own normal values. The following values obtained with healthy subjects are indicative only.

Age	N	Min.	Max.	Median	2.5 th percentile	97.5 th percentile
ng/mL						
Males (20 to 67 years)	98	0.07	0.86	0.32	0.11	0.77
Females						
Follicular phase	339	<0.03	1.63	0.31	0.04	0.93
Preovulatory peak	69	0.06	1.79	0.76	0.26	1.37
Luteal phase	306	0.50	32.2	6.79	0.94	22.5
Postmenopausal	48	0.08	0.62	0.22	0.09	0.59

(For more details, see APPENDIX)

QUALITY CONTROL

Good laboratory practices imply that control samples be used regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended that their results be analyzed using appropriate statistical methods.

In case of packaging deterioration or if data obtained show some performance alteration, please contact your local distributor or use the following e-mail address: imunochem@beckman.com

PERFORMANCE CHARACTERISTICS

(for more details, see the data sheet "APPENDIX")

Representative data are provided for illustration only. Performance obtained in individual laboratories may vary.

Sensitivity

Limit of detection (LoD): 0.04 ng/mL

The LoD for Progesterone is 0.04 ng/mL, determined consistent with guidelines in CLSI document EP17-A2 [1] based on the proportions of false positives (α) less than 5% and false negatives (β) less than 5%; using determinations, with 120 blank and 120 low level samples; and Limit of Blank (LoB) of 0.02 ng/mL.

Specificity

The antibody used in the immunoassay is highly specific for progesterone.

Precision

Repeatability and within-laboratory precision

The precision of the assay was determined consistent with guidelines in CLSI document EP05-A3 [2]. For repeatability the coefficients of variation were found below or equal to 9.48 % for serum samples. For within-laboratory precision the coefficients of variation were found below or equal to 16.85 % for serum samples.

Accuracy

Linearity

The assay demonstrated to be linear from 0.02 ng/mL to 76.02 ng/mL using serum samples (determined consistent with guidelines in CLSI document EP06-A [3]).

Dilution test

High-concentration serum samples were serially diluted in a sample with a low concentration of progesterone (e.g. male serum with concentration ≤ 0.1 ng/mL). The recovery percentages obtained were between 82.2% and 110%.

Recovery test

Low-concentration serum samples were spiked with known quantities of progesterone. The recovery percentages obtained were between 80.5% and 98.8%.

Measurement range (from LoD to highest calibrator): from 0.04 to approximately 50 ng/mL

LIMITATIONS

The non-respect of the instructions in this package insert may affect results significantly. Results should be interpreted in the light of the total clinical presentation of the patient, including clinical history, data from additional tests and other appropriate information.

Do not use hemolyzed, lipemic or icteric samples. For more details, see Appendix, § Interference.

For assays employing antibodies, the possibility exists for interference by heterophile antibodies in the patient sample. Patients who have been regularly exposed to animals or have received immunotherapy or diagnostic procedures utilizing immunoglobulins or immunoglobulin fragments may produce antibodies, e.g. HAMA, that interfere with immunoassays.

Such interfering antibodies may cause erroneous results. Carefully evaluate the results of patients suspected of having these antibodies.

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Interference

Serum samples containing Progesterone concentrations (low and high) were spiked with multiple concentrations of the substances below and assayed using Progesterone RIA kit. Values were calculated as described in CLSI EP07, 3rd ed. [4]. Interference was determined by testing controls (no interfering substance added) and matched test samples (with interfering substance added). No interference (defined as a shift in dose >15%) was found for addition of interferent up to concentration stated in the table below.

In spite of hemoglobin, bilirubin (conjugated, unconjugated) and triglyceride interference data in the table, we advise to avoid using hemolyzed, lipemic or icteric samples.

Interferent	Test concentration
Acetylsalicylic acid	39.21 µg/mL
Ascorbic acid	63.39 µg/mL
Biotin	1,590 ng/mL
Conjugated bilirubin	458.8 µg/mL
Hemoglobin	10,391 µg/mL
Heparin	7,790 ng/mL
Cholesterol	2.34 mg/mL
Ibuprofen	158.8 µg/mL
Prednisone	128.3 ng/mL
Prednisolone	784.0 ng/mL
Rheumatoid factor	34.47 IU/mL
Triglycerides	17.77 mg/mL
Unconjugated bilirubin	453.1 µg/mL

Specificity:

Data on cross-reactivity with several endogenous and pharmaceutically employed steroids are presented in the following table:

$$\text{Cross-reactivity (\%)} = \frac{\text{Progesterone concentration}}{\text{Steroid concentration}} \times 100 \text{ at 50\% binding of the zero calibrator}$$

Steroid	Cross-reactivity (%)
Progesterone	100
5α-Pregnanedione	14.11
5β-Pregnanedione	6.26
6β-Hydroxyprogesterone	5.16
Corticosterone	3.87
11-Deoxycorticosterone	1.56
20α-Dihydroprogesterone	0.84
16α-Hydroxyprogesterone	0.83
Pregnenolone	0.75
Pregnanolone	0.39
17-α-hydroxyprogesterone	0.66
Testosterone	0.24
Pregnanolone sulfate	0.24
Pregnenolone sulfate	0.17
Δ4-Androstenedione	0.05
Medroxyprogesterone acetate	0.03
19-Nortestosterone	0.02
Cortisol	0.01
11-Deoxycortisol	ND
19-Norethisterone	ND
Androstenediol	ND
Danazol	ND
DHEA sulfate	ND
Prednisolone	ND
Estradiol	ND
Spirolactone	ND

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Six serum samples and six EDTA plasma samples were assayed for 20 days, 2 runs per day, triplicates per run. Assays were performed by two lab technicians, by two reagent lots. There were 120 individual measurements per sample with no invalid results.

Serum	Mean (ng/mL)	Repeatability		Within-laboratory precision	
		SD (ng/mL)	C.V. (%)	SD (ng/mL)	C.V. (%)
S1	24.20	1.97	8.12	3.32	13.70
S2	11.55	0.71	6.14	1.34	11.64
S3	5.19	0.38	7.23	0.56	10.70
S4	1.52	0.09	6.04	0.16	10.27
S5	0.77	0.07	8.61	0.10	12.47
S6	0.32	0.03	9.48	0.05	16.85

EDTA-plasma	Mean (ng/mL)	Repeatability		Within-laboratory precision	
		SD (ng/mL)	C.V. (%)	SD (ng/mL)	C.V. (%)
P1	35.74	2.94	8.23	4.85	13.57
P2	11.95	0.87	7.27	1.30	10.85
P3	6.03	0.39	6.40	0.67	11.12
P4	2.26	0.12	5.14	0.23	10.08
P5	1.56	0.10	6.62	0.17	10.76
P6	0.42	0.03	7.74	0.06	15.00

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The assay demonstrated to be linear from 0.01 to 86.83 ng/mL using EDTA plasma samples (determined consistent with guidelines in CLSI document EP06-A [3]).

Dilution test

Three serum/EDTA plasma samples were serially diluted by human serum with low concentration of progesterone and assayed according to the procedure of the kit.

Serum	Dilution factor	Progesterone (ng/mL)		Ratio (%) Measured/Expected
		Measured	Expected	
S1	-	47.29	-	-
	1:2	25.96	23.65	109.8
	1:4	11.18	11.82	94.57
	1:8	5.47	5.91	92.54
	1:16	2.91	2.96	98.46
	1:32	1.49	1.48	100.8
S2	-	50.34	-	-
	1:2	25.65	25.17	101.9
	1:4	11.76	12.59	93.44
	1:8	5.17	6.29	82.16
	1:16	2.92	3.15	92.81
	1:32	1.52	1.57	96.62
S3	-	47.43	-	-
	1:2	25.83	23.72	108.9
	1:4	11.26	11.86	94.96
	1:8	5.28	5.93	89.06
	1:16	2.76	2.96	93.11
	1:32	1.54	1.48	103.9

EDTA-plasma	Dilution factor	Progesterone (ng/mL)		Ratio (%) Measured/Expected
		Measured	Expected	
P1	-	24.84	-	-
	1:2	14.04	12.42	113.0
	1:4	6.48	6.21	104.3
	1:8	2.87	3.11	92.43
	1:16	1.38	1.55	88.89
	1:32	0.63	0.78	81.16
P2	-	31.95	-	-
	1:2	17.29	15.98	108.2
	1:4	8.28	7.99	103.7
	1:8	3.29	3.99	82.38
	1:16	1.71	2.00	85.63
	1:32	0.80	1.00	80.13
P3	-	39.10	-	-
	1:2	21.34	19.55	109.2
	1:4	10.02	9.78	102.5
	1:8	4.33	4.89	88.59
	1:16	2.00	2.44	81.84
	1:32	0.99	1.22	81.02

Recovery test

Three serum/EDTA plasma samples were spiked by known amount of progesterone and samples were then assayed according to the procedure of the kit.

Serum	Endogen. conc. (ng/mL)	Added conc. (ng/mL)	Expected conc. (ng/mL)	Measured conc. (ng/mL)	Ratio (%) Measured/Expected
S1	0.27	0.18	0.45	0.37	81.40
	0.29	0.29	0.58	0.48	83.02
	0.28	0.71	0.99	0.81	81.77
S2	0.29	0.18	0.47	0.41	86.73
	0.31	0.40	0.71	0.58	81.36
	0.30	0.71	1.01	0.87	86.19
S3	1.70	0.78	2.48	2.12	85.55
	1.67	0.96	2.63	2.12	80.54
	1.72	3.10	4.82	4.76	98.77

EDTA plasma	Endogen. conc. (ng/mL)	Added conc. (ng/mL)	Expected conc. (ng/mL)	Measured conc. (ng/mL)	Ratio (%) Measured/Expected
P1	0.28	0.18	0.46	0.43	92.75
	0.30	0.29	0.59	0.52	88.44
	0.29	0.71	1.00	0.89	89.00
P2	0.54	0.18	0.72	0.60	83.54
	0.57	0.40	0.97	0.82	84.17
	0.56	0.71	1.26	1.06	83.85
P3	1.77	0.78	2.55	2.19	85.78
	1.74	0.96	2.71	2.38	87.95
	1.80	3.10	4.90	4.43	90.50

Expected values

Values of progesterone in fertile females

Phase	N	Min.	Max.	Median	2.5 th	97.5 th
					percentile	percentile
Follicular phase						
- Early	103	<0.03	1.63	0.44	0.05	1.08
- Middle	134	<0.03	1.02	0.28	<0.03	0.86
- Late	102	0.05	0.91	0.30	0.05	0.88
Preovulatory peak	69	0.06	1.79	0.76	0.26	1.37
Luteal phase						
- Early	108	0.50	17.2	2.57	0.79	10.5
- Middle	161	1.71	32.2	10.5	2.15	25.3
- Late	37	0.50	15.2	7.24	0.64	14.6

¹²⁵I Characteristics

T_{1/2} (¹²⁵I) = 1443 h = 60.14 d

¹²⁵ I	E (MeV)	%
γ	0.035	
X	0.027	114
	0.032	25

Symbols Key

- REF** Product Reference / Référence du produit / Produktreferenz / Riferimento prodotto / Número de referencia del producto / Referência do produto / Produktreferens / Κωδικός αναφοράς προϊόντος / 产品参考 / Gaminio nuoroda / Termékszám / Dane referencyjne produktu / Reference k produktu / Referenčné označenie výrobku / 제품 참조 자료 / Ürün Referansı / Ссылка на продукт / Референция за продукт / 產品參考
- IVD** In Vitro Diagnostic / Diagnostic in vitro / In-vitro-Diagnostikum / Diagnostica in vitro / Para diagnóstico in vitro / Diagnóstico in vitro / InVitro-diagnostik / Για διάγνωση in vitro / 体外诊断 / In vitro diagnostika / In vitro diagnosztikai felhasználásra / Diagnostyka in vitro / Diagnostika in vitro / 체외 진단 / In Vitro Diagnostik / Diagnostika in vitro / За ин витро диагностика / 體外診斷
- CONTENTS** Contents / Contenu / Inhalt / Contenuto / Contenido / Conteúdo / Περιεχόμενο / 组成 / Rinkinio sudėtis / Tartalom / Zawartość / Obsah / Obsah / 내용물 / İçindekiler / Содержание / Съдържание / 目錄
-  Manufactured by / Fabriqué par / Hergestellt von / Prodotto da / Fabricado por / Tillverkas av / Κατασκευαστής / 制造商 / Gamintojas / Gyártó: / Producent / Výrobce / Výrobca / 제조 / Üretici / Изготовлено / Произведено от / 製造商
-  Contains sufficient for <n> tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos / Conteúdo suficiente para "n" ensaios / Räcker till "n" antal tester / Περιεχόμενο επαρκές για "n" εξετάσεις / 含量足够 <n> 次测试 / Turinio pakanka <n> tyrim / <n> teszthez elegendő mennyiséget tartalmaz / Zawartość wystarcza na <n> testów / Lze použít pro <n> testů / Obsah vystačí na <n> testov / <n> 테스트에 대해 충분한 양 포함 / <n> sayida test için yeterlidir / Содержит достаточно для количества тестов: <n> / Съдържа достатъчно за <n> теста / 内容物足夠執行 <n> 次測試
- CE** CE Mark / Marquage CE / CE-Kennzeichnung / Marchio CE / Marcado CE / Marcação CE / CE-märkning / Σήμανση CE / CE 标志 / CE ženklas / CE jelzés / Znak CE / Značka CE / Označenie CE / CE 표시 / CE İşareti / Маркировка CE / CE маркировка / CE 標識
- SDS** Safety Data Sheet / Fiche technique santé-sécurité / Sicherheitsdatenblatt / Scheda dati di sicurezza / Hoja de datos de seguridad / Ficha de Dados de Segurança / Säkerhetsdatablad / Φύλλο Δεδομένων Ασφάλειας / 安全数据单 / Saugos duomenų lapas / Biztonsági adatlap / Karta Charakterystyki Bezpieczeństwa / Bezpečnostní list / Bezpečnostný list / 안전보건자료 / Güvenlik Bilgi Formu / Паспорт безопасности / Информационен Лист За Безопасност / 安全性資料表
-  Consult Instructions for Use / Consultez le mode d'emploi / Siehe Gebrauchsanweisung / Consultare le istruzioni per l'uso / Consulte las Instrucciones de uso / Instruções de utilização / Konsultera bruksanvisning / Συμβουλευτείτε τις οδηγίες χρήσης / 请参阅使用说明 / Skaitykite naudojimo instrukciją / Olvassa el a használati utasítást / Zapoznać się z instrukcją użycia / Postupujte podle návodu k použití / Prečítajte si návod na použitie / 사용 안내 문외 / Kullanna Talimatına Başvurun / Обратитесь к инструкциям / Вижте Инструкциите за употреба / 請參閱使用說明
-  Temperature range(s) / Plage(s) de température / Temperaturbereich(e) / Intervallo/i di temperatura / Intervalo(s) de temperatura / Intervalo(s) de temperatura / Temperaturintervall / Εύρος(-η) θερμοκρασίας / 溫度範圍 / Temperatūros diapazonas (-ai) / Hőmérséklet-tartomány(ok) / Zakres(y) temperatury / Rozsahy teplot / Rozsah(y) teploty / 온도 범위 / Sicaklık aralıkları / Диапазон(-ы) температуры / Температурен(ни) диапазон(и) / 溫度範圍
-  Caution / Précaution / Achtung / Attenzione / Precaución / Atenção / Försiktighet / Προσοχή / 注意事项 / Įspėjimas / Figyelem / Uwaga / Urozornění / Urozornenie / 주의 / Dikkat / Внимание / 注意
-  Expiration Date / Date D'expiration / Verfallsdatum, Verw. bis: / Data Di Scadenza / Fecha De Caducidad / Data de validade / Utgångsdatum / Ημερομηνία λήξης / 失效日期 / Galiojimo data / Lejárati idő / Data ważności / Datum expirace / Datum expirácie / 만료 날짜 / Son Kullanma Tarihi / Срок годности / Срок на годност / 到期日
- LOT** Lot Number / Numéro de lot / Chargennummer / Numero di lotto / Lote número / Número de lote / Satsnummer / Αριθ. παρτίδας / 批次号 / partijos numeris / Tételszám / Numer serii / Číslo šarže / 로트 번호 / Lot Numarası / Номер партии / Номер на партида / 批號
-  Date of Manufacture / Date de Fabrication / Herstellungsdatum / Data di Fabbricazione / Fecha de Fabricación / Data de Fabrico / Produktionsdatum / Ημερομηνία Παραγωγής / 生产日期 / Pagaminimo Data / Gyártás Dátuma / Data Produkcji / Datum Výroby / Datum Výroby / 제조 일자 / Üretim Tarihi / Дата Производства / Дата на Производство / 製造日期

Biohazard / Risque biologique / Biogefährdung / Rischio biologico / Riesgo biológico / Risco biológico / Biologisk fara / Βιολογικός κίνδυνος / 生物危害 / Biologisk fara / Veszélyes biológiai anyag / Zagrożenie biologiczne / Biologické riziko / Biologické riziko / 생물학적 위험 / Biolojik tehlike / Биологическая опасность / Биологична опасност / 生物危害

Radioactive / Radioactif / Radioaktiv / Radioattivo / Radiactivo / Radioactivo / Radioaktivt / Ραδιενεργό / 放射性 / Radioaktyvioji medžiaga / Radioaktiv / Radioaktywny / Radioaktivní / Rádiaoaktívny / 방사성 / Radyoaktif / Радиоактивный / Радиоактивен / 具放射性

Tracer / Traceur / Tracer / Marcato / Trazador / Marcador / / Ανιχνευτής / 追踪剂 / Atsekamoji medžiaga / Nyomjelző / Znacznik / Radioindikátor / Indikátor (tracer) / 트레이서 / Tracer'lar / метка / Индикатор / 追蹤劑

Calibrator / Calibrateur / Kalibrator / Calibratore / Calibrador / Calibrador / Kalibrator / Βαθμονομητής / 校准品 / Kalibravimo medžiaga / Kalibrátor / Kalibrator / kalibrátor / Kalibrátor / 보정 물질 / Kalibratör / Калибратор / Калибратор / 校正液

Control / Contrôle / Kontrolle / Controlla / Control / Controllo / Kontrolle / Μάρτυρας / 质控品 / Kontrolliné / Kontroll / Kontrola / Kontrola / Kontrola / 청도관리 / Kontrol / Контроль / Контролна / 質控品

Tubes / tubes / Röhrchen / provette / tubos / Tubos de amostra / Provrör / σωληνάρια / 试管 / Mégintuvéltai / Csövek / Probówki / Zkumavky / Skúmavky / 튜브 / Túpler / пробирки / Епруветки / 試管

Instruction for Use / Mode d'emploi / Gebrauchsanweisung / Istruzioni per l'uso / Instrucciones de uso / Instruções de utilização / Bruksanvisning / Οδηγίες χρήσης / 使用说明 / Naudojimo instrukcija / Használati utasítás / Instrukcja użycia / Návod k použití / Návod na použitie / 사용 안내 / Kullanma Talimatı / Инструкции / Инструкции за употреба / 使用说明

18 September 2020

REFERENCES

1. Approved Guideline - Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures, EP17-A2. June 2012. Clinical and Laboratory Standards Institute.
2. Approved Guideline – Evaluation of Precision of Quantitative Measurement Procedures, EP05-A3. October 2014. Clinical and Laboratory Standards Institute.
3. Approved Guideline - Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach, EP6-A. April 2003. Clinical and Laboratory Standards Institute.
4. Approved Guideline - Interference Testing in Clinical Chemistry, EP07, 3rd.ed. April 2018. Clinical and Laboratory Standards Institute.