

Pregnancy Test for Camels

Our quantitative **Progesterone camel ELISA** is specially developed for measuring progesterone in camel samples. With the highest sensitivity and precision, it supports veterinarians, researchers and breeders in **diagnosing pregnancy**.

- Kit Size: 96-well
- Sample Type: serum, plasma
- Sample volume: 10 µl
- Calibrator Range: 0.25 – 20 ng/ml
- Total Assay Time: 2 h
- Analytical sensitivity: 0.13 ng/ml
- No extraction needed
- Ready-to-use reagents
- Species-specific
- Validated with *Camelus dromedarius* samples


PLEASE NOTE:

A **Camel control kit** (Cat. No. DEV88CC) is available for this ELISA and can be used for internal quality control.

Contact our product specialist for your individual offer!

Dr. Manfred Czapp

Tel. +49(0)431/71922-32
czapp@demeditec.de



**demeditec
diagnostics gmbh**

Lise-Meitner-Straße 2
D-24145 Kiel
Germany

Tel. +49 (0)431 / 71922-0
info@demeditec.com



www.demeditec.com